

E3-BSA

Description

Product name	E3-BSA
Source	Synthetic

Specifications

Format	Liquid
---------------	--------

Stability and Storage

Storage	Store at -20°C, avoid repeated freeze/ thaw cycle.
Buffer	0.02M PBS, PH7.2
Stability.	Aliquot and store at -20°C for 24 months.

General information

Function	<p>Estriol (E3) is a metabolite of estradiol and estrone. E3 plays an important role in regulating intrauterine fetal development and can affect the sensitivity of pregnant uteri to oxytocin. Is an important index for clinically monitor placental function and judging fetal development.</p>
-----------------	--

For research use only.