




Catalogue of IVD antibody products



PTM BIO is a high-tech enterprise driven by an innovative platform, founded by renowned scientists from both domestic and international backgrounds. With over ten years of experience in antibody development, the company has established a platform for the development of mouse monoclonal antibodies and highly efficient rabbit monoclonal antibodies. We have successfully developed and launched over 150 raw materials for antibodies and antigens used in in vitro diagnostics. These materials cover a wide range of applications, including Alzheimer's disease, myocardial series, tumor series, thyroid series, hormone series, and veterinary antibodies and antigens. The products are suitable for various technology platforms such as chemiluminescence, immunoturbidimetry, and lateral flow chromatography.

CONTENTS |

Hormone	01	Cardiac markers	11
Infectious diseases	02	Veterinary	14
Thyroid diseases	04	Liver diseases	15
Kidney diseases	05	Diabetes and Metabolism	16
Inflammation	06	Bone metabolism	16
Tumor Markers	08	Hypertension, rheumatoid, anemia	17
Preeclampsia	10	Alzheimer's disease	17
Serology	11	Blood clotting	17

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Anti-müllerian hormone (AMH)	ABT-105008		FLF, CLIA
	ABT-105009		
	ABT-105010	ABT-105055/ABT-105009	
	ABT-105011	ABT-105009/ABT-105010	
	ABT-105012	ABT-105008/ABT-105012	
	ABT-105013	ABT-105012/ABT-105011	
	ABT-105052	ABT-105053/ABT-105052	
	ABT-105053	ABT-105053/ABT-105054	
	ABT-105054		
Parathyroid hormone (PTH)	ABT-105033	ABT-105033/ABT-105034 ABT-105036/ABT-105033	LF, CLIA
	ABT-105034		
	ABT-105035		
	ABT-105036		
	* Recombinant antigen ABT-9008		
Progesterone	ABT-105023		LF, CLIA
	ABT-105024		
	* Progesterone-BSA antigen ABT-9048		
Testosterone	ABT-105040		LF, CLIA
	ABT-105041		
	* Testosterone-BSA antigen ABT-9049		
Hydrocortisone	ABT-105042		LF, CLIA
Sex hormone binding globuli (SHBG)	ABT-105014	ABT-105014/ABT-105015	LF, CLIA
	ABT-105015		
Growth hormone (GH)	ABT-105003	ABT-105003/ABT-105004	LF, CLIA
	ABT-105004	ABT-105004/ABT-105025	
	ABT-105025	ABT-105003/ABT-105026	
	ABT-105026		
Human chorionic gonadotropin β -subunit (β -HCG)	ABT-105043	ABT-105043/ABT-105046 ABT-105049/ABT-105044 ABT-105048/ABT-105047 ABT-105051/ABT-105043 ABT-105050/ABT-105049	LF, CLIA
	ABT-105044		
	ABT-105045		
	ABT-105046		
	ABT-105047		
	ABT-105048		
	ABT-105049		
	ABT-105050		
ABT-105051			
Follicle stimulating hormone (FSH)	ABT-105028	ABT-105028/ABT-105029	LF, CLIA
	ABT-105029	ABT-105028/ABT-105030	
	ABT-105030	ABT-105067/ABT-105030	
	ABT-105067		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Luteinizing hormone (LH)	ABT-105016		LF, CLIA
	ABT-105017	ABT-105016/ABT-105017	
	ABT-105018	ABT-105016/ABT-105018	
	ABT-105056	ABT-105056/ABT-105057	
	ABT-105057	ABT-105056/ABT-105058	
	ABT-105058		
Prolactin (PRL)	ABT-105037	ABT-105038/ABT-105039	FLF, CLIA
	ABT-105038	ABT-105038/ABT-105037	
	ABT-105068	ABT-105068/ABT-105038	
	ABT-105039		
Insulin-like growth factor 1 (IGF-1)	ABT-105031		FLF, CLIA
	ABT-105032	ABT-105031/ABT-105032	
* Recombinant antigen	ABT-9014		
Estradiol (E2)			
* E2-BSA antigen	ABT-9047		
Insulin like growth factor binding protein 1 (IGFBP-1)	ABT-105019		FLF, CLIA
	ABT-105020	ABT-105022/ABT-105020	
	ABT-105021	ABT-105022/ABT-105021	
	ABT-105022		
Insulin like growth factor binding protein 3 (IGFBP-3)	ABT-105001		FLF, CLIA
	ABT-105002	ABT-105001/ABT-105002	
	ABT-105059	ABT-105059/ABT-105060	
	ABT-105060	ABT-105061/ABT-105060	
	ABT-105061		
	* Recombinant antigen	ABT-9007	
Insulin like growth factor binding protein 4 (IGFBP-4)	ABT-105005	ABT-105005/ABT-105006	FLF, CLIA
	ABT-105006	ABT-105006/ABT-105007	
	ABT-105007	ABT-105007/ABT-105005	

Infectious diseases

Product name	Catalog#	Recommended antibody pairs	Applicable platform
SARS-CoV-2 Nucleocapsid	ABT-113001		LF, ELISA, FLF, CLIA
	ABT-113002	ABT-113001/ABT-113002	
	ABT-113062	ABT-113063/ABT-113065	
	ABT-113063	ABT-113062/ABT-113063	
	ABT-113064	ABT-113063/ABT-113064	
	ABT-113065		
Monkeypox Virus A44R antigen	ABT-9044		
Monkeypox Virus B14R antigen	ABT-9045		
Monkeypox Virus B10R antigen	ABT-9046		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Yellow Fever Virus	ABT-113019	ABT-113019/ABT-113020	LF, ELISA, FLF, CLIA
	ABT-113020		
	ABT-113021	ABT-113020/ABT-113021	
	ABT-113022	ABT-113021/ABT-113022	
Adenovirus Hexon	ABT-113047	ABT-113047/ABT-113048	LF, ELISA, CLIA
	ABT-113048	ABT-113047/ABT-113049	
	ABT-113049	ABT-113050/ABT-113049	
	ABT-113050		
Entamoeba histolytica	ABT-113027	ABT-113027/ABT-113028	LF, ELISA
	ABT-113028		
Clostridium difficile	ABT-113036	ABT-113036/ABT-113037	LF, ELISA
	ABT-113037		
Zika Virus	ABT-113038	ABT-113038/ABT-113039	LF, ELISA
	ABT-113039	ABT-113079/ABT-113038	
	ABT-113079		
M.Pneumonia	ABT-113066	ABT-113066/ABT-113067	LF, ELISA
	ABT-113067	ABT-113066/ABT-113068	
	ABT-113068		
Human immunodeficiency virus P24 (HIV P24)	ABT-113056	ABT-113061/ABT-113059	LF, ELISA
	ABT-113057		
	ABT-113058		
	ABT-113059		
	ABT-113060		
	ABT-113061		
Hepatitis C virus (HCV)	ABT-113053	ABT-103053/ABT-103054	LF, ELISA
	ABT-113054	ABT-103053/ABT-103055	
	ABT-113055		
Influenza A Virus (Flu A)	ABT-113005	ABT-113005/ABT-113006	LF, ELISA, FLF, CLIA
	ABT-113006		
	ABT-113007		
	ABT-113029		
	ABT-113030		
	ABT-113031		
	ABT-113032		
	ABT-113033		
Hepatitis B surface antigen (HBsAg)	ABT-113008	ABT-113008/ABT-113010	LF, ELISA, FLF, CLIA
	ABT-113009	ABT-113009/ABT-113008	
	ABT-113010	ABT-113010/ABT-113011	
	ABT-113011	ABT-113011/ABT-113013	
	ABT-113012	ABT-113012/ABT-113014	
	ABT-113013	ABT-113015/ABT-113016	

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Hepatitis B surface antigen (HBsAg)	ABT-113014		
	ABT-113015		
	ABT-113016		
	ABT-113023		
Hepatitis B core antigen (HBcAg)	ABT-113017	ABT-113018/ABT-113017	LF, ELISA, FLF, CLIA
	ABT-113018		
Influenza B virus (Flu B)	ABT-113003		
	ABT-113024	ABT-113003/ABT-113025	
	ABT-113025	ABT-113004/ABT-113026	LF, ELISA, FLF, CLIA
	ABT-113026	ABT-113004/ABT-113025	
	ABT-113004		

Thyroid diseases

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Thyroid stimulating hormone (TSH)	ABT-104001		
	ABT-104002		
	ABT-104003	ABT-104001/ABT-104003	
	ABT-104004	ABT-104002/ABT-104004	FLF, CLIA
	ABT-104010	ABT-104003/ABT-104002	
	ABT-104011	ABT-104004/ABT-104012	
	ABT-104012		
	ABT-104013		
	ABT-104014		
Thyroglobulin	ABT-104015	ABT-104015/ABT-104016	FLF, CLIA
	ABT-104016		
	ABT-104017	ABT-104016/ABT-104017	
	ABT-104023	ABT-104023/ABT-104025	
	ABT-104024	ABT-104024/ABT-104023	
ABT-104025			
Triiodothyronine (T3)	ABT-104008		
	ABT-104009		
	ABT-104026		FLF, CLIA
	ABT-104027		
★ T3-BSA antigen	ABT-9051		
Thyroxine (T4)	ABT-104018		FLF, CLIA
	ABT-104019		
	ABT-104020		
	ABT-104021		
	★ T4-BSA antigen	ABT-9052	

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Sex hormone binding globulin (SHBG)	ABT-105014	ABT-105014/ABT-105015	FLF, CLIA
	ABT-105015		
Thyroid peroxidase (TPO)	ABT-104005	ABT-104005/ABT-104007	FLF, CLIA
	ABT-104006		
	ABT-104007	ABT-104006/ABT-104022	
	ABT-104022	ABT-104007/ABT-104006	
★ Recombinant antigen	ABT-9036		

Kidney diseases

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Cystatin C (CG)	ABT-107001	ABT-107002/ABT-107001 ABT-107002/ABT-107003 ABT-107002/ABT-107018 ABT-107019/ABT-107025 ABT-107020/ABT-107025 ABT-107020/ABT-107026	FLF, CLIA
	ABT-107002		
	ABT-107003		
	ABT-107018		
	ABT-107019		
	ABT-107020		
	ABT-107021		
	ABT-107022		
	ABT-107023		
	ABT-107024		
	ABT-107025		
	ABT-107026		
Neutrophil gelatinase-associated lipocalin (NGAL)	ABT-107010	ABT-107012/ABT-107014 ABT-107012/ABT-107015 ABT-107012/ABT-107016 ABT-107012/ABT-107017 ABT-107010/ABT-107013 ABT-107031/ABT-107032	FLF, CLIA
	ABT-107011		
	ABT-107012		
	ABT-107013		
	ABT-107014		
	ABT-107015		
	ABT-107016		
	ABT-107017		
ABT-107031			
ABT-107032			
Retinol binding protein4 (RBP4)	ABT-107006	ABT-107006/ABT-107009	FLF, CLIA, IT
	ABT-107007	ABT-107007/ABT-107009	
	ABT-107008	ABT-107008/ABT-107007	
	ABT-107009		
β2-microglobulin (β2-MG)	ABT-107027	ABT-107027/ABT-107028	FLF, CLIA, IT
	ABT-107028		
Kidney Injury Molecule-1 (KIM-1)	ABT-107004	ABT-107004/ABT-107005	FLF, CLIA
	ABT-107005	ABT-107005/ABT-107004	

Product name	Catalog#	Recommended antibody pairs	Applicable platform
C-reactive protein (CRP)	ABT-106001		
	ABT-106002		
	ABT-106003		
	ABT-106004	ABT-106001/ABT-106002	
	ABT-106005	ABT-106003/ABT-106005	
	ABT-106037	ABT-106004/ABT-106005	LF, FLF, CLIA, IT
	ABT-106038	ABT-106037/ABT-106038	
	ABT-106039	ABT-106040/ABT-106042	
	ABT-106040	ABT-106040/ABT-106043	
	ABT-106041		
	ABT-106042		
Interferon- γ (γ -IFN)	ABT-106013		
	ABT-106014	ABT-106013/ABT-106014	
	ABT-106017	ABT-106019/ABT-106017	FLF, CLIA
	ABT-106018	ABT-106019/ABT-106018	
	ABT-106019	ABT-106019/ABT-106049	
	ABT-106049		
Interferon- α (IFN- α)	ABT-106015		
	ABT-106016	ABT-106016/ABT-106015	FLF, CLIA
	★ Recombinant antigen ABT-9028		
Tumor Necrosis Factor α (TNF- α)	ABT-106026		
	ABT-106027		
	ABT-106050	ABT-106026/ABT-106027	
	ABT-106051	ABT-106050/ABT-106051	FLF, CLIA
	ABT-106052	ABT-106053/ABT-106052	
	ABT-106053		
★ Recombinant antigen ABT-9024			
Interleukin 1 β (IL-1 β)	ABT-106055		
	ABT-106056		
	ABT-106057	ABT-106055/ABT-106057	FLF, CLIA
	ABT-106073	ABT-106073/ABT-106055	
★ Recombinant antigen ABT-9029			
Interleukin 2RA (IL-2RA)	ABT-106058		
	ABT-106059	ABT-106058/ABT-106059	FLF, CLIA
	★ Recombinant antigen ABT-9042		
Interleukin 4 (IL-4)	ABT-106033		
	ABT-106034		
	ABT-106070	ABT-106034/ABT-106033	FLF, CLIA
	ABT-106071	ABT-106070/ABT-106072	
	ABT-106072		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Interleukin 5 (IL-5)	ABT-106028		
	ABT-106029		
	ABT-106030	ABT-106029/ABT-106028	
	ABT-106031	ABT-106029/ABT-106032	FLF, CLIA
	ABT-106032	ABT-106031/ABT-106028	
★ Recombinant antigen	ABT-9034		
Interleukin 6 (IL-6)	ABT-106011		
	ABT-106020		
	ABT-106021	ABT-106021/ABT-106011	
	ABT-106022	ABT-106025/ABT-106022	FLF, CLIA
	ABT-106023	ABT-106021/ABT-106025	
	ABT-106024		
Interleukin 8 (IL-8)	ABT-106009		
	ABT-106010	ABT-106009/ABT-106010	FLF, CLIA
	ABT-106065	ABT-106065/ABT-106010	
Interleukin 10 (IL-10)	ABT-106060		
	ABT-106061	ABT-106060/ABT-106061	
	ABT-106062	ABT-106062/ABT-106064	FLF, CLIA
	ABT-106063	ABT-106063/ABT-106061	
	ABT-106064		
★ Recombinant antigen	ABT-9035		
Interleukin 15 (IL-15)	ABT-106069		
	ABT-106068	ABT-106069/ABT-106067	FLF, CLIA
	ABT-106067	ABT-106069/ABT-106066	
	ABT-106066		
Calprotectin	ABT-102081		
	ABT-102082		
	ABT-102083		
	ABT-102084	ABT-102083/ABT-102087	
	ABT-102085	ABT-102084/ABT-102087	
	ABT-102086	ABT-102081/ABT-102082	LF
	ABT-102112	ABT-102112/ABT-102114	
	ABT-102113	ABT-102113/ABT-102115	
	ABT-102114		
	ABT-102115		
Procalcitonin (PCT)	ABT-102087		
	ABT-106006	ABT-106006/ABT-106008	FLF, CLIA, IT
	ABT-106007	ABT-106007/ABT-106008	
	ABT-106008		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Carcinoembryonic antigen (CEA)	ABT-102021		
	ABT-102022		
	ABT-102023	ABT-102021/ABT-102022	
	ABT-102024	ABT-102023/ABT-102021	
	ABT-102025	ABT-102024/ABT-102025	FLF, CLIA, IT
	ABT-102026	ABT-102027/ABT-102024	
	ABT-102101	ABT-102024/ABT-102023	
	ABT-102102	ABT-102101/ABT-102102	
* Recombinant antigen	ABT-9018		
Neuron-specific enolase (NSE)	ABT-102066		
	ABT-102067		
	ABT-102068	ABT-102070/ABT-102068	
	ABT-102069	ABT-102070/ABT-102071	FLF, CLIA
	ABT-102070	ABT-102066/ABT-102071	
	ABT-102071		
Human epidermal growth factor receptor-2 (HER2)	ABT-102088		
	ABT-102089	ABT-102088/ABT-102089	FLF, CLIA
	ABT-9033		
* Recombinant antigen	ABT-9022		
Carbohydrate protein CA50	ABT-102078		
	ABT-102079	ABT-102078/ABT-102080	FLF, CLIA
	ABT-102080	ABT-102078/ABT-102079	
Carbohydrate protein CA12-5	ABT-102046		
	ABT-102047	ABT-102047/ABT-102048	
	ABT-102048	ABT-102047/ABT-102049	FLF, CLIA
	ABT-102049	ABT-102047/ABT-102050	
	ABT-102050	ABT-102046/ABT-102047	
Carbohydrate protein CA72-4	ABT-102098		
	ABT-102099	ABT-102099/ABT-102100	FLF, CLIA
	ABT-102100	ABT-102100/ABT-102098	
Carbohydrate protein CA19-9	ABT-102109		
	ABT-102110	ABT-102109/ABT-102110	FLF, CLIA
Carbohydrate protein CA24-2	ABT-102103		
	ABT-102104	ABT-102103/ABT-102104	FLF, CLIA
	ABT-102111	ABT-102111/ABT-102103	
Carbohydrate protein CA15-3	ABT-102038		
	ABT-102039		
	ABT-102040	ABT-102038/ABT-102040	
	ABT-102041	ABT-102040/ABT-102042	FLF, CLIA
	ABT-102042	ABT-102043/ABT-102021	
	ABT-102043		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Glial fibrillary acidic protein (GFAP)	ABT-102034	ABT-102034/ABT-102035	FLF, CLIA
	ABT-102035	ABT-102035/ABT-102037	
	ABT-102036	ABT-102036/ABT-102034	
	ABT-102037	ABT-102037/ABT-102036	
Fibroblast growth factor receptor 3 (FGF R3) ★ Recombinant antigen	ABT-102072	ABT-102072/ABT-102073	FLF, CLIA
	ABT-102073		
	ABT-9040		
Prostate-specific antigen (PSA)	ABT-102001	ABT-102002/ABT-102001 ABT-102001/ABT-102058 ABT-102059/ABT-102057	FLF, CLIA, IT
	ABT-102002		
	ABT-102057		
	ABT-102058		
	ABT-102059		
Pepsinogen I (PG I) ★ Recombinant antigen	ABT-102012	ABT-102012/ABT-102014 ABT-102013/ABT-102015 ABT-102014/ABT-102013 ABT-102015/ABT-102014 ABT-102105/ABT-102107 ABT-102108/ABT-102106	FLF, CLIA, IT
	ABT-102013		
	ABT-102014		
	ABT-102015		
	ABT-102105		
	ABT-102106		
	ABT-102107		
	ABT-102108		
	ABT-9017		
Pepsinogen II (PG II) ★ Recombinant antigen	ABT-102016	ABT-102016/ABT-102017 ABT-102017/ABT-102019 ABT-102017/ABT-102020 ABT-102020/ABT-102017	FLF, CLIA, IT
	ABT-102017		
	ABT-102018		
	ABT-102019		
	ABT-102020		
Gastrin-17	ABT-102090	ABT-102090/ABT-102093 ABT-102091/ABT-102097 ABT-102093/ABT-102094 ABT-102093/ABT-102095	FLF, CLIA
	ABT-102091		
	ABT-102092		
	ABT-102093		
	ABT-102094		
	ABT-102095		
	ABT-102096		
	ABT-102097		
Pro-gastrin-releasing peptide (proGRP)	ABT-102060	ABT-102061/ABT-102063 ABT-102062/ABT-102063 ABT-102060/ABT-102061 ABT-102064/ABT-102065	FLF, CLIA
	ABT-102061		
	ABT-102062		
	ABT-102063		
	ABT-102064		
	ABT-102065		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Squamous Cell Carcinoma Antigen (SCCA)	ABT-102051	ABT-102052/ABT-102053	FLF, CLIA
	ABT-102052		
	ABT-102053	ABT-102054/ABT-102052	
	ABT-102054		
Cytokeratin 19 fragment (Cyfra 21-1)	ABT-102044	ABT-102044/ABT-102027	FLF, CLIA
	ABT-102045	ABT-102044/ABT-102045	
	ABT-102027		
Epididymal prote 4 (HE4)	ABT-102028	ABT-102028/ABT-102029	FLF, CLIA
	ABT-102029	ABT-102029/ABT-102031	
	ABT-102030	ABT-102029/ABT-102033	
	ABT-102031	ABT-102030/ABT-102031	
	ABT-102032	ABT-102032/ABT-102033	
	ABT-102033		
Alpha-fetoprotein (AFP)	ABT-102005	ABT-102005/ABT-102007	FLF, CLIA
	ABT-102006		
	ABT-102007	ABT-102005/ABT-102009	
	ABT-102008	ABT-102008/ABT-102010	
	ABT-102009	ABT-102009/ABT-102011	
	ABT-102010	ABT-102010/ABT-102005	
	ABT-102011		
* Recombinant antigen	ABT-9005		
Abnormal prothrombin protein Induced by vitamin K absence or antagonist-II (PIVKA-II)	ABT-102003	ABT-102004/ABT-102003 ABT-102075/ABT-102077 ABT-102076/ABT-102077	LF, CLIA
	ABT-102004		
	ABT-102075		
	ABT-102076		
	ABT-102077		
* Recombinant antigen	ABT-9006		
Vascular endothelial growth factor			LF, CLIA
* Recombinant antigen VEGF121	ABT-9038		
* Recombinant antigen VEGF165	ABT-9039		

Preeclampsia

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Soluble vascular endothelial growth factor receptor 1 (sFlt-1)	ABT-111001	ABT-111003/ABT-111004 ABT-111003/ABT-111001 ABT-111001/ABT-111002	LF, CLIA
	ABT-111002		
	ABT-111003		
	ABT-111004		
* Recombinant antigen	ABT-9013		
Placental growth factor (PLGF)	ABT-111005	ABT-111005/ABT-111006	LF, CLIA
	ABT-111006		


Serology

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Immunoglobulin A (IgA)	ABT-116006	ABT-116007/ABT-116006	FLF, CLIA
	ABT-116007		
Immunoglobulin G (IgG)	ABT-116009		FLF, CLIA
	ABT-116010		
Immunoglobulin M (IgM)	ABT-116011		FLF, CLIA
	ABT-116012		
	ABT-116013		
Immunoglobulin E (IgE)	ABT-116001	ABT-116001/ABT-116003	FLF, CLIA
	ABT-116002	ABT-116002/ABT-116003	
	ABT-116003	ABT-116004/ABT-116005	
	ABT-116004	ABT-116005/ABT-116001	
	ABT-116005		


Cardiac markers

Product name	Catalog#	Recommended antibody pairs	Applicable platform
B type natriuretic peptide (BNP)	ABT-101027	ABT-101027/ABT-101028	FLF, CLIA
	ABT-101028		
	ABT-101029		
	ABT-101030		
	ABT-101094		
	ABT-101095		
	ABT-101096		
	ABT-101097		
	ABT-101098		
	ABT-101099		
Creatine Kinase MB (CK-MB)	ABT-101082	ABT-101083/ABT-101082	LF, CLIA
	ABT-101083	ABT-101088/ABT-101083	
	ABT-101088		
N terminal pro B type natriuretic peptide (NT-proBNP)	ABT-101040	ABT-101044/ABT-101041	LF, CLIA
	ABT-101041		
	ABT-101042		
	ABT-101043		
	ABT-101044		
	ABT-101045		
	ABT-101110		
	ABT-101102		
	ABT-101103		
	ABT-101104		
ABT-101105			

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Growth differentiation factor-15 (GDF-15)	ABT-101049	ABT-101049/ABT-101050 ABT-101049/ABT-101051	LF, CLIA
	ABT-101050		
	ABT-101051		
Procalcitonin (PCT)	ABT-106006	ABT-106007/ABT-106008 ABT-106007/ABT-106045 ABT-106044/ABT-106007	LF, CLIA
	ABT-106007		
	ABT-106008		
	ABT-106044		
	ABT-106045		
	ABT-106046		
	ABT-106047		
ABT-106048			
Fibrinogen (FIB)	ABT-101072	ABT-101072/ABT-101073	LF, CLIA
	ABT-101073		
D-Dimer	ABT-101011	ABT-101011/ABT-101012 ABT-101014/ABT-101012 ABT-101015/ABT-101023 ABT-101011/ABT-101013 ABT-101047/ABT-101023 ABT-101047/ABT-101090 ABT-101090/ABT-101091 ABT-101106/ABT-101012	LF, FLF, CLIA, IT
	ABT-101012		
	ABT-101013		
	ABT-101014		
	ABT-101015		
	ABT-101023		
	ABT-101046		
	ABT-101047		
	ABT-101048		
	ABT-101106		
	ABT-101089		
	ABT-101090		
	ABT-101091		
Soluble growth stimulating expression of gene 2 protein (sST2)	ABT-101052	ABT-101053/ABT-101055 ABT-101056/ABT-101057 ABT-101061/ABT-101060	LF, CLIA
	ABT-101053		
	ABT-101054		
	ABT-101055		
	ABT-101056		
	ABT-101057		
	ABT-101058		
	ABT-101059		
	ABT-101060		
ABT-101061			
Galectin-3	ABT-101018	ABT-101018/ABT-101020	FLF, CLIA
	ABT-101019	ABT-101019/ABT-101018	
	ABT-101020	ABT-101020/ABT-101021	
	ABT-101021	ABT-101022/ABT-101020	
	ABT-101022		

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Calcium-binding protein B (S100B)	ABT-101001		
	ABT-101002	ABT-101036/ABT-101002	
	ABT-101035	ABT-101035/ABT-101002	FLF, CLIA
	ABT-101036	ABT-101001/ABT-101002	
	* Recombinant antigen	ABT-9032	
Myoglobin (MYO)	ABT-101003		
	ABT-101004	ABT-101004/ABT-101003	FLF, CLIA, IT
	* Recombinant antigen	ABT-9004	
Myeloperoxidase (MPO)	ABT-101005	ABT-101005/ABT-101006	
	ABT-101006	ABT-101005/ABT-101007	
	ABT-101007	ABT-101007/ABT-101005	FLF, CLIA
	ABT-101074	ABT-101074/ABT-101101	
	ABT-101101		
Cardiac troponin I (cTnI)	ABT-101008		
	ABT-101009		
	ABT-101010		
	ABT-101107		
	ABT-101062	ABT-101009/ABT-101065	
	ABT-101063	ABT-101010/ABT-101066	
	ABT-101064	ABT-101010/ABT-101008	LF, FLF, CLIA, IT
	ABT-101065	ABT-101107/ABT-101109	
	ABT-101066	ABT-101067/ABT-101062	
	ABT-101067		
ABT-101109			
ABT-101108			
Cardiac troponin T (cTnT)	ABT-101068	ABT-101069/ABT-101068	
	ABT-101069	ABT-101068/ABT-101070	FLF, CLIA
	ABT-101070	ABT-101070/ABT-101071	
	ABT-101071		
Lipoprotein-associated phospholipase A2 (Lp-PLA2)	ABT-101031		
	ABT-101032	ABT-101031/ABT-101076	
	ABT-101033	ABT-101032/ABT-101076	
	ABT-101075	ABT-101076/ABT-101077	FLF, CLIA
	ABT-101076	ABT-101076/ABT-101078	
	ABT-101077	ABT-101033/ABT-101076	
	ABT-101078		
Fatty acid-binding protein (FABP)	ABT-101024		
	ABT-101025	ABT-101024/ABT-101026	
	ABT-101026	ABT-101025/ABT-101026	FLF, CLIA, IT
	* Recombinant antigen	ABT-9003	

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Adiponectin	ABT-101037		FLF, CLIA
	ABT-101038	ABT-101039/ABT-101037	
	ABT-101039	ABT-101038/ABT-101039	
	★ Recombinant antigen ABT-9037		

Veterinary

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Canine distemper virus (CDV)	ABT-117004		LF, CLIA
	ABT-117005	ABT-117004/ABT-117006	
	ABT-117006	ABT-117005/ABT-117006	
	★ Recombinant antigen ABT-9021		
Canine N-terminal pro-B type natriuretic peptide (CNT-proBNP)	ABT-117012		LF, CLIA
	ABT-117013	ABT-117013/ABT-117012	
	ABT-117014	ABT-117014/ABT-117013	
	ABT-117015	ABT-117014/ABT-117015	
	ABT-117036	ABT-117036/ABT-117037	
	ABT-117037		
Canine C-reactive protein (CCRP)	ABT-117007		LF, CLIA
	ABT-117008	ABT-117007/ABT-117008	
	ABT-117009	ABT-117009/ABT-117010	
	ABT-117010	ABT-117007/ABT-117009	
	ABT-117011		
	★ Recombinant antigen ABT-9010		
Canine parvovirus (CPV)	ABT-117025		LF, CLIA
	ABT-117026	ABT-117025/ABT-117027	
	ABT-117027	ABT-117026/ABT-117028	
	ABT-117028	ABT-117041/ABT-117042	
	ABT-117041	ABT-117043/ABT-117042	
	ABT-117042		
	ABT-117043		
Feline serum amyloid A (FSAA)	ABT-117016		LF, CLIA
	ABT-117017		
	ABT-117018	ABT-117016/ABT-117018	
	ABT-117019	ABT-117016/ABT-117019	
	ABT-117020	ABT-117016/ABT-117020	
	ABT-117032	ABT-117020/ABT-117033	
	ABT-117033		
★ Recombinant antigen ABT-9027			

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Feline leukemia virus (FeLV)	ABT-117034	ABT-117034/ABT-117035	LF, CLIA
	ABT-117035		
	★ Recombinant antigen ABT-9043		
Feline parvovirus (FPV)	ABT-117022	ABT-117023/ABT-117024	LF, CLIA
	ABT-117023	ABT-117022/ABT-117023	
	ABT-117024	ABT-117038/ABT-117039	
	ABT-117038		
	ABT-117039		
Feline Cystatin C (fCys-C)	ABT-9026		LF, CLIA
★ Recombinant antigen			
Canine Coronavirus (CCV)	ABT-117045	ABT-117045/ABT-117046	LF, CLIA
	ABT-117046		
	★ Recombinant antigen ABT-9025		
African Swine Fever Virus (ASFV)	ABT-117001	ABT-117001/ABT-117002	LF, ELISA, CLIA
	ABT-117002	ABT-117002/ABT-117003	
	ABT-117003		
	★ Recombinant antigen ABT-9011		

Liver diseases

Product name	Catalog#	Recommended antibody pairs	Applicable platform
N-terminal peptide of type III Procollagen (PIIINP)	ABT-109003	ABT-109003/ABT-109004	LF, CLIA
	ABT-109004		
	★ Recombinant antigen ABT-9019		
	★ Recombinant antigen ABT-9020		
Collagen IV (C-IV)	ABT-109009	ABT-109009/ABT-109010	LF, CLIA
	ABT-109010		
	★ Recombinant antigen ABT-9031		
Laminin (LN)	ABT-109011	ABT-109012/ABT-109011	LF, CLIA
	ABT-109012		
Serum cholyglycin (CG)	ABT-109005		LF, CLIA
	ABT-109006		
	ABT-109007		
	ABT-109008		
	ABT-109016		
Chitinase 3-like protein 1 (CHI3L1)	ABT-9041		LF, CLIA
★ Recombinant antigen			
Prealbumin	ABT-109013	ABT-109013/ABT-109015	LF, CLIA
	ABT-109014	ABT-109013/ABT-109014	
	ABT-109015	ABT-109014/ABT-109017	
	ABT-109017		

Diabetes and Metabolism

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Ferritin	ABT-103007		
	ABT-103008		
	ABT-103015	ABT-103008/ABT-103007	
	ABT-103016	ABT-103007/ABT-103015	FLF, CLIA
	ABT-103029	ABT-103029/ABT-103030	
	ABT-103030		
Transferrin	ABT-103001		
	ABT-103002	ABT-103001/ABT-103028	FLF, CLIA, IT
	ABT-103028		
Transferrin receptor (sTR)	ABT-103003		
	ABT-103004	ABT-103003/ABT-103005	
	ABT-103005	ABT-103004/ABT-103006	FLF, CLIA
	ABT-103006		
Glycated hemoglobin (HbA1c)	ABT-103021		
	ABT-103022	ABT-103022/ABT-103024	LF, CLIA
	ABT-103023		
Hemoglobin	ABT-103024		
	ABT-103025	ABT-103025/ABT-103027	
	ABT-103026	ABT-103027/ABT-103026	LF, CLIA
	ABT-103027		
C-Peptide	ABT-103018		
	ABT-103019	ABT-103018/ABT-103019	
	ABT-103020	ABT-103018/ABT-103020	LF, CLIA
Albumin	ABT-103009		
	ABT-103010	ABT-103009/ABT-103010	
	ABT-103011	ABT-103010/ABT-103011	LF, CLIA, IT
	ABT-103012	ABT-103010/ABT-103012	
Lactoferrin	ABT-103013		
	ABT-103014	ABT-103013/ABT-103014	FLF, CLIA, IT

Bone metabolism

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Osteocalcin (OST)	ABT-108003		
	ABT-108004	ABT-108003/ABT-108004	LF, CLIA
	ABT-108005		
	ABT-108006		
N-terminal peptide of type I Procollagen (PINP)	ABT-108001		
	ABT-108002	ABT-108002/ABT-108001	LF, CLIA
	ABT-9030		

* Recombinant antigen

Hypertension, rheumatoid, anemia

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Renin	ABT-110001		
	ABT-110002	ABT-110002/ABT-110001	
	ABT-110003	ABT-110001/ABT-110003	FLF, CLIA
	ABT-110004	ABT-110004/ABT-110005	
	ABT-110005		
* Recombinant antigen	ABT-9009		
Folic acid (FA)	ABT-110009		FLF, CLIA
Aldosterone	ABT-110010		FLF, CLIA
Adrenocorticotrophic hormone (ACTH)	ABT-110007	ABT-110007/ABT-110008	FLF, CLIA
	ABT-110008		
Matrix metalloproteinase-3 (MMP-3)	ABT-110011	ABT-110011/ABT-110013	
	ABT-110012	ABT-110012/ABT-110014	FLF, CLIA
	ABT-110013	ABT-110013/ABT-110014	
	ABT-110014		
Angiotensin II (All)	ABT-110006		FLF, CLIA
* All-BSA antigen	ABT-9050		

Alzheimer's disease

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Beta Amyloid 1-42 (A β 1-42)	ABT-112002		
	ABT-112003		
	ABT-112004	ABT-112002/ABT-112005	
	ABT-112005	ABT-112003/ABT-112005	CLIA
	ABT-112010	ABT-112010/ABT-112011	
	ABT-112011		
	ABT-112012		
Brain-derived neurotrophic factor (BDNF)	ABT-112006	ABT-112006/ABT-112007	LF, CLIA
	ABT-112007		
Phosphorylated tau 217 (p-Tau 217)	ABT-112013		CLIA
Phosphorylated tau 181 (p-Tau 181)	ABT-112001	ABT-112001/ABT-112009	CLIA
Total Tau Protein (T-Tau)	ABT-112008		CLIA
	ABT-112009		

Blood clotting

Product name	Catalog#	Recommended antibody pairs	Applicable platform
Thrombomodulin	ABT-119001		
	ABT-119002	ABT-119002/ABT-119003	
	ABT-119003	ABT-119001/ABT-119002	LF, CLIA
	ABT-9015		
* Recombinant antigen			



PTM BIO LLC

Telephone: +1-312-802-6843

Email: info@ptmbio.com

Address: 2201 West Campbell Park Drive, Suite 305, Chicago, IL 60612

Website: www.ptmbio.com