




PTMab Antibody Products

PTMab Antibodies

Experts in PTM-Specific Antibodies

Driven by Innovation, Dedicated to PTM Excellence

2024	D-lactylation	「Nature Chemical Biology」
2021	Methacrylation	「Cell Discovery」
2019	L-lactylation 🔥	「Nature」
2018	Benzoylation	「Nature Communications」
2016	β-Hydroxybutyrylation	「Molecular Cell」
2014	2-Hydroxyisobutyrylation	「Nature Chemical Biology」
2014	Glutarylation	「Cell Metabolism」
2012	Malonylation	「Cell」
2011	Crotonylation	「Molecular & Cellular Proteomics」
2010	Succinylation	「Nature Chemical Biology」
2007	Propionylation	「Molecular & Cellular Proteomics」
2007	Butyrylation	「Molecular & Cellular Proteomics」

PTM BIO PTM-Specific Antibody Advantages



Broad PTM Coverage and Comprehensive Site Selection

50+ pan-PTM antibodies
330+ histone site-specific PTM antibodies



Dedicated to PTM Antibodies, Committed to Quality

Each PTM-specific antibody is validated by at least two specificity assays

Validation methods include dot blot, blocking peptide assay, drug treatment, and IP-MS



Broad PTM Coverage with Rigorous Validation

Validated for WB, IHC, IP, ICC/IF, FC, and ChIP

Verified across multiple sample types, species, and negative controls

Application and species recommendations are based on experimental validation data



PTM BIO PTM-Specific Antibody Advantages

High specificity

High affinity

High batch-to-batch consistency

Contents

Pan-PTM Antibodies

- P01 Pan-PTM Antibodies
- P03 Pan-PTM Antibody Kits (PTM-6680/6681)
- P05 Pan-PTM Antibody Resins

Histone Site-Specific PTM Antibodies

- P05 Histone Site-Specific PTM Antibodies
- P16 Histone Lactylation / Acetylation Antibody Kits (PTM-7093/7094)
- P17 Histone H3K18 Site-Specific PTM Antibody Kit (PTM-7477)
- P18 Popular Histone H3 Modification / Site Antibody Kit (PTM-7480)

Epigenetics Antibodies

- | | |
|----------------------------------|---|
| P19 Histone Antibodies | P35 Cancer Research Antibodies |
| P19 Histone Acyltransferases | P55 Immunology Research Antibodies |
| P19 Histone Deacylases | P70 Cardiovascular Research Antibodies |
| P20 Histone Methyltransferases | P77 Metabolism Research Antibodies |
| P21 Histone Demethylases | P88 Developmental & Stem Cell Research Antibodies |
| P21 DNA Methyltransferases | P93 Neuroscience Research Antibodies |
| P21 RNA Methyltransferases | P106 Cell Biology Research Antibodies |
| P21 RNA Demethylases | P125 Loading Control Antibodies |
| P21 Chromatin Remodeling Factors | P126 Isotype Controls |
| P22 Transcription Factors | P126 Secondary Antibodies |
| P26 Others | P127 Tag Antibodies |

Contents

Pan-PTM Antibodies

Cat.#	Product	Recommended Applications
PTM-3201RM	Anti-Pyruvyl Lysine Rabbit mAb	WB
PTM-1429	Anti-D-Lactyl Lysine Rabbit pAb	WB
PTM-1429RM	Anti-D-Lactyl Lysine Rabbit mAb	WB/ICC/IF/IP
PTM-1401RM	Anti-L-Lactyl Lysine Rabbit mAb	WB/IHC-P/IP/ICC/IF/FC/ChIP
PTM-1401	Anti-L-Lactyl Lysine Rabbit pAb	WB/FC
PTM-1425	Anti-L-Lactyl Lysine Rabbit pAb-BSA and Azide free	WB/IHC-P
PTM-102	Anti-Acetyllysine Mouse mAb	WB/ICC/IF
PTM-101	Anti-Acetyllysine Mouse mAb	WB/IHC-P
PTM-105RM	Anti-Acetyllysine Rabbit mAb	WB/IHC-P
PTM-132	Anti-Acetyllysine Rabbit pAb-BSA and Azide free	WB/FC
PTM-502	Anti-Crotonyllysine Mouse mAb	WB
PTM-501	Anti-Crotonyllysine Rabbit pAb	WB
PTM-419	Anti-Succinyllysine Mouse mAb	WB
PTM-401	Anti-Succinyllysine Rabbit pAb	WB
PTM-802	Anti-2-Hydroxyisobutyryllysine Mouse mAb	WB
PTM-801	Anti-2-Hydroxyisobutyryllysine Rabbit pAb	WB
PTM-1201RM	Anti-β-Hydroxybutyryllysine Rabbit mAb	WB/ICC/IF/IHC-P
PTM-329	Anti-Butyryllysine Mouse mAb	WB/ICC/IF
PTM-301RM	Anti-Butyryllysine Rabbit mAb	WB/IHC-P/ICC/IF
PTM-301	Anti-Butyryllysine Rabbit pAb	WB
PTM-337RM	Anti-Butyryllysine Rabbit mAb-BSA and Azide free	WB
PTM-902	Anti-Malonyllysine Mouse mAb	WB
PTM-901	Anti-Malonyllysine Rabbit pAb	WB
PTM-203	Anti-Propionyllysine Mouse mAb	WB/IHC-P
PTM-201	Anti-Propionyllysine Rabbit pAb	WB
PTM-1152	Anti-Glutaryllysine Mouse mAb	WB
PTM-1151	Anti-Glutaryllysine Rabbit pAb	WB
PTM-762	Anti-Benzoyllysine Mouse mAb	WB/IHC-P
PTM-1501	Anti-Methacryllysine Mouse mAb	WB/IHC-P
PTM-1601	Anti-Acetylmethionine Rabbit pAb	WB/ICC/IF/FC
PTM-1701RM	Anti-Carboxyethyllysine Rabbit mAb	WB
PTM-606	Anti-Di-Methyllysine Rabbit pAb	WB
PTM-606F	Anti-Di-Methyllysine Rabbit mAb-BSA and Azide free	WB
PTM-602	Anti-Mono/Di-Methyllysine Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-601	Anti-Tri-Methyllysine Rabbit pAb	WB
PTM-5798	Anti-Ubiquitin Mouse mAb	WB/ICC/IF
PTM-7205	Anti-Ubiquitin (linkage-specific K48) Rabbit mAb	WB/IHC-P/FC
PTM-1107	Anti-Ubiquitin Mouse mAb (NT)	WB/IHC-P
PTM-7228	Anti-Ubiquitin (linkage-specific K63) Rabbit mAb	WB
PTM-7311	Anti-Ubiquitin Rabbit mAb	WB/IHC-P
PTM-1124RM	Anti-Ubiquitin Rabbit mAb-BSA and Azide free	WB
PTM-1106RM	Anti-Ubiquitin Rabbit mAb	WB/ICC/IF
PTM-1109	Anti-SUMO1/2/3 Mouse mAb	WB/IHC-P
PTM-617RM	Anti-Symmetric Di-Methyl Arginine Rabbit mAb	WB/IHC-P/ICC/IF
PTM-617	Anti-Symmetric Di-Methyl Arginine Rabbit pAb	WB/ICC/IF
PTM-605RM	Anti-Asymmetric Di-Methyl Arginine Rabbit mAb	WB
PTM-701	Anti-Phospho-Tyrosine Mouse mAb	WB
PTM-702RM	Anti-Phospho-Tyrosine Rabbit mAb	WB/IHC-P/IP
PTM-729RM	Anti-Phospho-Tyrosine Rabbit mAb-BSA and Azide free	WB
PTM-702	Anti-Phospho-Tyrosine Rabbit pAb	WB
PTM-705RM	Anti-Phospho-Threonine Rabbit mAb	WB/ICC/IF
PTM-730RM	Anti-Phospho-Threonine Rabbit mAb-BSA and Azide free	WB
PTM-951RM	Anti-O-Linked N-Acetylglucosamine Rabbit mAb	WB/ICC/IF/IHC-P
PTM-955RM	Anti-O-Linked N-Acetylglucosamine Rabbit mAb-BSA and Azide free	WB
PTM-752	Anti-3-Nitrotyrosine Mouse mAb	WB/ICC/IF
PTM-751	Anti-3-Nitrotyrosine Rabbit pAb	WB

All-in-One Acylation Antibody Kits

「 Start Smart, Validate with Confidence 」

「 Test More, Waste Less 」

「 Two Kits, One Flexible Solution 」

「 From Initial Screening to Deeper Discovery 」

Acylation Antibody Kit Plus

Cat.No.: PTM-6680

Product	Cat.#	Product Size
Anti-Acetyllysine Rabbit mAb	PTM-105RM	50 µl
Anti-L-Lactyl Lysine Rabbit mAb	PTM-1401RM	20 µl
Anti-Succinyllysine Mouse mAb	PTM-419	20 µl
Anti-Crotonyllysine Mouse mAb	PTM-502	20 µl
Anti-β-Hydroxybutyryllysine Rabbit mAb	PTM-1201RM	20 µl
Anti-Malonyllysine Mouse mAb	PTM-902	20 µl
Anti-2-Hydroxyisobutyryllysine Rabbit pAb	PTM-802	20 µl
Anti-Propionyllysine Mouse mAb	PTM-203	20 µl
Anti-Butyryllysine Rabbit mAb	PTM-301RM	20 µl
Anti-Benzoyllysine Mouse mAb	PTM-762	20 µl
Anti-Glutaryllysine Mouse mAb	PTM-1152	20 µl

Acylation Antibody Kit

Cat. No.: PTM-6681

Product	Cat.#	Product Size
Anti-Acetyllysine Rabbit mAb	PTM-105RM	50 µl
Anti-L-Lactyl Lysine Rabbit mAb	PTM-1401RM	20 µl
Anti-Succinyllysine Mouse mAb	PTM-419	20 µl
Anti-Crotonyllysine Mouse mAb	PTM-502	20 µl
Anti-β-Hydroxybutyryllysine Rabbit mAb	PTM-1201RM	20 µl
Anti-Malonyllysine Mouse mAb	PTM-902	20 µl
Anti-2-Hydroxyisobutyryllysine Rabbit pAb	PTM-802	20 µl

PTM Type	Biological Relevance
Lactylation First reported in <i>Nature</i> , 2019	Closely associated with glycolysis and mitochondrial oxidative metabolism. Involved in key biological processes such as tumor development, inflammation, metabolism, and hypoxic stress.
Succinylation First reported in 2010	Widely involved in metabolic regulation. Many enzymes in central metabolism are succinylated, and this modification is closely associated with cancer, cardiovascular disease, and inflammation.
Crotonylation Recognized as a major highlight in epigenetics by <i>Cell</i> , 2011	Regulates gene expression, stem cell differentiation, reproductive development, tumor metabolism, DNA damage repair, and other biological processes.
β-Hydroxybutyrylation First reported in 2016	Closely linked to core cellular metabolism, including fatty acid oxidation. Involved in energy metabolism, DNA damage repair, and tumorigenesis.
2-Hydroxyisobutyrylation First reported in 2014	Associated with amino acid biosynthesis and glycolysis. Involved in gene expression regulation, germ cell differentiation, and pathogenicity of plant-pathogenic fungi.
Phosphorylation A classic and widely studied PTM	Regulates enzymatic activity and ligand–receptor interactions, playing important roles in signal transduction, development, differentiation, cancer progression, and stress responses.
Acetylation A classic PTM with a long research history	Affects chromatin structure and gene expression, and participates in epigenetic regulation, signal transduction, metabolic regulation, immune response, tumorigenesis, and apoptosis.
Methylation A classic PTM with a long research history	Regulates cancer, aging, neurodegenerative diseases, and other processes through epigenetic mechanisms, including gene expression control, protein function regulation, and RNA processing.
N-Acetylation Acetylation of the first residue at the protein N-terminus	Regulates protein stability, protein–protein interactions, and membrane association, and plays an important role in autophagy.
Malonylation First reported in 2012	Involved in cellular energy metabolism, metabolic disorders, immune regulation, and plant stress responses.
Benzoylation First reported in 2018	Participates in the regulation of signaling pathways, hormone secretion, epigenetic regulation, and metabolic control.
Glutarylation First reported in 2014	Regulates chromatin-related functions, including nucleosome assembly, gene expression, DNA damage repair, and cell cycle progression.
Methacrylation First reported in 2021	An isomeric modification of crotonylation, with potential regulatory roles in epigenetics and cellular detoxification. It may also be associated with neurological and liver diseases.
Carboxyethylation Non-enzymatic chemical modification	Involved in oxidative damage regulation and associated with atherosclerosis, diabetes, uremia, heart failure, neurodegenerative diseases, and other conditions.
O-GlcNAcylation O-linked glycosylation of serine/threonine residues	Affects protein folding, conformation, localization, stability, and activity. It participates in multiple signaling pathways and regulates development, hormone response, cellular stress, transcriptional control, signal transduction, metabolic diseases, and tumor progression.
Ubiquitination A classic PTM involved in protein degradation	Mediates protein degradation and regulates substrate function by altering protein localization and activity. It plays important roles in cell cycle control, apoptosis, aging, neurodegenerative diseases, and signal transduction.
SUMOylation A common ubiquitin-like modification	Participates in key physiological and biochemical processes, including the regulation of protein–protein and protein–DNA interactions. It is involved in cardiovascular disease, cancer, neurodegenerative diseases, immune regulation, and other pathological processes.


Application Strategy: Immunoprecipitation Followed by Western Blotting for Single-Protein PTM Profiling

Case Study: Post-Translational Modification of Viral Proteins

> *PLoS Pathog.* 2023 Apr 6;19(4):e1011324. doi: 10.1371/journal.ppat.1011324.
eCollection 2023 Apr.

Immune evasion strategy involving propionylation by the KSHV interferon regulatory factor 1 (vIRF1)

Research Context | Viruses rely on host-cell translation machinery and may acquire host-mediated post-translational modifications. In this study, cells were transfected with the KSHV vIRF1 gene. The vIRF1 protein was enriched by immunoprecipitation, and modification signals were detected using a panel of pan-PTM antibodies. Strong acetylation and propionylation signals were observed, demonstrating the utility of pan-PTM antibodies for profiling PTMs on viral proteins.



Pan-PTM Antibodies Resins

Cat.#	Product	Recommended Applications
PTM-804	Anti-2-hydroxyisobutyryllysine antibody conjugated agarose beads	IAC
PTM-104	Anti-acetyllysine antibody conjugated agarose beads	IAC
PTM-503	Anti-crotonyllysine antibody conjugated agarose beads	IAC
PTM-1104	Anti-diglycine lysine antibody conjugated agarose beads	IAC
PTM-1154	Anti-glutaryllysine antibody conjugated agarose beads	IAC
PTM-1404	Anti-L-lactyllysine antibody conjugated agarose beads	IAC
PTM-1434	Anti-D-Lactyl Lysine Antibody Conjugated Agarose beads	IAC
PTM-904	Anti-malonyllysine antibody conjugated agarose beads	IAC
PTM-703	Anti-phosphotyrosine antibody conjugated agarose beads	IAC
PTM-202	Anti-propionyllysine antibody conjugated agarose beads	IAC
PTM-402	Anti-succinyllysine antibody conjugated agarose beads	IAC
PTM-1204	Anti-β-hydroxybutyryllysine antibody conjugated agarose beads	IAC
/	beads buffer kit	For use with beads

Histone Site-Specific PTM Antibodies

Histone Lysine Lactylation

PTM-1422RM	Anti-L-Lactyl-Histone H2A.Z (Lys11) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP
PTM-1426RM	Anti-L-Lactyl-Histone H2B (Lys15) Rabbit mAb	WB/ChIP
PTM-1424RM	Anti-L-Lactyl-Histone H2B (Lys16) Rabbit mAb	WB
PTM-1423	Anti-L-Lactyl-Histone H2B (Lys120) Rabbit pAb	WB
PTM-1419RM	Anti-L-Lactyl-Histone H3 (Lys9) Rabbit mAb	WB/IHC-P/ChIP
PTM-1414RM	Anti-L-Lactyl-Histone H3 (Lys14) Rabbit mAb	WB/ChIP/ICC/IF
PTM-1414	Anti-L-Lactyl-Histone H3 (Lys14) Rabbit pAb	WB/ChIP
PTM-1406RM	Anti-L-Lactyl-Histone H3 (Lys18) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-1427RM	Anti-L-Lactyl-Histone H3 (Lys18) Rabbit mAb-ChIP Grade	WB/ChIP/CUT&Tag
PTM-1413RM	Anti-L-Lactyl-Histone H3 (Lys23) Rabbit mAb	WB/IHC-P/ChIP
PTM-1428	Anti-L-Lactyl-Histone H3 (Lys27) Rabbit pAb	WB/ChIP
PTM-1421RM	Anti-L-Lactyl-Histone H3 (Lys56) Rabbit mAb	WB/FC
PTM-1409	Anti-L-Lactyl-Histone H4 (Lys5) Mouse mAb	WB/FC
PTM-1407RM	Anti-L-Lactyl-Histone H4 (Lys5) Rabbit mAb	WB/ICC/IF/FC/IP/CUT&Tag
PTM-1407	Anti-L-Lactyl-Histone H4 (Lys5) Rabbit pAb	WB
PTM-1415RM	Anti-L-Lactyl-Histone H4 (Lys8) Rabbit mAb	WB/ChIP/CUT&Tag

Cat.#	Product	Recommended Applications
PTM-1411RM	Anti-L-Lactyl-Histone H4 (Lys12) Rabbit mAb	WB/IHC-P/IP/ChIP
PTM-1417RM	Anti-L-Lactyl-Histone H4 (Lys16) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP
Histone Lysine Acetylation		
PTM-177	Anti-Acetyl-Histone H1.4 (Lys25) Rabbit pAb	WB
PTM-106	Anti-Acetyl-Histone H2A (Lys5) Rabbit pAb	WB
PTM-106RM	Anti-Acetyl-Histone H2A (Lys5) Rabbit mAb	WB/ChIP
PTM-193	Anti-Acetyl-Histone H2A (Lys9) Mouse mAb	WB
PTM-194	Anti-Acetyl-Histone H2A (Lys9) Rabbit mAb	WB/ChIP/IP
PTM-175	Anti-Acetyl-Histone H2A (Lys15) Mouse mAb	WB/IHC-P
PTM-172	Anti-Acetyl-Histone H2A (Lys15) Rabbit pAb	WB
PTM-173	Anti-Acetyl-Histone H2A (Lys118) Rabbit pAb	WB
PTM-152	Anti-Acetyl-Histone H2B (Lys5) Mouse mAb	WB/IHC-P/ChIP/IP
PTM-107	Anti-Acetyl-Histone H2B (Lys5) Rabbit pAb	WB/IHC-P
PTM-176	Anti-Acetyl-Histone H2B (Lys11) Mouse mAb	WB/IHC-P/IP/ChIP/ICC/IF
PTM-130	Anti-Acetyl-Histone H2B (Lys11) Rabbit pAb	WB
PTM-153	Anti-Acetyl-Histone H2B (Lys12) Mouse mAb	WB/IHC-P/ChIP
PTM-108	Anti-Acetyl-Histone H2B (Lys12) Rabbit pAb	WB/ChIP
PTM-181	Anti-Acetyl-Histone H2B (Lys15) Mouse mAb	WB/IHC-P
PTM-109	Anti-Acetyl-Histone H2B (Lys15) Rabbit pAb	WB/IHC-P
PTM-167	Anti-Acetyl-Histone H2B (Lys16) Mouse mAb	WB/IHC-P/ICC/IF
PTM-123	Anti-Acetyl-Histone H2B (Lys16) Rabbit pAb	WB/IHC-P/ChIP
PTM-155	Anti-Acetyl-Histone H2B (Lys20) Mouse mAb	WB/IHC-P/ChIP/ICC/IF
PTM-110	Anti-Acetyl-Histone H2B (Lys20) Rabbit pAb	WB/IHC-P/ChIP
PTM-174	Anti-Acetyl-Histone H2B (Lys23) Mouse mAb	WB
PTM-171	Anti-Acetyl-Histone H2B (Lys23) Rabbit pAb	WB/ChIP
PTM-169	Anti-Acetyl-Histone H2B (Lys24) Mouse mAb	WB/ChIP/ICC/IF
PTM-126	Anti-Acetyl-Histone H2B (Lys24) Rabbit pAb	WB
PTM-182	Anti-Acetyl-Histone H2B (Lys46) Rabbit pAb	WB
PTM-182RM	Anti-Acetyl-Histone H2B (Lys46) Rabbit mAb	WB/ChIP/ICC/IF
PTM-131	Anti-Acetyl-Histone H2B (Lys116) Rabbit pAb	WB/ChIP
PTM-111	Anti-Acetyl-Histone H2B (Lys120) Rabbit pAb	WB/ChIP
PTM-168	Anti-Acetyl-Histone H3 (Lys4) Mouse mAb	WB/ICC/IF
PTM-188	Anti-Acetyl-Histone H3 (Lys4) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-124	Anti-Acetyl-Histone H3 (Lys4) Rabbit pAb	WB/ChIP

Cat.#	Product	Recommended Applications
PTM-156	Anti-Acetyl-Histone H3 (Lys9) Mouse mAb	WB/ChIP/IP
PTM-112RM	Anti-Acetyl-Histone H3 (Lys9) Rabbit mAb	WB/IHC-P/ChIP/CUT&Tag
PTM-112	Anti-Acetyl-Histone H3 (Lys9) Rabbit pAb	WB/IHC-P/IP
PTM-157	Anti-Acetyl-Histone H3 (Lys14) Mouse mAb	WB/IHC-P/ICC/IF/ChIP/IP
PTM-113RM	Anti-Acetyl-Histone H3 (Lys14) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP
PTM-158	Anti-Acetyl-Histone H3 (Lys18) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-114RM	Anti-Acetyl-Histone H3 (Lys18) Rabbit mAb	WB/IHC-P/ChIP
PTM-115	Anti-Acetyl-Histone H3 (Lys23) Rabbit pAb	WB
PTM-115RM	Anti-Acetyl-Histone H3 (Lys23) Rabbit mAb	WB/IHC-P/ChIP/IP
PTM-160	Anti-Acetyl-Histone H3 (Lys27) Mouse mAb	WB/IHC-P/ICC/IF
PTM-116RM	Anti-Acetyl-Histone H3 (Lys27) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP/IP
PTM-117	Anti-Acetyl-Histone H3 (Lys36) Rabbit pAb	WB/IHC-P/ICC/IF/ChIP
PTM-117RM	Anti-Acetyl-Histone H3 (Lys36) Rabbit mAb	WB/IHC-P/ChIP/ICC/IF
PTM-128	Anti-Acetyl-Histone H3 (Lys37) Rabbit pAb	WB
PTM-162	Anti-Acetyl-Histone H3 (Lys56) Mouse mAb	WB/ICC/IF/ChIP
PTM-129RM	Anti-Acetyl-Histone H3 (Lys64) Rabbit mAb	WB/IHC-P/ChIP
PTM-129	Anti-Acetyl-Histone H3 (Lys64) Rabbit pAb	WB/ChIP
PTM-170	Anti-Acetyl-Histone H3 (Lys115) Rabbit pAb	WB
PTM-184RM	Anti-Acetyl-Histone H3 (Lys122) Rabbit mAb	WB/IHC-P
PTM-184	Anti-Acetyl-Histone H3 (Lys122) Rabbit pAb	WB
PTM-163	Anti-Acetyl-Histone H4 (Lys5) Mouse mAb	WB/IHC-P/ICC/IF
PTM-119	Anti-Acetyl-Histone H4 (Lys5) Rabbit pAb	WB/IHC-P/ICC/IF/ChIP
PTM-189	Anti-Acetyl-Histone H4 (Lys5/8/12/16) Mouse mAb	WB/IHC-P/ICC/IF
PTM-190	Anti-Acetyl-Histone H4 (Lys8) Mouse mAb	WB/IHC-P/ChIP
PTM-164	Anti-Acetyl-Histone H4 (Lys8) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-120	Anti-Acetyl-Histone H4 (Lys8) Rabbit pAb	WB/IHC-P/ICC/IF
PTM-192	Anti-Acetyl-Histone H4 (Lys12) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-121RM	Anti-Acetyl-Histone H4 (Lys12) Rabbit mAb	ELISA/WB/IHC-P
PTM-121	Anti-Acetyl-Histone H4 (Lys12) Rabbit pAb	WB/IHC-P/ICC/IF
PTM-187	Anti-Acetyl-Histone H4 (Lys16) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-122	Anti-Acetyl-Histone H4 (Lys16) Rabbit pAb	WB/IHC-P/ICC/IF/ChIP
PTM-127RM	Anti-Acetyl-Histone H4 (Lys77) Rabbit mAb	WB/IHC-P/IP/ChIP
Histone Lysine Crotonylation		
PTM-543	Anti-Crotonyl-Histone H2A (Lys119) Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-505	Anti-Crotonyl-Histone H2A (Lys119) Rabbit pAb	WB/ChIP
PTM-508RM	Anti-Crotonyl-Histone H2B (Lys11) Rabbit mAb	WB/IHC-P
PTM-508	Anti-Crotonyl-Histone H2B (Lys11) Rabbit pAb	WB
PTM-529	Anti-Crotonyl-Histone H2B (Lys11) Mouse mAb	WB
PTM-528	Anti-Crotonyl-Histone H2B (Lys12) Mouse mAb	WB/IHC-P/ChIP
PTM-509	Anti-Crotonyl-Histone H2B (Lys12) Rabbit pAb	WB/ChIP
PTM-546RM	Anti-Crotonyl-Histone H2B (Lys15/Lys16/Lys20) Rabbit mAb	WB/ChIP
PTM-533	Anti-Crotonyl-Histone H2B (Lys16) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-534	Anti-Crotonyl-Histone H2B (Lys20) Mouse mAb	WB/ChIP
PTM-512	Anti-Crotonyl-Histone H2B (Lys20) Rabbit pAb	WB/IP
PTM-514	Anti-Crotonyl-Histone H2B (Lys34) Rabbit pAb	WB/ChIP
PTM-527	Anti-Crotonyl-Histone H3 (Lys4) Mouse mAb	WB/ChIP
PTM-515RM	Anti-Crotonyl-Histone H3 (Lys4) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-547RM	Anti-Crotonyl-Histone H3 (Lys4) Rabbit mAb	WB/IHC-P/ChIP/IP
PTM-539	Anti-Crotonyl-Histone H3 (Lys9) Mouse mAb	WB/ChIP
PTM-516RM	Anti-Crotonyl-histone H3 (Lys9) Rabbit mAb	WB/IHC-P/ICC/IF/IP/ChIP
PTM-516	Anti-Crotonyl-Histone H3 (Lys9) Rabbit pAb	WB/ChIP
PTM-537	Anti-Crotonyl-Histone H3 (Lys14) Mouse mAb	WB/ICC/IF/ChIP/IP
PTM-535RM	Anti-Crotonyl-Histone H3 (Lys14) Rabbit mAb	WB/IHC-P/ChIP
PTM-535	Anti-Crotonyl-Histone H3 (Lys14) Rabbit pAb	WB/IP
PTM-540	Anti-Crotonyl-Histone H3 (Lys18) Mouse mAb	WB/IHC-P/ChIP
PTM-517	Anti-Crotonyl-Histone H3 (Lys18) Rabbit pAb	WB/IHC-P
PTM-517RM	Anti-Crotonyl-Histone H3 (Lys18) Rabbit mAb	WB/ICC/IF
PTM-544RM	Anti-Crotonyl-histone H3 (Lys18/ Lys23) Rabbit mAb	WB/ICC/IF/IP
PTM-519	Anti-Crotonyl-Histone H3 (Lys23) Mouse mAb	WB/IHC-P/ChIP
PTM-545RM	Anti-Crotonyl-Histone H3 (Lys27) Rabbit mAb	WB/ICC/IF/ChIP
PTM-536	Anti-Crotonyl-Histone H3 (Lys36) Rabbit pAb	WB/ChIP
PTM-521	Anti-Crotonyl-Histone H4 (Lys5) Rabbit pAb	WB
PTM-522RM	Anti-Crotonyl-Histone H4 (Lys8) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP
PTM-542	Anti-Crotonyl-Histone H4 (Lys8) Rabbit mAb	WB/IHC-P
PTM-530	Anti-Crotonyl-Histone H4 (Lys12) Mouse mAb	WB/ChIP
PTM-523	Anti-Crotonyl-Histone H4 (Lys12) Rabbit pAb	WB
PTM-524	Anti-Crotonyl-Histone H4 (Lys16) Mouse mAb	WB/IP/IHC-P
PTM-532	Anti-Crotonyl-Histone H4 (Lys16) Rabbit pAb	WB
PTM-541	Anti-Crotonyl-Histone H4 (Ser1) Rabbit pAb	WB/ICC/IF/ChIP

Cat.#	Product	Recommended Applications
-------	---------	--------------------------

Histone Lysine 2-hydroxyisobutyrylation

PTM-820	Anti-2-Hydroxyisobutyryl-Histone H2A (Lys5) Rabbit pAb	WB
PTM-871	Anti-2-Hydroxyisobutyryl-Histone H2B (Lys12) Mouse mAb	WB/ChIP
PTM-831	Anti-2-Hydroxyisobutyryl-Histone H2B (Lys12) Rabbit pAb	WB
PTM-881	Anti-2-Hydroxyisobutyryl-Histone H3 (Lys14) Mouse mAb	WB
PTM-841RM	Anti-2-Hydroxyisobutyryl-Histone H3 (Lys14) Rabbit mAb	WB
PTM-882	Anti-2-Hydroxyisobutyryl-Histone H3 (Lys18) Mouse mAb	WB/IHC-P/ChIP/ICC/IF
PTM-843RM	Anti-2-Hydroxyisobutyryl-Histone H3 (Lys23) Rabbit mAb	WB/ICC/IF/ChIP
PTM-845	Anti-2-Hydroxyisobutyryl-Histone H3 (Lys79) Rabbit pAb	WB
PTM-854	Anti-2-Hydroxyisobutyryl-Histone H4 (Lys5) Mouse mAb	WB/ChIP
PTM-850	Anti-2-Hydroxyisobutyryl-Histone H4 (Lys5/8/12) Rabbit pAb	WB
PTM-805	Anti-2-Hydroxyisobutyryl-Histone H4 (Lys8) Rabbit pAb	WB/ICC/IF
PTM-808	Anti-2-Hydroxyisobutyryl-Histone H4 (Lys20) Rabbit pAb	WB

Histone Lysine β -hydroxybutyrylation

PTM-1270	Anti- β -Hydroxybutyryl-Histone H2A (Lys5) Mouse mAb	WB
PTM-1220	Anti- β -Hydroxybutyryl-Histone H2A (Lys5) Rabbit pAb	WB
PTM-1259	Anti- β -Hydroxybutyryl-Histone H2A/H2A.X (Lys9) Rabbit pAb	WB
PTM-1224	Anti- β -Hydroxybutyryl-Histone H2A (Lys118) Rabbit pAb	WB
PTM-1225	Anti- β -Hydroxybutyryl-Histone H2A.X (Lys134) Rabbit pAb	WB
PTM-1230	Anti- β -Hydroxybutyryl-Histone H2B (Lys5) Rabbit pAb	WB
PTM-1231	Anti- β -Hydroxybutyryl-Histone H2B (Lys11) Rabbit pAb	WB
PTM-1232	Anti- β -Hydroxybutyryl-Histone H2B (Lys12) Rabbit pAb	WB
PTM-1233RM	Anti- β -Hydroxybutyryl-Histone H2B (Lys15) Rabbit mAb	WB
PTM-1234	Anti- β -Hydroxybutyryl-Histone H2B (Lys16) Rabbit pAb	WB
PTM-1235	Anti- β -Hydroxybutyryl-Histone H2B (Lys20) Rabbit pAb	WB
PTM-1236	Anti- β -Hydroxybutyryl-Histone H2B (Lys23) Rabbit pAb	WB
PTM-1238	Anti- β -Hydroxybutyryl-Histone H2B (Lys34) Rabbit pAb	WB
PTM-1296	Anti- β -Hydroxybutyryl-Histone H3 (Lys4) Mouse mAb	WB/IP
PTM-1258	Anti- β -Hydroxybutyryl-Histone H3 (Lys4) Rabbit pAb	WB/IP
PTM-1260	Anti- β -Hydroxybutyryl-Histone H3 (Lys4) Rabbit pAb	WB/ChIP
PTM-1250RM	Anti- β -Hydroxybutyryl-Histone H3 (Lys9) Rabbit mAb	WB/IHC-P
PTM-1250	Anti- β -Hydroxybutyryl-Histone H3 (Lys9) Rabbit pAb	WB
PTM-1292	Anti- β -Hydroxybutyryl-Histone H3 (Lys18) Mouse mAb	WB
PTM-1252	Anti- β -Hydroxybutyryl-Histone H3 (Lys18) Rabbit pAb	WB
PTM-1300	Anti- β -Hydroxybutyryl-Histone H3 (Lys23) Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-1293	Anti- β -Hydroxybutyryl-Histone H3 (Lys27) Mouse mAb	WB
PTM-1257	Anti- β -Hydroxybutyryl-Histone H3 (Lys27) Rabbit pAb	WB
PTM-1205	Anti- β -Hydroxybutyryl-Histone H4 (Lys5) Rabbit pAb	WB
PTM-1253RM	Anti- β -Hydroxybutyryl-Histone H4 (Lys8) Rabbit mAb	WB/ICC/IF
PTM-1206RM	Anti- β -Hydroxybutyryl-Histone H4 (Lys12) Rabbit mAb	WB/ICC/IF/ChIP
PTM-1206	Anti- β -Hydroxybutyryl-Histone H4 (Lys12) Rabbit pAb	WB/ChIP
PTM-1262	Anti- β -Hydroxybutyryl-Histone H4 (Lys16) Mouse mAb	WB/ChIP

Histone Lysine Succinylation

PTM-409	Anti-Succinyl-Histone H2B (Lys120) Rabbit pAb	WB
PTM-421	Anti-Succinyl-Histone H3 (Lys14) Mouse mAb	WB/IHC-P
PTM-422	Anti-Succinyl-Histone H3 (Lys23) Mouse mAb	WB/IHC-P
PTM-7479	Anti-Succinyl-Histone H3 (Lys23) Mouse mAb	WB/IHC-P
PTM-412	Anti-Succinyl-Histone H3 (Lys79) Rabbit pAb	WB
PTM-413	Anti-Succinyl-Histone H3 (Lys122) Rabbit pAb	WB

Histone Lysine Butyrylation

PTM-327	Anti-Butyryl-Histone H2A (Lys5) Mouse mAb	WB/ChIP
PTM-317	Anti-Butyryl-Histone H2A (Lys5) Rabbit pAb	WB/ICC/IF/IHC-P
PTM-320	Anti-Butyryl-Histone H2B (Lys5) Mouse mAb	WB/IHC-P/ChIP
PTM-303	Anti-Butyryl-Histone H2B (Lys5) Rabbit pAb	WB/ICC/IF/IHC-P
PTM-323	Anti-Butyryl-Histone H2B (Lys12) Mouse mAb	WB
PTM-318	Anti-Butyryl-Histone H2B (Lys12) Rabbit pAb	WB
PTM-324	Anti-Butyryl-Histone H2B (Lys15) Mouse mAb	WB/ICC/IF
PTM-333	Anti-Butyryl-Histone H2B (Lys15) Rabbit pAb	WB/ICC/IF/ChIP
PTM-322	Anti-Butyryl-Histone H2B (Lys16) Mouse mAb	WB/ICC/IF/ChIP
PTM-319	Anti-Butyryl-Histone H2B (Lys16) Rabbit pAb	WB/ICC/IF/IHC-P
PTM-316	Anti-Butyryl-Histone H2B (Lys20) Mouse mAb	WB/IHC-P
PTM-304	Anti-Butyryl-Histone H2B (Lys20) Rabbit pAb	WB/ICC/IF/IHC-P
PTM-335RM	Anti-Butyryl-Histone H2B (Lys23) Rabbit mAb	WB/ChIP
PTM-312	Anti-Butyryl-Histone H3 (Lys9) Mouse mAb	WB/ChIP/ICC/IF
PTM-305	Anti-Butyryl-Histone H3 (Lys9) Rabbit pAb	WB/ICC/IF/ChIP
PTM-334RM	Anti-Butyryl-Histone H3 (Lys14) Rabbit mAb	WB/ICC/IF/ChIP
PTM-331	Anti-Butyryl-Histone H3 (Lys18) Mouse mAb	WB/ChIP
PTM-306RM	Anti-Butyryl-Histone H3 (Lys18) Rabbit mAb	WB/ICC/IF
PTM-338RM	Anti-Butyryl-Histone H3 (Lys18) Rabbit mAb	WB/ChIP/IP

Cat.#	Product	Recommended Applications
PTM-307	Anti-Butyryl-Histone H3 (Lys23) Mouse mAb	WB/IHC-P
PTM-326	Anti-Butyryl-Histone H3 (Lys27) Mouse mAb	WB/ICC/IF
PTM-315	Anti-Butyryl-Histone H3 (Lys27) Rabbit pAb	WB/ChIP
PTM-336	Anti-Butyryl-Histone H3 (Lys115) Rabbit pAb	WB
PTM-310	Anti-Butyryl-Histone H4 (Lys5) Mouse mAb	WB/IHC-P/ChIP
PTM-313	Anti-Butyryl-Histone H4 (Lys5) Rabbit pAb	WB
PTM-311	Anti-Butyryl-Histone H4 (Lys8) Rabbit pAb	WB/ICC/IF/ChIP
PTM-314	Anti-Butyryl-Histone H4 (Lys12) Mouse mAb	WB

Histone Lysine Malonylation

PTM-911RM	Anti-Malonyl-Histone H3 (Lys23) Rabbit mAb	WB/ChIP/ICC/IF
PTM-910RM	Anti-Malonyl-Histone H3 (Lys122) Rabbit mAb	WB/IHC-P
PTM-910	Anti-Malonyl-Histone H3 (Lys122) Rabbit pAb	WB

Histone Lysine Propionylation

PTM-217	Anti-Propionyl-Histone H2A (Lys125) Rabbit pAb	WB
PTM-217RM	Anti-Propionyl-Histone H2A (Lys125) Rabbit mAb	WB
PTM-215	Anti-Propionyl-Histone H3 (Lys14) Mouse mAb	WB
PTM-211	Anti-Propionyl-Histone H3 (Lys14) Rabbit pAb	WB
PTM-214	Anti-Propionyl-Histone H3 (Lys18) Mouse mAb	WB
PTM-213	Anti-Propionyl-Histone H3 (Lys18) Rabbit pAb	WB
PTM-208	Anti-Propionyl-Histone H3 (Lys23) Mouse mAb	WB/IHC-P/ChIP
PTM-205	Anti-Propionyl-Histone H3 (Lys23) Rabbit pAb	WB/IHC-P/ChIP
PTM-220RM	Anti-Propionyl-Histone H3 (Lys56) Rabbit mAb	WB/IHC-P/IP/ChIP
PTM-216	Anti-Propionyl-Histone H4 (Lys5) Mouse mAb	WB/IHC-P
PTM-219RM	Anti-Propionyl-Histone H4 (Lys5) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-218	Anti-Propionyl-Histone H4 (Lys8) Mouse mAb	WB/IHC-P/ChIP
PTM-209	Anti-Propionyl-Histone H4 (Lys12) Mouse mAb	WB/IHC-P
PTM-206	Anti-Propionyl-Histone H4 (Lys12) Rabbit pAb	WB/ICC/IF/ChIP
PTM-210	Anti-Propionyl-Histone H4 (Lys16) Mouse mAb	WB/IHC-P/ChIP

Histone Lysine Ubiquitination/SUMOylation

PTM-1121	Anti-Ubiquityl-Histone H2A (Lys119) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP/IP
PTM-1119	Anti-Ubiquityl-Histone H2A (Lys119) Rabbit pAb	WB
PTM-7449	Anti-Ubiquityl-Histone H2A (Lys119) Rabbit mAb	WB/ICC/IF/IP/ChIP
PTM-1120	Anti-Ubiquityl-Histone H2B (Lys120) Mouse mAb	WB
PTM-1122	Anti-Ubiquityl-Histone H2B (Lys120) Mouse mAb	WB/IHC-P
PTM-1123RM	Anti-Ubiquityl-Histone H2B (Lys120) Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
-------	---------	--------------------------

Histone Arginine Citrullination

PTM-1305	Anti-Citrulline-Histone H2A (Arg3) Rabbit pAb	WB
PTM-1306RM	Anti-Citrulline-Histone H3 (Arg2/8/17) Rabbit mAb	WB/IHC-P

Histone Lysine Methacrylation

PTM-1503	Anti-Methacryl-Histone H3 (Lys18) Mouse mAb	WB
----------	---	----

Histone Lysine Mono-Methylation

PTM-698	Anti-Mono-Methyl-Histone H2B (Lys12) Rabbit pAb	WB/ChIP
PTM-661	Anti-Mono-Methyl-Histone H2B (Lys15) Mouse mAb	WB/IHC-P/IP/ChIP
PTM-5158	Anti-Mono-Methyl-Histone H3 (Lys4) Mouse mAb	WB/IHC-P/ChIP
PTM-611RM	Anti-Mono-Methyl-Histone H3 (Lys4) Rabbit mAb	WB/ICC/IF/IP
PTM-611	Anti-Mono-Methyl-Histone H3 (Lys4) Rabbit pAb	WB/IHC-P/ICC/IF
PTM-614	Anti-Mono-Methyl-Histone H3 (Lys9) Rabbit pAb	WB/ICC/IF/IP
PTM-7287	Anti-Mono-Methyl-Histone H3 (Lys9) Rabbit mAb	WB/ChIP
PTM-670	Anti-Mono-Methyl-Histone H3 (Lys14) Rabbit pAb	WB
PTM-618	Anti-Mono-Methyl-Histone H3 (Lys18) Rabbit pAb	WB/ChIP
PTM-618RM	Anti-Mono-Methyl-Histone H3 (Lys18) Rabbit mAb	WB
PTM-619RM	Anti-Mono-Methyl-Histone H3 (Lys23) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP/IP
PTM-619	Anti-Mono-Methyl-Histone H3 (Lys23) Rabbit pAb	WB/IP/ChIP
PTM-5110	Anti-Mono-Methyl-Histone H3 (Lys27) Mouse mAb	WB/ChIP
PTM-649	Anti-Mono-Methyl-Histone H3 (Lys27) Mouse mAb	WB/ChIP
PTM-620	Anti-Mono-Methyl-Histone H3 (Lys27) Rabbit pAb	WB
PTM-620RM	Anti-Mono-Methyl-Histone H3 (Lys27) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP
PTM-623RM	Anti-Mono-methyl-histone H3 (Lys36) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-623	Anti-Mono-Methyl-Histone H3 (Lys36) Rabbit pAb	WB/ICC/IF
PTM-655	Anti-Mono-Methyl-Histone H3 (Lys56) Mouse mAb	WB/IP
PTM-626	Anti-Mono-Methyl-Histone H3 (Lys56) Rabbit pAb	WB/ICC/IF/ChIP
PTM-5112	Anti-Mono-Methyl-Histone H3 (Lys79) Rabbit mAb	WB/IHC-P
PTM-628	Anti-Mono-Methyl-Histone H3 (Lys79) Rabbit pAb	WB/IHC-P/ICC/IF/IP
PTM-628RM	Anti-Mono-Methyl-Histone H3 (Lys79) Rabbit mAb	WB/ChIP/IP
PTM-631RM	Anti-Mono-Methyl-Histone H3 (Lys122) Rabbit mAb	WB
PTM-631	Anti-Mono-Methyl-Histone H3 (Lys122) Rabbit pAb	WB/ChIP
PTM-685RM	Anti-Mono-Methyl-Histone H4 (Lys12) Rabbit mAb	WB/IHC-P/ICC/IF/FC/ChIP/IP
PTM-5111	Anti-Mono-Methyl-Histone H4 (Lys20) Mouse mAb	WB/ICC/IF/ChIP
PTM-634	Anti-Mono-Methyl-Histone H4 (Lys20) Rabbit pAb	WB/ChIP
PTM-637	Anti-Mono-Methyl-Histone H4 (Lys31) Rabbit pAb	WB

Cat.#	Product	Recommended Applications
-------	---------	--------------------------

Histone Lysine Di-Methylation

PTM-641	Anti-Di-Methyl-Histone H3 (Lys4) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-6032	Anti-Di-Methyl-Histone H3 (Lys4) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-612	Anti-Di-Methyl-Histone H3 (Lys4) Rabbit pAb	WB/ICC/IF/ChIP
PTM-644	Anti-Di-Methyl-Histone H3 (Lys9) Mouse mAb	WB/IHC-P/ICC/IF/IP
PTM-5003	Anti-Di-Methyl-Histone H3 (Lys9) Mouse mAb	WB/IHC-P/ICC/IF
PTM-615RM	Anti-Di-Methyl-Histone H3 (Lys9) Rabbit mAb	WB/ICC/IF/ChIP/IP
PTM-640RM	Anti-Di-Methyl-Histone H3 (Lys14) Rabbit mAb	WB/IHC-P/ICC/IF/ChIP
PTM-676	Anti-Di-Methyl-Histone H3 (Lys23) Mouse mAb	WB/IP/IHC-P/ChIP/ICC/IF
PTM-5312	Anti-Di-Methyl-Histone H3 (Lys27) Mouse mAb	WB/IHC-P/ChIP
PTM-5010	Anti-Di-Methyl-Histone H3 (Lys27) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5156	Anti-Di-Methyl-Histone H3 (Lys36) Mouse mAb	WB/IHC-P
PTM-662	Anti-Di-Methyl-Histone H3 (Lys36) Mouse mAb	WB/IHC-P/ChIP
PTM-624	Anti-Di-Methyl-Histone H3 (Lys36) Rabbit pAb	WB/IHC-P/ICC/IF
PTM-687	Anti-Di-Methyl-Histone H3 (Lys37) Rabbit pAb	WB
PTM-627RM	Anti-Di-Methyl-Histone H3 (Lys56) Rabbit mAb	WB/IHC-P/ICC/IF/IP/ChIP
PTM-627	Anti-Di-Methyl-Histone H3 (Lys56) Rabbit pAb	WB/ICC/IF
PTM-5159	Anti-Di-Methyl-Histone H3 (Lys79) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-658	Anti-Di-Methyl-Histone H3 (Lys79) Mouse mAb	WB/IHC-P/ICC/IF/IP
PTM-629	Anti-Di-Methyl-Histone H3 (Lys79) Rabbit pAb	WB/IHC-P/IP
PTM-645RM	Anti-Di-Methyl-Histone H3 (Lys122) Rabbit mAb	WB/ICC/IF/IP/ChIP

Histone Lysine Tri-Methylation

PTM-5019	Anti-Tri-Methyl-Histone H3 (Lys4) Mouse mAb	WB/ChIP
PTM-613	Anti-Tri-Methyl-Histone H3 (Lys4) Rabbit pAb	WB/IHC-P/ChIP
PTM-616	Anti-Tri-Methyl-Histone H3 (Lys9) Rabbit pAb	WB/IHC-P/ICC/IF
PTM-616RM	Anti-Tri-Methyl-Histone H3 (Lys9) Rabbit mAb	WB/IHC-P/ChIP/IP
PTM-646	Anti-Tri-Methyl-Histone H3 (Lys14) Mouse mAb	WB
PTM-632	Anti-Tri-Methyl-Histone H3 (Lys18) Rabbit pAb	WB/ChIP
PTM-648	Anti-Tri-Methyl-Histone H3 (Lys23) Mouse mAb	WB
PTM-638	Anti-Tri-Methyl-Histone H3 (Lys23) Rabbit pAb	WB/ChIP
PTM-651	Anti-Tri-Methyl-Histone H3 (Lys27) Mouse mAb	WB/IHC-P/ICC/IF
PTM-5002	Anti-Tri-Methyl-Histone H3 (Lys27) Mouse mAb	WB/IHC-P/ICC/IF/ChIP
PTM-647RM	Anti-Tri-Methyl-Histone H3 (Lys27) Rabbit mAb	WB/ChIP/IP
PTM-625RM	Anti-Tri-Methyl-Histone H3 (Lys36) Rabbit mAb	WB/IHC-P
PTM-625	Anti-Tri-Methyl-Histone H3 (Lys36) Rabbit pAb	WB/IHC-P/ICC/IF/IP

Cat.#	Product	Recommended Applications
PTM-653	Anti-Tri-Methyl-Histone H3 (Lys36) Rabbit pAb	WB/ChIP
PTM-653RM	Anti-Tri-methyl-Histone H3 (Lys36) Rabbit mAb	WB/IHC-P/ChIP/IP
PTM-690	Anti-Tri-Methyl-Histone H3 (Lys64) Rabbit pAb	WB/ChIP
PTM-5157	Anti-Tri-Methyl-Histone H3 (Lys79) Mouse mAb	WB/IHC-P
PTM-689	Anti-Tri-Methyl-Histone H3 (Lys79) Mouse mAb	WB/ChIP
PTM-630	Anti-Tri-Methyl-Histone H3 (Lys79) Rabbit pAb	WB/IHC-P/IP
PTM-665	Anti-Tri-Methyl-Histone H4 (Lys20) Mouse mAb	WB
PTM-5005	Anti-Tri-Methyl-Histone H4 (Lys20) Mouse mAb	WB
PTM-636	Anti-Tri-Methyl-Histone H4 (Lys20) Rabbit pAb	WB/IHC-P

Histone Lysine Mono/Di/Tri-Methylation

PTM-677	Anti-Di/Tri-Methyl-Histone H3 (Lys23) Rabbit pAb	WB/IP
PTM-5006	Anti-Di/Tri-Methyl-Histone H4 (Lys20) Mouse mAb	WB/IHC-P
PTM-669	Anti-Mono/Di-Methyl-Histone H3 (Lys56) Mouse mAb	WB
PTM-657RM	Anti-Mono/Di-Methyl-Histone H3 (Lys79) Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP/ChIP
PTM-7435	Anti-Mono/Di/Tri-Methyl-Histone H3 (Lys79) Rabbit mAb	WB/IHC-P

Histone Lysine Tri-Methylation/ Histone Glutamine Serotonylation

PTM-1420	Anti-Histone H3 (tri methyl K4, serotonyl Q5) Rabbit pAb	WB/ChIP
----------	--	---------

Histone Lysine Tri-Methylation/ Histone Serine Phosphorylation

PTM-743	Anti-Histone H3 (tri methyl K9, phospho S10) Mouse mAb	WB
PTM-723	Anti-Histone H3 (tri methyl K9, phospho S10) Rabbit pAb	WB/IHC-P

Histone Arginine Mono-Methylation

PTM-633	Anti-Mono-Methyl-Histone H3 (Arg2) Rabbit pAb	WB/IP
PTM-700	Anti-Mono-Methyl-Histone H3 (Arg2) Rabbit mAb	WB/ICC/IF
PTM-695	Anti-Mono-Methyl-Histone H3 (Arg8) Mouse mAb	WB/IHC-P/IP/ICC/IF
PTM-610	Anti-Mono-Methyl-Histone H3 (Arg42) Rabbit pAb	WB
PTM-691	Anti-Mono-Methyl-Histone H4 (Arg3) Mouse mAb	WB/IHC-P/ChIP/ICC/IF
PTM-688	Anti-Mono-Methyl-Histone H4/H2A (Arg3) Rabbit pAb	WB
PTM-697	Anti-Mono-Methyl-Histone H4 (Arg17) Rabbit pAb	WB/ChIP
PTM-699	Anti-Mono-Methyl-Histone H4 (Arg19) Rabbit pAb	WB
PTM-694	Anti-Mono-Methyl-Histone H4 (Arg23) Rabbit pAb	WB

Histone Arginine Symmetric Di-Methylation

PTM-673	Anti-Symmetric Di-Methyl-Histone H2A (Arg3) Rabbit pAb	WB
PTM-693	Anti-Symmetric Di-Methyl-Histone H3 (Arg2) Mouse mAb	WB
PTM-652RM	Anti-Symmetric Di-Methyl-Histone H3 (Arg2) Rabbit mAb	WB/ChIP
PTM-672	Anti-Symmetric Di-Methyl-Histone H3 (Arg8) Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-679	Anti-Symmetric Di-Methyl-Histone H3 (Arg17) Rabbit pAb	WB
PTM-639	Anti-Symmetric Di-Methyl-Histone H4 (Arg3) Rabbit pAb	WB
Histone Arginine Asymmetric Di-Methylation		
PTM-660	Anti-Asymmetric Di-Methyl-Histone H2A/H4 (Arg3) Mouse mAb	WB
PTM-668	Anti-Asymmetric Di-Methyl-Histone H3 (Arg2) Mouse mAb	WB/IHC-P/ChIP/ICC/IF
PTM-671	Anti-Asymmetric Di-Methyl-Histone H3 (Arg2) Rabbit pAb	WB/ChIP
PTM-667	Anti-Asymmetric Di-Methyl-Histone H4 (Arg3) Mouse mAb	WB/IHC-P/ChIP
PTM-650RM	Anti-Asymmetric Di-Methyl-Histone H3 (Arg17) Rabbit mAb	WB/ChIP
Histone Threonine Phosphorylation		
PTM-722	Anti-Phospho-Histone H1.4 (Thr146) Rabbit pAb	WB/IP
PTM-718	Anti-Phospho-Histone H2A (Thr120) Rabbit pAb	WB/IP
PTM-730	Anti-Phospho-Histone H3 (Thr3) Mouse mAb	WB
PTM-710	Anti-Phospho-Histone H3 (Thr3) Rabbit pAb	WB/IP
PTM-711	Anti-Phospho-Histone H3 (Thr11) Rabbit pAb	WB
PTM-720	Anti-Phospho-Histone H3 (Thr32) Rabbit pAb	WB
PTM-724RM	Anti-Phospho-Histone H3 (Thr45) Rabbit mAb	WB
PTM-724	Anti-Phospho-Histone H3 (Thr45) Rabbit pAb	WB
PTM-725	Anti-Phospho-Histone H3 (Thr80) Rabbit pAb	WB
Histone Tyrosine Phosphorylation		
PTM-754	Anti-Phospho-Histone H2B (Tyr37) Mouse mAb	WB/IHC-P
Histone Serine Phosphorylation		
PTM-726	Anti-Phospho-sHistone H2A (Ser129) Rabbit pAb	WB
PTM-736	Anti-Phospho-Histone H2A.X (Ser139) Mouse mAb	WB/IP
PTM-753	Anti-Phospho-Histone H2A.X (Ser139) Mouse mAb	WB/IHC-P/ICC/IF
PTM-716	Anti-Phospho-Histone H2A.X (Ser139) Rabbit pAb	WB/IP
PTM-727	Anti-Phospho-Histone H2A.X (Ser139) Rabbit mAb	WB/IHC-P
PTM-7224	Anti-Phospho-Histone H2B (Ser14) Rabbit mAb	WB
PTM-757	Anti-Phospho-Histone H3 (Ser10) Mouse mAb	WB/IHC-P/ICC/IF
PTM-732	Anti-Phospho-Histone H3 (Ser10) Mouse mAb	WB/IHC-P
PTM-712	Anti-Phospho-Histone H3 (Ser10) Rabbit pAb	WB/IHC-P/IP
PTM-712RM	Anti-Phospho-Histone H3 (Ser10) Rabbit mAb	WB/ChIP/IP
PTM-756	Anti-Phospho-Histone H3 (Ser28) Mouse mAb	WB/ChIP
PTM-714	Anti-Phospho-Histone H3 (Ser28) Rabbit pAb	WB/IHC-P
PTM-728	Anti-Phospho-Histone H3 (Ser28) Rabbit pAb	WB

Histone Acylation Antibody Kits

Histone Lactylation & Acetylation Antibodies



• Innovative Histone Lactylation

Histone lactylation is an emerging epigenetic mark first reported in *Nature* in 2019

Tumor immunity

Inflammatory response

Hypoxic stress

• Classic Histone Acetylation First reported in 1964, histone acetylation remains one of the most extensively studied acylation modifications.

A classic epigenetic mark with enduring research value.

Epigenetic regulation

Gene expression

Metabolic regulation

Lactyl-Histone Antibody Kit Catalog


No.: PTM-7093

Product	Cat.#	Product Size
Anti-Lactyl-Histone H2B (Lys16) Rabbit mAb	PTM-1424RM	20 µL
Anti-Lactyl-Histone H3 (Lys9) Rabbit mAb	PTM-1419RM	20 µL
Anti-Lactyl-Histone H3 (Lys14) Rabbit mAb	PTM-1414RM	20 µL
Anti-Lactyl-Histone H3 (Lys18) Rabbit mAb	PTM-1406RM	20 µL
Anti-Lactyl-Histone H4 (Lys5) Rabbit mAb	PTM-1407RM	20 µL
Anti-L-Lactyl-Histone H4 (Lys8) Rabbit mAb	PTM-1415RM	20 µL
Anti-Lactyl-Histone H4 (Lys12) Rabbit mAb	PTM-1411RM	20 µL
Anti-Lactyl-Histone H4 (Lys16) Rabbit mAb	PTM-1417RM	20 µL
Anti-Histone H3 Rabbit mAb (NT)	PTM-1001RM	100 µL
Anti-Histone H4 Rabbit mAb (CT)	PTM-1015RM	100 µL

Product	Cat.#	Product Size
Anti-Acetyl-Histone H3 (Lys4) Rabbit mAb	PTM-188	20 µL
Anti-Acetyl-Histone H3 (Lys9) Rabbit pAb	PTM-112	20 µL
Anti-Acetyl-Histone H3(Lys14) Rabbit mAb	PTM-113RM	20 µL
Anti-Acetyl-Histone H3 (Lys18) Mouse mAb	PTM-158	20 µL
Anti-Acetyl-Histone H3 (Lys27) Mouse mAb	PTM-160	20 µL
Anti-Acetyl-Histone H3 (Lys36) Rabbit pAb	PTM-117	20 µL
Anti-Acetyl-Histone H3 (Lys56) Mouse mAb	PTM-162	20 µL
Anti-Acetyl-Histone H3 (Lys122) Rabbit mAb	PTM-184RM	20 µL

H3K18 Modification Antibody Kit

H3



Histone H3K18 harbors diverse modification types that can influence chromatin structure and accessibility, thereby regulating gene expression and transcriptional activity. Altered H3K18 modification patterns have been implicated in cancer, neurodegenerative disorders, and other disease-related processes.

Broad H3K18 PTM Coverage


Smart Bundle, Better Value

Histone H3(Lys18) modification Antibody Kit

Cat.No.: PTM-7477

Product	Cat. No	Product Size
Anti-L-Lactyl-Histone H3 (Lys18) Rabbit mAb-ChIP Grade	PTM-1427RM	20 µl
Anti-Acetyl-Histone H3 (Lys18) Mouse mAb	PTM-158	20 µl
Anti-Butyryl-Histone H3 (Lys18) Mouse mAb	PTM-331	20 µl
Anti-Crotonyl-Histone H3 (Lys18) Mouse mAb	PTM-540	20 µl
Anti-Propionyl-Histone H3 (Lys18) Mouse mAb	PTM-214	20 µl
Anti-2-Hydroxyisobutyryl-Histone H3 (Lys18) Mouse mAb	PTM-882	20 µl
Anti-β-Hydroxybutyryl-Histone H3 (Lys18) Mouse mAb	PTM-1292	20 µl
Anti-Mono-Methyl-Histone H3 (Lys18) Rabbit mAb	PTM-618RM	20 µl
Anti-Tri-Methyl-Histone H3 (Lys18) Rabbit pAb	PTM-632	20 µl
Anti-Histone H3 Rabbit mAb	PTM-1002RM	100 µl

High-Demand Histone PTM Antibody Kit



Popular Histone PTM Antibodies

Includes high-demand histone modification types and site-specific antibodies, supporting diverse epigenetics research needs.

ChIP-Validated Performance

All included PTM antibodies are suitable for ChIP applications, providing reliable support for chromatin and epigenetic regulation studies.

Histone H3(Ac/Me/La) modification Antibody Kit

Cat.No.: PTM-7480

Product	Cat. No	Product Size
Anti-L-Lactyl-Histone H3 (Lys14) Rabbit mAb	PTM-1414RM	20 µL
Anti-L-Lactyl-Histone H3 (Lys18) Rabbit mAb-ChIP Grade	PTM-1427RM	20 µL
Anti-L-Lactyl-Histone H3 (Lys27) Rabbit pAb	PTM-1428	20 µL
Anti-L-Lactyl-Histone H3 (Lys9) Rabbit mAb	PTM-1419RM	20 µL
Anti-Acetyl-Histone H3 (Lys14) Rabbit mAb	PTM-113RM	20 µL
Anti-Acetyl-Histone H3 (Lys27) Rabbit mAb	PTM-116RM	20 µL
Anti-Acetyl-Histone H3 (Lys4) Rabbit pAb	PTM-124	20 µL
Anti-Acetyl-Histone H3 (Lys9) Rabbit mAb	PTM-112RM	20 µL
Anti-Tri-Methyl-Histone H3 (Lys27) Rabbit mAb	PTM-647RM	20 µL
Anti-Tri-Methyl-Histone H3 (Lys4) Rabbit pAb	PTM-613	20 µL
Anti-Histone H3 Rabbit mAb (CT)	PTM-1002RM	100 µL

Epigenetics Antibodies

Cat.#	Product	Recommended Applications
Histone Antibodies		
PTM-6054	Anti-Histone H1 Rabbit mAb (CT)	WB/IHC-P/ICC/IF/FC
PTM-20151	Anti-Histone H1.2 Rabbit mAb	WB
PTM-1008	Anti-Histone H2A Rabbit pAb	WB/IHC-P
PTM-1005	Anti-Histone H2A.X Mouse mAb (CT)	WB
PTM-7034	Anti-Histone H2A.Z Rabbit mAb	WB/IHC-P/ICC/IF
PTM-1007	Anti-Histone H2B Rabbit pAb (CT)	WB/IHC-P
PTM-1006	Anti-Histone H2B Mouse mAb (NT)	WB/IHC-P
PTM-1007RM	Anti-Histone H2B Rabbit mAb	WB/IHC-P
PTM-1014	Anti-Histone H3.3 (mutated K36M) Rabbit pAb	WB
PTM-7294	Anti-Histone H3.3 Rabbit mAb	WB/IHC-P
PTM-6600	Anti-Histone H3 Rabbit mAb	WB/IHC-P/ChIP
PTM-6613	Anti-Histone H3 Rabbit mAb	WB/IHC-P
PTM-1002RM	Anti-Histone H3 Rabbit mAb (CT)	WB/IHC-P/ICC/IF
PTM-1001RM	Anti-Histone H3 Rabbit mAb (NT)	WB/IHC-P/ICC/IF/FC/ChIP
PTM-1002	Anti-Histone H3 Rabbit pAb (CT)	WB/IHC-P
PTM-6621	Anti-Histone H3 (HRP) Rabbit mAb (NT)	WB/IHC-P
PTM-1009	Anti-Histone H4 Mouse mAb (NT)	WB/IHC-P
PTM-1015RM	Anti-Histone H4 Rabbit mAb (CT)	WB/IHC-P
Histone Acyltransferases		
PTM-5711	Anti-PCAF Rabbit mAb	WB/IP
PTM-5281	Anti-CBP Rabbit mAb	WB/ICC/IF
PTM-6980	Anti-CBP Rabbit pAb	WB
PTM-5115	Anti-GCN5 Mouse mAb	WB
PTM-5195	Anti-HAT1 Mouse mAb	WB/IHC-P
PTM-7389	Anti-MYST1 Mouse mAb	WB/ICC/IF
PTM-20296	Anti-ACSS2 Rabbit pAb	WB
Histone Deacylases		
PTM-5377	Anti-HDAC1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6350	Anti-HDAC1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6409	Anti-HDAC1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6006	Anti-HDAC10 Mouse mAb	WB/ICC/IF
PTM-7219	Anti-HDAC2 Mouse mAb	WB
PTM-5183	Anti-HDAC3 Mouse mAb	WB
PTM-5551	Anti-HDAC4 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-7412	Anti-HDAC6 Mouse mAb	WB
PTM-6323	Anti-HDAC6 Rabbit mAb	WB/IHC-P
PTM-6669	Anti-HDAC6 Rabbit mAb	WB/IHC-P
PTM-5172	Anti-HDAC8 Mouse mAb	WB
PTM-755	Anti-Phospho-HDAC1 (Ser406) Mouse mAb	WB
PTM-6945	Anti-Phospho-HDAC4(Ser632)/HDAC5(Ser661)/HDAC7(Ser486) Rabbit pAb	WB
PTM-7418	Anti-Phospho-HDAC5 (Ser498) Rabbit mAb	WB
PTM-6601	Anti-Phospho-SIRT1(Ser27) Rabbit mAb	WB/ICC/IF
PTM-5021	Anti-SIRT1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6043	Anti-SIRT2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6318	Anti-SIRT2 Rabbit mAb	WB/IHC-P
PTM-6462	Anti-SIRT3 Rabbit mAb	WB
PTM-7127	Anti-SIRT5 Mouse mAb	WB
PTM-7223	Anti-SIRT5 Rabbit mAb	WB/IHC-P
PTM-7128	Anti-SIRT5 Rabbit mAb	WB
PTM-6658	Anti-SIRT6 Rabbit mAb	WB
PTM-7231	Anti-SIRT7 Rabbit mAb	WB

Histone Methyltransferases

PTM-6205	Anti-EHMT1 Rabbit mAb	WB
PTM-5327	Anti-EHMT2 Mouse mAb	WB/ICC/IF
PTM-7414	Anti-EZH1 Rabbit mAb	WB
PTM-6310	Anti-EZH2 Mouse mAb	WB/IHC-P
PTM-6272	Anti-EZH2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5227	Anti-EZH2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7292	Anti-NSD3 Rabbit mAb	WB
PTM-5154	Anti-PRMT1 Mouse mAb	WB
PTM-5889	Anti-PRMT5 Rabbit mAb	WB/IHC-P/IP
PTM-7383	Anti-PRMT7 Rabbit mAb	WB
PTM-7305	Anti-SET7 Rabbit mAb	WB
PTM-7278	Anti-SET7 Rabbit mAb	WB
PTM-6112	Anti-SETD1A Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5365	Anti-SETDB1 Mouse mAb	WB
PTM-5333	Anti-SETDB1 Mouse mAb	WB
PTM-6668	Anti-SETDB2 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-6552	Anti-SUV39H2 Rabbit mAb	WB/IHC-P
PTM-5987	Anti-SUV39H2 Rabbit mAb	WB
PTM-5948	Anti-SUV39H1 Rabbit mAb	WB
Histone Demethylases		
PTM-7321	Anti-JMJD2A Rabbit mAb	WB
PTM-7306	Anti-JMJD2D Rabbit mAb	WB
PTM-7290	Anti-JMJD6 Mouse mAb	WB/IHC-P
PTM-7462	Anti-KDM1 Rabbit mAb	WB
PTM-6203	Anti-KDM2A Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5215	Anti-KDM5C/JARID1C/SmcX Mouse mAb	WB/ICC/IF
PTM-6654	Anti-KDM5C/JARID1C/SmcX Mouse mAb	WB
PTM-5938	Anti-LSD1 Mouse mAb	WB/IHC-P
PTM-5960	Anti-LSD1 Mouse mAb	WB/IHC-P
PTM-6278	Anti-LSD2/AOF1 Rabbit mAb	WB/ICC/IF/FC
PTM-7039	Anti-PHF2 Rabbit mAb	WB
PTM-6496	Anti-PHF8 Rabbit mAb	WB
DNA Methyltransferases		
PTM-6317	Anti-DNMT1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7448	Anti-DNMT3B Rabbit mAb	WB/ChIP
PTM-6074	Anti-DNMT3L Mouse mAb	WB
RNA Methyltransferases		
PTM-6001	Anti-NSUN2 Mouse mAb	WB/IHC-P
PTM-7028	Anti-METTTL3 Mouse mAb	WB/IHC-P
PTM-7032	Anti-METTTL3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6706	Anti-METTTL7A Mouse mAb	WB/IHC-P
PTM-7386	Anti-TRMT112 Mouse mAb	WB/ICC/IF
RNA Demethylases		
PTM-6922	Anti-FTO Rabbit mAb	WB/IHC-P
Chromatin Remodeling Factors		
PTM-6308	Anti-SMARCD1 Rabbit mAb	WB
PTM-6893	Anti-SMARCD2 Rabbit mAb	WB
PTM-5557	Anti-BAF60C Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5419	Anti-BRD8 Mouse mAb	WB/IHC-P
PTM-5316	Anti-BRD2 (phospho S37) Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-7331	Anti-BRD2 Rabbit mAb	WB/IHC-P/ChIP
PTM-7344	Anti-BAF53B Rabbit mAb	WB
PTM-5826	Anti-BRG1 Mouse mAb	WB
PTM-6990	Anti-CHD1 Rabbit mAb	WB
PTM-6377	Anti-CHD3 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6982	Anti-CHD4 Rabbit mAb	WB/IHC-P
PTM-6762	Anti-CHD4 Rabbit pAb	WB/IHC-P
PTM-6213	Anti-HP1 alpha Rabbit mAb	WB
PTM-7322	Anti-HP1 beta Rabbit mAb	WB
PTM-6383	Anti-HP1 gamma Mouse mAb	WB/IHC-P/ICC/IF
PTM-7111	Anti-PBRM1/BAF180 Rabbit mAb	WB
PTM-6889	Anti-SMARCA2/BRM Rabbit mAb	WB/IHC-P
PTM-5144	Anti-SMARCA5 Mouse mAb	WB/ICC/IF
PTM-7336	Anti-SMARCB1/BAF47 Rabbit mAb	WB
PTM-7011	Anti-SMARCC1/BAF155 Rabbit pAb	WB/IHC-P
PTM-7472	Anti-SMARCC2/BAF170 Rabbit mAb	WB/IHC-P/FC/IP
PTM-7036	Anti-SMARCC2/BAF170 Rabbit pAb	WB/IHC-P

Transcription Factors

PTM-5035	Anti-Androgen Receptor Rabbit mAb	WB/IHC-P
PTM-6606	Anti-ATF1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7396	Anti-ATF3 Rabbit mAb	WB/IHC-P
PTM-5041	Anti-BCL6 Rabbit mAb	WB/IHC-P
PTM-6653	Anti-c-Fos Mouse mAb	WB
PTM-7474	Anti-c-Jun Rabbit mAb	WB/IP
PTM-7317	Anti-c-Myc Rabbit mAb	WB/ChIP/IP
PTM-7314	Anti-C/EBP beta Mouse mAb	WB/IHC-P
PTM-7216	Anti-C/EBP delta Mouse mAb	WB
PTM-7478	Anti-CDX1 Rabbit mAb	WB/IHC-P/IP
PTM-7113	Anti-COUP-TFII Rabbit mAb	WB/IHC-P/ChIP
PTM-6548	Anti-CREB3L1 Mouse mAb	WB
PTM-6772	Anti-CTCF Mouse mAb	WB/IHC-P
PTM-7455	Anti-CTCF L Rabbit mAb	WB/IHC-P/FC
PTM-6156	Anti-Ctip2 Rabbit mAb	WB/IHC-P/FC
PTM-7323	Anti-Dlx4/BP1 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-7415	Anti-Dlx5 Mouse mAb	WB
PTM-7106	Anti-DMRT1 Rabbit mAb	WB
PTM-6630	Anti-E2F1 Mouse mAb	WB
PTM-7264	Anti-E2F6 Rabbit mAb	WB
PTM-7199	Anti-EGR2 Rabbit mAb	WB
PTM-7374	Anti-EN1/Engrailed 1 Rabbit mAb	WB
PTM-6265	Anti-ERG Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6967	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6838	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6792	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6816	Anti-Estrogen Related Receptor alpha Rabbit mAb	WB
PTM-5395	Anti-ETV5 Mouse mAb	WB
PTM-7295	Anti-FOXA2 Rabbit mAb	WB/ChIP
PTM-7002	Anti-FOXC1 Rabbit mAb	WB
PTM-5324	Anti-FOXA1 Rabbit mAb	IHC-P/ICC/IF
PTM-6083	Anti-FOXO3A Mouse mAb	WB/IHC-P/ICC/IF
PTM-6857	Anti-FOXP2 Mouse mAb	WB
PTM-5742	Anti-FOXL2 Rabbit mAb	WB/IHC-P
PTM-5135	Anti-FOXP1 Rabbit mAb	WB/IHC-P
PTM-5016	Anti-FOXP3 Rabbit mAb	IHC-P
PTM-7347	Anti-FUBP1 Rabbit mAb	WB
PTM-7444	Anti-GATA1 Rabbit mAb	WB/IP/ChIP
PTM-6994	Anti-GATA1 Rabbit mAb	WB/IHC-P
PTM-6465	Anti-GATA2 Mouse mAb	WB
PTM-6077	Anti-GATA3 Mouse mAb	WB/ICC/IF
PTM-5197	Anti-GATA3 Mouse mAb	WB/IHC-P
PTM-6429	Anti-GATA3 Rabbit mAb	WB/IHC-P
PTM-7446	Anti-Glucocorticoid receptor Rabbit mAb	WB/ICC/IF
PTM-6720	Anti-HES5 Rabbit mAb	WB
PTM-6320	Anti-HNF-1B Rabbit mAb	WB
PTM-6937	Anti-HOXB13 Rabbit mAb	WB
PTM-6814	Anti-HOXB4 Rabbit mAb	WB
PTM-7301	Anti-HSF1 Mouse mAb	WB
PTM-6237	Anti-IRF6 Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6442	Anti-IRF8 Mouse mAb	WB/ICC/IF
PTM-7463	Anti-IRF9 Rabbit mAb	WB
PTM-5672	Anti-JunB Rabbit mAb	WB
PTM-5238	Anti-KLF4 Rabbit mAb	IHC-P/ICC/IF
PTM-5057	Anti-LMO2 Rabbit mAb	IHC-P
PTM-6651	Anti-LMO4 Rabbit mAb	WB
PTM-7175	Anti-MASH1 Rabbit mAb	WB
PTM-7343	Anti-MAX Rabbit mAb	WB/ChIP
PTM-6877	Anti-MAZ Rabbit pAb	WB/IHC-P
PTM-7341	Anti-MEF2C Rabbit mAb	WB/IHC-P
PTM-7356	Anti-N-Myc Rabbit mAb	WB
PTM-6010	Anti-Nanog Mouse mAb	WB
PTM-7312	Anti-NeuroD1 Rabbit mAb	WB
PTM-5706	Anti-NF-κB1 p105/p50 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6593	Anti-NF-κB2 p100/p52 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7042	Anti-NFIB Rabbit mAb	WB/IHC-P
PTM-7324	Anti-Notch1 Rabbit mAb	WB/IHC-P
PTM-7225	Anti-Notch4 Rabbit mAb	WB
PTM-6529	Anti-NR2F6/EAR-2 Mouse mAb	WB
PTM-5593	Anti-OCT4 Mouse mAb	WB
PTM-6466	Anti-p53 Mouse mAb	WB
PTM-6517	Anti-p53 Rabbit mAb	WB
PTM-6319	Anti-p53 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7442	Anti-PAX5 Rabbit mAb	WB/ICC/IF/FC
PTM-5338	Anti-PAX8 Rabbit mAb	IHC-P
PTM-6248	Anti-PAX9 Rabbit mAb	WB/IHC-P
PTM-7445	Anti-PBX1 Rabbit mAb	WB/IHC-P/IP/ChIP
PTM-6907	Anti-PDX1 Rabbit pAb	WB/IHC-P
PTM-7209	Anti-Phospho-c-Myb (Ser11) Rabbit mAb	WB
PTM-7425	Anti-Phospho-p53 (Thr377) Rabbit mAb	WB
PTM-7419	Anti-Phospho-STAT1 (Ser727) Rabbit mAb	WB/IHC-P
PTM-7024	Anti-Phospho-Stat3 (Ser727) Rabbit mAb	WB
PTM-7293	Anti-Pim-1 Rabbit mAb	WB
PTM-7010	Anti-PLAG1 Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6523	Anti-PRAME Rabbit mAb	WB
PTM-5204	Anti-PRDM1 Mouse mAb	WB/ICC/IF
PTM-6337	Anti-PRDM1/Blimp1 Rabbit mAb	WB
PTM-7338	Anti-PRDM4/PFM1 Rabbit mAb	WB/IHC-P
PTM-5085	Anti-Progesterone Receptor Rabbit mAb	IHC-P
PTM-7297	Anti-PROX1 Mouse mAb	WB/IHC-P
PTM-5270	Anti-PU.1 Rabbit mAb	WB/IHC-P
PTM-7303	Anti-RUNX1/AML1 Rabbit mAb	WB
PTM-6840	Anti-RUNX1/AML1 Rabbit pAb	WB
PTM-5253	Anti-RUNX2 Mouse mAb	WB
PTM-6617	Anti-RUNX2 Rabbit mAb	WB
PTM-6547	Anti-RUNX3 Rabbit mAb	WB/IP
PTM-7025	Anti-RXRA Rabbit mAb	WB/ICC/IF
PTM-6638	Anti-SALL2 Mouse mAb	WB
PTM-5293	Anti-SATB1 Rabbit mAb	IHC-P
PTM-5295	Anti-SATB2 Rabbit mAb	IHC-P
PTM-5799	Anti-SF1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6929	Anti-Smad2 (phospho S255) Rabbit mAb	WB
PTM-6688	Anti-Snail Mouse mAb	WB
PTM-7065	Anti-SOX17 Rabbit mAb	WB/IHC-P
PTM-6390	Anti-SOX2 Rabbit mAb	WB/IHC-P
PTM-5072	Anti-SOX2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5076	Anti-SOX10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7053	Anti-SP1 Rabbit mAb	WB/IHC-P
PTM-6599	Anti-STAT1 Mouse mAb	WB
PTM-5754	Anti-STAT1 Rabbit mAb	WB/IHC-P/IP
PTM-5090	Anti-STAT3 Mouse mAb	WB/IHC-P
PTM-5318	Anti-STAT5b Mouse mAb	WB/ICC/IF
PTM-5828	Anti-STAT5b Rabbit mAb	WB
PTM-6057	Anti-STAT6 Rabbit mAb	WB/ICC/IF/FC
PTM-6805	Anti-STAT6 Rabbit pAb	WB
PTM-5244	Anti-T-bet Mouse mAb	WB
PTM-7461	Anti-TATA binding protein TBP antibody	WB/IP
PTM-6042	Anti-TBP Mouse mAb	WB/IHC-P/FC

Cat.#	Product	Recommended Applications
PTM-7422	Anti-TBX1 Rabbit mAb	WB
PTM-6288	Anti-TCF-3/TCF7L1 Mouse mAb	WB
PTM-6518	Anti-TCF3/E2A Mouse mAb	WB/IHC-P
PTM-6451	Anti-TCF3/E2A Mouse mAb	WB
PTM-7280	Anti-TCF4 Rabbit mAb	WB
PTM-6998	Anti-USF1 Rabbit mAb	WB
PTM-7460	Anti-Vitamin D Receptor Rabbit mAb	WB/IP
PTM-7298	Anti-WT1 Rabbit mAb	WB/IHC-P
PTM-7035	Anti-WT1 Rabbit pAb	WB
Others		
PTM-6766	Anti-Aconitase 1/ACO1 Rabbit mAb	WB
PTM-7349	Anti-ADA Rabbit mAb	WB
PTM-7283	Anti-ADA2 Rabbit mAb	WB/IHC-P
PTM-7253	Anti-AIF Mouse mAb	WB/IHC-P
PTM-6735	Anti-AKT (phospho T308) Rabbit mAb	WB
PTM-7179	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-7271	Anti-APE1 Rabbit mAb	WB/IHC-P
PTM-6430	Anti-APE1 Rabbit mAb	WB/IHC-P/IP
PTM-5149	Anti-ARD1A Mouse mAb	WB
PTM-5560	Anti-Arp5/ACTR5 Mouse mAb	WB
PTM-5206	Anti-ASH2L Mouse mAb	WB
PTM-7207	Anti-Ataxin-1 Rabbit mAb	WB
PTM-5946	Anti-ATM Rabbit mAb	WB
PTM-5709	Anti-BMI1 Mouse mAb	WB
PTM-7387	Anti-BRD1 Mouse mAb	WB/ICC/IF
PTM-7331	Anti-BRD2 Rabbit mAb	WB/IHC-P/ChIP
PTM-6487	Anti-c-Cbl Mouse mAb	WB/ICC/IF
PTM-6562	Anti-CBX4 Rabbit mAb	WB/IHC-P
PTM-6627	Anti-CD82 Mouse mAb	FC
PTM-7410	Anti-Cdc25B Rabbit mAb	WB/IHC-P
PTM-6585	Anti-Cdc25C Mouse mAb	WB
PTM-6850	Anti-Cdc25C Rabbit mAb	WB
PTM-5476	Anti-Cdc37 Mouse mAb	WB
PTM-6705	Anti-Cdc42 Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6521	Anti-CDK1/Cdc2 Mouse mAb	WB/ICC/IF
PTM-6396	Anti-CDK2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7154	Anti-CDKN2A/p16INK4a Rabbit mAb	WB/IHC-P
PTM-6704	Anti-CDKN2A/p16INK4a Rabbit mAb	WB
PTM-5152	Anti-CECR2 Mouse mAb	WB
PTM-6089	Anti-CENP-A Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6536	Anti-CEP63 Mouse mAb	WB/ICC/IF
PTM-6276	Anti-Chk2 Mouse mAb	WB/IHC-P
PTM-5953	Anti-Chk2 Mouse mAb	WB
PTM-7275	Anti-CIRBP Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7073	Anti-Citrulline-NPM1 (Arg196) Rabbit pAb	WB
PTM-7451	Anti-Cleaved PARP1 p25 Rabbit mAb	WB/IHC-P/ICC/IF/C/IP/CUT&Tag
PTM-7466	Anti-Cleaved PARP1 p85 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6182	Anti-CNBP Mouse mAb	WB/IHC-P
PTM-5174	Anti-CNOT1 Mouse mAb	WB
PTM-6098	Anti-CNOT1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6748	Anti-CNOT6 Rabbit pAb	WB/IHC-P
PTM-5494	Anti-COMMD1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7385	Anti-CRP2 Rabbit mAb	WB
PTM-5193	Anti-CRTC2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7473	Anti-CSDE1 Rabbit mAb	WB/IHC-P/FC/IP
PTM-6987	Anti-CSN2 Mouse mAb	WB
PTM-7430	Anti-CTNNBIP1 Rabbit mAb	WB
PTM-6659	Anti-Cyclin B1 Rabbit mAb	WB/ICC/IF/IP
PTM-6664	Anti-Cyclin D2 Rabbit mAb	WB/IHC-P
PTM-7181	Anti-Cytoplasmic Dynein Intermediate chain Mouse mAb	WB
PTM-5413	Anti-DBR1 Mouse mAb	WB/IHC-P
PTM-7183	Anti-DcR1 Rabbit mAb	WB/ICC/IF
PTM-6963	Anti-DDB2 Rabbit mAb	WB/IHC-P
PTM-6795	Anti-DDX3 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6370	Anti-DDX4 Rabbit mAb	IHC-P
PTM-6904	Anti-DDX4 Rabbit pAb	WB/IHC-P
PTM-7453	Anti-DGCR8 Rabbit mAb	WB/IHC-P/IP
PTM-6913	Anti-DHX29 Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-6636	Anti-DHX9 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5335	Anti-Dicer Mouse mAb	WB/ICC/IF
PTM-6757	Anti-DJ-1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7262	Anti-DNA Polymerase beta Rabbit mAb	WB/IHC-P
PTM-7352	Anti-EED Rabbit mAb	WB
PTM-6530	Anti-eEF1A1 Rabbit mAb	WB/ICC/IF
PTM-7114	Anti-eEF1A2 Rabbit mAb	WB/IHC-P
PTM-6843	Anti-eEF1A2 Rabbit pAb	WB/IHC-P
PTM-6236	Anti-eEF1B2 Mouse mAb	WB/ICC/IF
PTM-5719	Anti-eEF1G Rabbit mAb	WB/IHC-P
PTM-7267	Anti-eEF2 Mouse mAb	WB
PTM-7360	Anti-EID1 Rabbit mAb	WB/IHC-P
PTM-5609	Anti-eIF3B Rabbit mAb	WB/IHC-P
PTM-6079	Anti-eIF3E Mouse mAb	WB
PTM-6549	Anti-eIF3F Rat mAb	WB
PTM-7101	Anti-eIF3G Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7468	Anti-eIF4A2 Rabbit mAb	WB/IHC-P/IP
PTM-6802	Anti-eIF4E Rabbit mAb	WB/IHC-P
PTM-7033	Anti-eIF4G1 Rabbit mAb	WB/IHC-P
PTM-5961	Anti-eIF5A Mouse mAb	WB/ICC/IF
PTM-7355	Anti-ELL Mouse mAb	WB
PTM-5980	Anti-ELP3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6385	Anti-ELP4 Rabbit mAb	WB/ICC/IF/IP
PTM-7348	Anti-ENO1 Rabbit mAb	WB
PTM-6780	Anti-ENO1 Rabbit pAb	WB/IHC-P
PTM-7052	Anti-EWSR1/EWS Rabbit pAb	WB/IHC-P
PTM-7388	Anti-FARSA Mouse mAb	WB/IHC-P
PTM-6930	Anti-FARSB Mouse mAb	WB/IHC-P
PTM-6960	Anti-FEN1 Rabbit mAb	WB/IHC-P
PTM-6941	Anti-FKBP52 Rabbit mAb	WB
PTM-6519	Anti-FKBP52 Rabbit mAb	WB
PTM-6737	Anti-GART Mouse mAb	WB/ICC/IF
PTM-7450	Anti-GTF2E2 Rabbit mAb	WB/ICC/IF/FC
PTM-6556	Anti-GTF3C4 Mouse mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7212	Anti-HACE1 Rabbit mAb	WB
PTM-7193	Anti-Hamartin Rabbit mAb	WB
PTM-6380	Anti-HEC1 Mouse mAb	WB
PTM-6303	Anti-HMGA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6123	Anti-HMGB1 Rabbit mAb	WB/IHC-P/FC/IP
PTM-7137	Anti-HMGB1 Rat mAb	WB
PTM-6611	Anti-HMGB2 Rabbit mAb	WB/IHC-P/IC/IF/IP
PTM-7092	Anti-hnRNP A1 Rabbit pAb	WB
PTM-6555	Anti-hnRNP A2/B1 Mouse mAb	WB/ICC/IF
PTM-7255	Anti-hnRNP E2/PCBP2 Rabbit mAb	WB
PTM-6026	Anti-hnRNP F Mouse mAb	WB/ICC/IF
PTM-6426	Anti-hnRNP H Mouse mAb	WB/IHC-P/ICC/IF
PTM-7334	Anti-hnRNP K Mouse mAb	WB
PTM-7320	Anti-hnRNP K Mouse mAb	WB/IHC-P
PTM-5382	Anti-HRPT2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5432	Anti-HuR Mouse mAb	WB/IHC-P/ICC/IF
PTM-6822	Anti-IBP160 Rabbit pAb	WB
PTM-7237	Anti-IFI16 Mouse mAb	WB
PTM-6682	Anti-IGF2BP3 Mouse mAb	WB
PTM-6546	Anti-IKK gamma Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6490	Anti-ILF2/NF45 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6366	Anti-ILF3 Rabbit mAb	WB/IHC-P/IC/IF/FC
PTM-6460	Anti-ING1 Rabbit mAb	WB
PTM-7300	Anti-ING3 Rabbit mAb	WB
PTM-6102	Anti-ING5 Rabbit mAb	WB/IHC-P
PTM-6927	Anti-ISG20 Rabbit mAb	WB/IHC-P
PTM-7135	Anti-JAK2 Rabbit mAb	WB
PTM-6456	Anti-KAP1 Mouse mAb	WB/IHC-P
PTM-7304	Anti-KEAP1 Rabbit mAb	WB
PTM-5298	Anti-Keratin 7 Rabbit mAb	WB/IHC-P/IC/IF/IP
PTM-7359	Anti-LIF Rat mAb	WB/IHC-P
PTM-7340	Anti-LIMK1 Rabbit mAb	WB
PTM-7420	Anti-LIN28A Rabbit mAb	WB/ICC/IF
PTM-5867	Anti-LONP1 Mouse mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6121	Anti-LysRS Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7102	Anti-MAVS Rabbit pAb	WB
PTM-7411	Anti-MBD2 Mouse mAb	WB
PTM-6892	Anti-MBD2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7447	Anti-MBD3 Rabbit mAb	WB/ChIP
PTM-6592	Anti-MBNL2 Mouse mAb	WB/IHC-P
PTM-5074	Anti-MCM2 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5903	Anti-MCM4 Mouse mAb	WB
PTM-5839	Anti-MCM5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6263	Anti-MCM7 Mouse mAb	WB
PTM-7469	Anti-MCM7 Rabbit mAb	WB/IP/ICC/IF
PTM-6363	Anti-MCM7 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7310	Anti-MDM2 Mouse mAb	WB
PTM-7302	Anti-MED1 Rabbit mAb	WB
PTM-7329	Anti-MED12 Mouse mAb	WB
PTM-7254	Anti-MED26 Mouse mAb	WB
PTM-6459	Anti-MED30 Mouse mAb	WB
PTM-6284	Anti-MED4 Rabbit mAb	WB/IHC-P
PTM-7222	Anti-MED7 Rabbit mAb	WB/IHC-P
PTM-7316	Anti-MEP50 Mouse mAb	WB/IHC-P
PTM-7177	Anti-Met Mouse mAb	WB
PTM-3102	Anti-Methylglyoxal Rabbit pAb	WB
PTM-6215	Anti-MIS12 Mouse mAb	WB
PTM-6512	Anti-MLH1 Rabbit mAb	WB
PTM-7404	Anti-MORF4L2 Rabbit mAb	WB/IHC-P
PTM-6851	Anti-Mre11 Mouse mAb	WB
PTM-5342	Anti-MSH2 Rabbit mAb	IHC-P
PTM-5034	Anti-MSH2 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5196	Anti-MSH6 Mouse mAb	WB/IHC-P/FC
PTM-5014	Anti-MSH6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6038	Anti-MTA2 Mouse mAb	WB/IHC-P
PTM-7351	Anti-MTF2 Rabbit mAb	WB
PTM-6993	Anti-mtTFA Rabbit mAb	WB/IHC-P
PTM-7345	Anti-NAP1L1 Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6035	Anti-NARS Mouse mAb	WB/ICC/IF
PTM-7299	Anti-NCK1 Rabbit mAb	WB
PTM-7353	Anti-NCOA4 Mouse mAb	WB
PTM-7328	Anti-NCoR1 Rabbit mAb	WB
PTM-7471	Anti-NELF-B Rabbit mAb	WB/HC-P/FC/IP
PTM-5330	Anti-NO66 Mouse mAb	WB/HC-P
PTM-5352	Anti-NPM1 Mouse mAb	WB
PTM-5310	Anti-NPM1 Mouse mAb	WB/ICC/IF
PTM-6129	Anti-NPM1 Rabbit mAb	WB/HC-P/IC/IF/FC
PTM-6903	Anti-NPM1 Rabbit pAb	WB
PTM-7171	Anti-Nucleolin Rabbit mAb	WB/HC-P
PTM-6012	Anti-NUDT21 Mouse mAb	WB/HC-P/ICC/IF
PTM-7279	Anti-NUP107 Rabbit mAb	WB
PTM-7392	Anti-OAS2 Mouse mAb	WB/ICC/IF
PTM-7125	Anti-p21 Rabbit mAb	WB
PTM-6510	Anti-p21 Rabbit mAb	WB/HC-P/IP
PTM-6471	Anti-p27 Kip1 Rabbit mAb	WB/HC-P
PTM-7260	Anti-p53R2 Rabbit mAb	WB
PTM-6338	Anti-p60 CAF1/MPP7 Rabbit mAb	WB/HC-P/ICC/IF/IP
PTM-5384	Anti-PABPN1 Mouse mAb	WB/HC-P/ICC/IF
PTM-5171	Anti-PARP1 Mouse mAb	WB
PTM-7475	Anti-PARP2 Rabbit mAb	WB
PTM-7431	Anti-PARP4 Rabbit mAb	WB
PTM-7337	Anti-PCGF5 Rabbit mAb	WB
PTM-7363	Anti-PCGF6 Rabbit mAb	WB
PTM-5187	Anti-PCNA Mouse mAb	WB/HC-P
PTM-7257	Anti-PCQAP/MED15 Rabbit mAb	WB
PTM-6410	Anti-PDCD2 Rabbit mAb	WB
PTM-6482	Anti-PDCD4 Mouse mAb	WB/HC-P/ICC/IF
PTM-7066	Anti-PEG10 Rabbit mAb	WB
PTM-7174	Anti-Perilipin-1 Rabbit mAb	WB
PTM-5914	Anti-PHB2 Rabbit mAb	WB/ICC/IF
PTM-5941	Anti-PHF20 Rabbit mAb	WB/ICC/IF
PTM-6885	Anti-Phospho S6 Ribosomal Protein (Ser235/236) Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6515	Anti-Phospho-4E-BP1 (Thr37/46) Rabbit mAb	WB
PTM-6649	Anti-Phospho-AKT1 (Ser473) Rabbit mAb	WB
PTM-6887	Anti-Phospho-CAD (Ser1859) Rabbit pAb	WB
PTM-7417	Anti-Phospho-Cdc37 (Ser13) Rabbit mAb	WB
PTM-6786	Anti-Phospho-Cdc37 (Ser13) Rabbit pAb	WB
PTM-7117	Anti-Phospho-Chk2 (Ser19) Rabbit pAb	WB/ICC/IF
PTM-7136	Anti-Phospho-Chk2 (Thr68) Rabbit mAb	WB
PTM-7123	Anti-Phospho-HSF1 (Ser326) Rabbit mAb	WB
PTM-6957	Anti-Phospho-HSF1 (Ser326) Rabbit pAb	WB
PTM-7236	Anti-Phospho-IkB alpha (Tyr42) Rabbit mAb	WB
PTM-7384	Anti-Phospho-MCM2 (Ser139) Rabbit mAb	WB/IHC-P/IP
PTM-7350	Anti-Phospho-NEDD4L (Ser448) Rabbit mAb	WB
PTM-6954	Anti-Phospho-NPM1 (Thr199) Rabbit pAb	WB
PTM-7381	Anti-Phospho-Rb(Thr826) Rabbit mAb	WB
PTM-7309	Anti-Phospho-RNA Polymerase II CTD Repeat YSPTSPS (Ser2) Rabbit mAb	WB/IHC-P/ChIP
PTM-7191	Anti-Phospho-RNA Polymerase II CTD Repeat YSPTSPS (Ser5) Rabbit mAb	WB
PTM-7122	Anti-Phospho-RNA polymerase II CTD repeat YSPTSPS (Ser5) Rabbit pAb	WB
PTM-5321	Anti-Phospho-S6 Ribosomal Protein (Ser235/Ser236) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5299	Anti-Phospho-S6 Ribosomal Protein (Ser240/244) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6714	Anti-Phospho-TBK1/NAK (Ser172) Rabbit mAb	WB/IP
PTM-7170	Anti-PKM2 Rabbit mAb	WB/IHC-P
PTM-7416	Anti-PLK1 Rabbit mAb	WB
PTM-7361	Anti-PNKP Rabbit mAb	WB/IHC-P
PTM-7409	Anti-PRDX5 Mouse mAb	WB/IHC-P
PTM-6865	Anti-pro Caspase-1 Rabbit mAb	WB/IHC-P
PTM-6616	Anti-PURA Mouse mAb	WB
PTM-6687	Anti-QKI Rabbit mAb	WB/ICC/IF
PTM-7443	Anti-Rad21 Rabbit mAb	WB/IHC-P/ICC/IF/FC/ChIP
PTM-6419	Anti-Rad51 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5911	Anti-RAG2 Rabbit mAb	WB
PTM-6914	Anti-RALY Rabbit pAb	WB
PTM-6773	Anti-Ras Mouse mAb	WB/IHC-P
PTM-5979	Anti-RBBP4 Rabbit mAb	WB/IHC-P
PTM-7459	Anti-RBBP7 Rabbit mAb	WB/IP

Cat.#	Product	Recommended Applications
PTM-7470	Anti-RBBP7/RBBP4 Rabbit mAb	WB/IHC-P/FC/ChIP/IP
PTM-6312	Anti-RBM15B Mouse mAb	WB
PTM-5421	Anti-RBM4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7452	Anti-RBMS2 Rabbit mAb	WB
PTM-6474	Anti-RBX1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6605	Anti-RING1B Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7009	Anti-RING2/RING1B/RNF2 Rabbit pAb	WB
PTM-5259	Anti-RMP/C19orf2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7476	Anti-RPA32/RPA2 Rabbit mAb	WB/IP
PTM-7003	Anti-RPA70 Rabbit mAb	WB
PTM-6433	Anti-RPL26L1 Rabbit mAb	WB/IHC-P
PTM-5496	Anti-RPL4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7068	Anti-RPS3 Rabbit mAb	WB
PTM-7421	Anti-RSK1 Rabbit mAb	WB/IHC-P
PTM-7096	Anti-S100B Rabbit pAb	WB/IHC-P
PTM-6933	Anti-SA2 Rabbit mAb	WB/IHC-P
PTM-7046	Anti-SAM68 Rabbit pAb	WB/IHC-P
PTM-5459	Anti-SBDS Mouse mAb	WB/IHC-P/ICC/IF
PTM-6326	Anti-SBDS Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7173	Anti-SCP3 Rabbit mAb	WB
PTM-7220	Anti-SEN1 Rabbit mAb	WB
PTM-7234	Anti-SET/TAF-I Rabbit mAb	WB/IHC-P
PTM-5347	Anti-SETBP1 Rabbit mAb	IHC-P
PTM-6972	Anti-SF2 Rabbit pAb	WB
PTM-6618	Anti-SF3B3 Rabbit mAb	WB/IHC-P/IP
PTM-7229	Anti-SFT Rabbit mAb	WB
PTM-5309	Anti-SGF29 Mouse mAb	WB
PTM-7458	Anti-SIN3A Rabbit mAb	WB/IP/IHC-P
PTM-5858	Anti-SIX2 Mouse mAb	WB
PTM-7160	Anti-SMC1A Rabbit mAb	WB
PTM-5786	Anti-SMYD2 Mouse mAb	WB/ICC/IF
PTM-7457	Anti-SMYD3 Rabbit mAb	WB
PTM-5406	Anti-SND1 Mouse mAb	WB
PTM-5393	Anti-SND1 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5902	Anti-SNRPN Mouse mAb	WB
PTM-7288	Anti-SOCS3 Rabbit mAb	WB
PTM-6434	Anti-SQSTM1/p62 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7215	Anti-SRSF3 Rabbit mAb	WB/IHC-P
PTM-6446	Anti-SSRP1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7211	Anti-StAR Rabbit mAb	WB/IHC-P
PTM-7210	Anti-STUB1 Rabbit mAb	WB/FC
PTM-6700	Anti-TAF1 Rabbit mAb	WB
PTM-7218	Anti-TAF15 Rabbit mAb	WB/IHC-P
PTM-5812	Anti-TAF6L Mouse mAb	WB/IHC-P
PTM-7186	Anti-Tau(phospho s396) Rabbit mAb	WB/IHC-P
PTM-6143	Anti-TAZ Mouse mAb	WB
PTM-7142	Anti-TBK1/ NAK Rabbit mAb	WB
PTM-6946	Anti-TBLR1/TBL1XR1 Rabbit mAb	WB/IHC-P
PTM-5469	Anti-TDP43 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6509	Anti-TDP43 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5998	Anti-TDRD3 Rabbit mAb	WB
PTM-5478	Anti-TELO2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7456	Anti-TFIIF beta Rabbit mAb	WB/IP/FC
PTM-5955	Anti-TH1L Mouse mAb	WB/ICC/IF
PTM-7196	Anti-Thymidine Kinase 1 Mouse mAb	WB
PTM-6589	Anti-Thymidylate Synthase Rabbit mAb	WB/ICC/IF/FC
PTM-6975	Anti-TIP49A Rabbit pAb	WB/IHC-P
PTM-7110	Anti-TOP3B Rabbit mAb	WB/ICC/IF
PTM-7467	Anti-Topoisomerase I Rabbit mAb	WB/IHC-P/FC/IP/ChIP
PTM-6388	Anti-Topoisomerase I Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7407	Anti-Topoisomerase II alpha Mouse mAb	WB/ICC/IF
PTM-7424	Anti-TORC1/CRTC1 Mouse mAb	WB/ICC/IF
PTM-7400	Anti-Tox Rabbit mAb	WB/IHC-P
PTM-6056	Anti-TPX2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6995	Anti-TREX1 Rabbit mAb	WB
PTM-7464	Anti-TRF2 Mouse mAb	WB/IHC-P/IP
PTM-7339	Anti-TRIAD1 Rabbit mAb	WB
PTM-5637	Anti-TRIM21 Mouse mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-7408	Anti-TRIM24 Rabbit mAb	WB/ICC/IF
PTM-6550	Anti-TRIM5 Mouse mAb	WB
PTM-7044	Anti-TRIM56 Rabbit pAb	WB
PTM-6754	Anti-TTF1 Rabbit mAb	WB/IHC-P
PTM-7030	Anti-TTK Mouse mAb	WB
PTM-5416	Anti-TYMS Mouse mAb	WB/IHC-P
PTM-5966	Anti-U2AF1 Mouse mAb	WB/IHC-P
PTM-7454	Anti-U2AF2 Rabbit mAb	WB/IHC-P/ChIP
PTM-6376	Anti-UBE2C Mouse mAb	WB/IHC-P/ICC/IF
PTM-7235	Anti-UBE2L3 Rabbit mAb	WB/IHC-P
PTM-7346	Anti-UBE3A Rabbit mAb	WB
PTM-7390	Anti-UHRF1 Mouse mAb	WB/IHC-P
PTM-5471	Anti-UPF1 Mouse mAb	WB/ICC/IF
PTM-7265	Anti-USP39 Rabbit mAb	WB/IHC-P
PTM-7465	Anti-USP7 Mouse mAb	WB/IP
PTM-6741	Anti-ValRS Mouse mAb	WB/IHC-P/ICC/IF
PTM-6839	Anti-WHSC1 Rabbit mAb	WB
PTM-7001	Anti-WSTF Rabbit mAb	WB
PTM-6351	Anti-XLF Mouse mAb	WB/IHC-P
PTM-5976	Anti-XRCC4 Mouse mAb	WB
PTM-5405	Anti-XRCC5 Mouse mAb	WB/ICC/IF
PTM-5968	Anti-XRN2 Mouse mAb	WB/IHC-P
PTM-7202	Anti-YAP1 Rabbit mAb	WB/IHC-P
PTM-6199	Anti-YTHDF1 Mouse mAb	WB
PTM-7291	Anti-ZPR1 Rabbit mAb	WB

Cancer Research Antibodies

PTM-7040	Anti-14-3-3 gamma Rabbit mAb	WB/IHC-P
PTM-6392	Anti-5T4 Rabbit mAb	WB/IHC-P/IP
PTM-6166	Anti-ABAT Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5996	Anti-ABCB11/BSEP Rabbit mAb	IHC-P
PTM-7060	Anti-ACADVL/LCAD Rabbit pAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6540	Anti-Acetyl-CoA Carboxylase 1 Mouse mAb	WB/IHC-P
PTM-185	Anti-Acetyl-p53 (Lys139) Mouse mAb	WB
PTM-6794	Anti-Acetyl-SOD2 (Lys68) Rabbit mAb	WB/IHC-P
PTM-7056	Anti-ACTH Rabbit mAb	IHC-P
PTM-7349	Anti-ADA Rabbit mAb	WB
PTM-6848	Anti-Adiponectin Rabbit mAb	WB/IHC-P
PTM-6692	Anti-ADK Rabbit mAb	WB/IHC-P
PTM-7208	Anti-ADRA1A Rabbit mAb	WB/FC/ICC/IF
PTM-6243	Anti-AFP Mouse mAb	WB/ICC/IF
PTM-6397	Anti-AFP Rabbit mAb	WB/IHC-P/FC
PTM-6322	Anti-AFP Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5229	Anti-AFP Rabbit mAb	WB/FC
PTM-6755	Anti-AGK Mouse mAb	WB
PTM-6853	Anti-AGR2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7253	Anti-AIF Mouse mAb	WB/IHC-P
PTM-6771	Anti-AIFM2/FSP1 Mouse mAb	WB/IHC-P
PTM-6735	Anti-AKT (phospho T308) Rabbit mAb	WB
PTM-6643	Anti-AKT1 + AKT2 + AKT3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6476	Anti-AKT1 + AKT2 + AKT3 Rabbit mAb	WB/IHC-P/FC
PTM-5191	Anti-AKT1 Mouse mAb	WB
PTM-6470	Anti-AKT1 Rabbit mAb	WB
PTM-6071	Anti-AKT1 Rabbit mAb	WB/IHC-P
PTM-6925	Anti-AKT2 Mouse mAb	WB
PTM-7179	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-5913	Anti-ALDH1L1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5790	Anti-ALDOA Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6729	Anti-Aldolase A Rabbit pAb	WB/IHC-P/ICC/IF
PTM-6730	Anti-Aldolase B Rabbit pAb	WB/IHC-P/ICC/IF
PTM-6078	Anti-Aldolase C Mouse mAb	WB/IHC-P/ICC/IF
PTM-5981	Anti-ALOX15 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-7397	Anti-alpha 1 Antitrypsin Mouse mAb	WB/IHC-P
PTM-5566	Anti-alpha Actinin 4 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5026	Anti-alpha Skeletal Muscle Actin Rabbit mAb	WB/IP
PTM-6511	Anti-alpha Skeletal Muscle Actin Rabbit mAb	WB/IHC-P/IP

Cat.#	Product	Recommended Applications
PTM-6721	Anti-Alpha-1-Antitrypsin Mouse mAb	WB/IHC-P/ICC/IF
PTM-5581	Anti-AMACR Rabbit mAb	WB/IHC-P
PTM-6739	Anti-AMPK alpha 1+AMPK alpha 2 Rabbit mAb	WB
PTM-6718	Anti-AMPK alpha 2 Mouse mAb	WB
PTM-6526	Anti-AMPK beta 2 Mouse mAb	WB/ICC/IF
PTM-6149	Anti-AMPK gamma 1 Mouse mAb	WB
PTM-5325	Anti-Amy2 Rabbit mAb	WB
PTM-5291	Anti-Androgen Receptor (AR) (C-term) Rabbit mAb	WB/IHC-P
PTM-5035	Anti-Androgen Receptor Rabbit mAb	WB/IHC-P
PTM-5423	Anti-Annexin A4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5424	Anti-Annexin V Mouse mAb	WB/IHC-P/ICC/IF
PTM-7149	Anti-APAF1 Rabbit mAb	WB/IHC-P
PTM-5730	Anti-Arginase-1 Rabbit mAb	WB/IHC-P
PTM-6812	Anti-Asparagine Synthetase Mouse mAb	WB/IHC-P
PTM-7396	Anti-ATF3 Rabbit mAb	WB/IHC-P
PTM-6650	Anti-ATG13 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7043	Anti-ATG16L1 Mouse mAb	WB
PTM-6296	Anti-ATG4B Rabbit mAb	WB/IHC-P/FC
PTM-6142	Anti-ATG5 Mouse mAb	WB/ICC/IF
PTM-5868	Anti-ATG5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6514	Anti-ATG9A Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5797	Anti-ATM Mouse mAb	WB
PTM-6394	Anti-ATM Rabbit mAb	WB/IHC-P
PTM-5946	Anti-ATM Rabbit mAb	WB
PTM-5260	Anti-ATM Rabbit mAb	IHC-P/ICC/IF
PTM-7372	Anti-ATP Citrate Lyase Rabbit mAb	WB/IHC-P
PTM-5170	Anti-ATR Mouse mAb	WB
PTM-6125	Anti-Aurora B Rabbit mAb	WB/IHC-P
PTM-5947	Anti-Aurora B Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6881	Anti-Axl Rabbit mAb	WB
PTM-5176	Anti-B7-H3 Mouse mAb	WB
PTM-5218	Anti-B7-H3/CD276 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6797	Anti-B7-H4 Rabbit mAb	WB
PTM-5420	Anti-BACH1 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5557	Anti-BAF60C Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6094	Anti-BAFF/TNFSF13B Mouse mAb	ICC/IF/FC
PTM-6371	Anti-BAG2 Rabbit mAb	WB
PTM-5781	Anti-BAG6 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5482	Anti-BAX Mouse mAb	WB/ICC/IF
PTM-5627	Anti-BCAR1 Mouse mAb	WB/IHC-P
PTM-6049	Anti-BCAS3 Mouse mAb	WB/ICC/IF
PTM-5777	Anti-Bcl-2 Rabbit mAb	WB/ICC/IF
PTM-5202	Anti-Bcl-2 Rabbit mAb	WB/IHC-P/IP
PTM-7008	Anti-Bcl-XL Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6793	Anti-Bcr Rabbit mAb	WB/IHC-P
PTM-6588	Anti-beta Catenin Rabbit pAb	WB/IHC-P
PTM-6608	Anti-beta-2-Microglobulin Rabbit mAb	WB/IHC-P
PTM-5353	Anti-BID Mouse mAb	WB/IHC-P
PTM-6111	Anti-Bid Rabbit mAb	WB/IHC-P/IP
PTM-7145	Anti-Bim Rabbit mAb	WB
PTM-5709	Anti-BMI1 Mouse mAb	WB
PTM-5205	Anti-BOB-1 Rabbit mAb	WB/IHC-P
PTM-6936	Anti-BTN2A2/BTF2 Rabbit mAb	WB
PTM-6487	Anti-c-Cbl Mouse mAb	WB/ICC/IF
PTM-6653	Anti-c-Fos Mouse mAb	WB
PTM-6457	Anti-c-IAP1 Mouse mAb	WB/ICC/IF
PTM-6686	Anti-c-IAP2 Rabbit mAb	WB/ICC/IF
PTM-6360	Anti-c-Maf Rabbit mAb	WB
PTM-5143	Anti-c-Myc Mouse mAb	WB
PTM-7317	Anti-c-Myc Rabbit mAb	WB/ChIP/IP
PTM-7216	Anti-C/EBP delta Mouse mAb	WB
PTM-6342	Anti-CAB39/MO25 Rabbit mAb	WB
PTM-6911	Anti-CAD Rabbit pAb	WB/IHC-P
PTM-5590	Anti-Cadherin-17 Rabbit mAb	IHC-P
PTM-5137	Anti-Cadherin-17 Rabbit mAb	IHC-P
PTM-7282	Anti-Cadherin-6 Rabbit mAb	WB
PTM-6403	Anti-Calcitonin Mouse mAb	WB
PTM-6690	Anti-Calcitonin Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-5378	Anti-Caldesmon/CDM Mouse mAb	WB/IHC-P
PTM-5261	Anti-Caldesmon/CDM Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6314	Anti-Calmodulin 1/2/3 Rabbit mAb	WB
PTM-6809	Anti-Calpain 2 Mouse mAb	WB
PTM-6806	Anti-Calretinin Rabbit pAb	WB/IHC-P
PTM-6299	Anti-CAMKII alpha Rabbit mAb	WB/IHC-P
PTM-6807	Anti-Carbonic Anhydrase 9 / CA9 Rabbit pAb	IHC-P
PTM-6661	Anti-Casein Kinase 1 gamma 2 Rabbit mAb	WB/ICC/IF/FC
PTM-7248	Anti-Caspase-2 Rabbit mAb	WB
PTM-6233	Anti-Caspase-6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6781	Anti-Caspase-7 Rabbit mAb	WB/IHC-P
PTM-6782	Anti-Caspase-7 Rabbit pAb	WB/IHC-P
PTM-6085	Anti-Caspase-8 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5630	Anti-Catalase Mouse mAb	WB/IHC-P
PTM-6372	Anti-Cathepsin B Rabbit mAb	WB
PTM-6369	Anti-Cathepsin D Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7436	Anti-Cathepsin K Mouse mAb	WB/IHC-P
PTM-6398	Anti-Caveolin-1 Rabbit mAb	WB
PTM-5056	Anti-Caveolin-1 Rabbit mAb	WB/IHC-P
PTM-5281	Anti-CBP Rabbit mAb	WB/ICC/IF
PTM-6980	Anti-CBP Rabbit pAb	WB
PTM-7203	Anti-CCR2 Rabbit mAb	IHC-P
PTM-5022	Anti-CCR4 Mouse mAb	FC
PTM-6489	Anti-Cdc23 Rabbit mAb	WB/ICC/IF
PTM-7410	Anti-Cdc25B Rabbit mAb	WB/IHC-P
PTM-6705	Anti-Cdc42 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5919	Anti-CDCP1 Mouse mAb	FC
PTM-6396	Anti-CDK2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5636	Anti-CDK4 Mouse mAb	WB
PTM-6131	Anti-CDK4 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5665	Anti-CDK6 Rabbit mAb	WB/ICC/IF/FC
PTM-5302	Anti-CDKN2A/p14ARF Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6704	Anti-CDKN2A/p16INK4a Rabbit mAb	WB
PTM-5924	Anti-CEACAM1 Mouse mAb	IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5446	Anti-CEACAM1 Rabbit mAb	WB/IHC-P/FC
PTM-7082	Anti-CGA/hCG alpha Rabbit mAb	WB/IHC-P
PTM-6276	Anti-Chk2 Mouse mAb	WB/IHC-P
PTM-5953	Anti-Chk2 Mouse mAb	WB
PTM-7204	Anti-Chromogranin A Rabbit mAb	IHC-P
PTM-5806	Anti-Chromogranin A Rabbit mAb	WB/IHC-P/FC
PTM-5586	Anti-Claudin-1 Rabbit mAb	WB/IHC-P
PTM-5679	Anti-Claudin-7 Rabbit mAb	WB/IHC-P
PTM-7246	Anti-Cleaved caspase-3 Rabbit mAb	WB
PTM-5272	Anti-Cleaved PARP1 p25 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6014	Anti-Cleaved PARP1 p85 Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6873	Anti-Clusterin beta Chain Rabbit mAb	WB
PTM-7261	Anti-Clusterin Rabbit mAb	WB/IHC-P
PTM-5729	Anti-Cofilin-1 Rabbit mAb	WB/ICC/IF/FC
PTM-7242	Anti-COL4A1 Rabbit mAb	WB
PTM-5217	Anti-Complement C4-B Rabbit mAb	WB/IHC-P
PTM-6528	Anti-COX IV Mouse mAb	WB/FC
PTM-5937	Anti-COX IV Mouse mAb	WB/IHC-P/ICC/IF
PTM-5908	Anti-COX IV Mouse mAb	WB/IHC-P/ICC/IF
PTM-7147	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-6951	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-6612	Anti-CPT1A Mouse mAb	WB/IHC-P
PTM-7069	Anti-CPT1B Rabbit mAb	WB/IHC-P
PTM-5502	Anti-CRABP2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7058	Anti-Creatine kinase B type Rabbit mAb	WB/IHC-P
PTM-6008	Anti-CSE1L Mouse mAb	WB/ICC/IF
PTM-5257	Anti-CSF-1R Rabbit mAb	IHC-P
PTM-5339	Anti-CTAG1B Rabbit mAb	IHC-P
PTM-5653	Anti-CTBP1 Mouse mAb	WB
PTM-5220	Anti-CTEN Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6452	Anti-Cullin 1 Mouse mAb	WB
PTM-6333	Anti-CX3CR1 Rabbit mAb	WB
PTM-6901	Anti-CXCR5 Rabbit pAb	WB
PTM-6659	Anti-Cyclin B1 Rabbit mAb	WB/ICC/IF/IP

Cat.#	Product	Recommended Applications
PTM-6029	Anti-Cyclin D1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5936	Anti-Cyclin D1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6664	Anti-Cyclin D2 Rabbit mAb	WB/IHC-P
PTM-5950	Anti-Cyclin D3 Mouse mAb	WB/ICC/IF
PTM-5138	Anti-Cyclin D3 Rabbit mAb	WB/IHC-P
PTM-5411	Anti-Cyclophilin B Mouse mAb	WB/IHC-P/ICC/IF
PTM-6685	Anti-Cyclophilin F Rabbit mAb	WB/IHC-P
PTM-5370	Anti-CYP11A1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7192	Anti-CYP17A1 Rabbit mAb	WB
PTM-5618	Anti-Cystatin A Rabbit mAb	WB/ICC/IF/FC
PTM-7181	Anti-Cytoplasmic Dynein Intermediate chain Mouse mAb	WB
PTM-7183	Anti-DcR1 Rabbit mAb	WB/ICC/IF
PTM-5059	Anti-delta 1 Catenin Rabbit mAb	IHC-P
PTM-6539	Anti-Desmin Rabbit mAb	WB/IHC-P
PTM-5188	Anti-Desmin Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6171	Anti-DHCR24 Mouse mAb	WB/ICC/IF
PTM-6150	Anti-DIO2 Mouse mAb	WB/IHC-P
PTM-6757	Anti-DJ-1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6256	Anti-DOPA Decarboxylase/DDC Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7109	Anti-Dopamine beta Hydroxylase Rabbit mAb	IHC-P
PTM-5673	Anti-DPEP1 Rabbit mAb	IHC-P
PTM-7057	Anti-DPP4/CD26 Rabbit mAb	WB/IHC-P
PTM-6695	Anti-DR5 Rabbit mAb	WB/IHC-P
PTM-5024	Anti-E Cadherin Mouse mAb	WB/IHC-P/ICC/IF
PTM-6222	Anti-E-Cadherin Rabbit mAb	WB/IP
PTM-6630	Anti-E2F1 Mouse mAb	WB
PTM-5428	Anti-ECHS1 Mouse mAb	WB/IHC-P
PTM-5789	Anti-eEF1A1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7114	Anti-eEF1A2 Rabbit mAb	WB/IHC-P
PTM-6843	Anti-eEF1A2 Rabbit pAb	WB/IHC-P
PTM-6745	Anti-EGFR (phospho S1166) Rabbit pAb	WB
PTM-5438	Anti-EGFR Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5042	Anti-EGFR Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5327	Anti-EHMT2 Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7360	Anti-EID1 Rabbit mAb	WB/IHC-P
PTM-5887	Anti-eIF4G2/DAP5 Mouse mAb	WB/ICC/IF
PTM-7348	Anti-ENO1 Rabbit mAb	WB
PTM-6780	Anti-ENO1 Rabbit pAb	WB/IHC-P
PTM-6021	Anti-ENO3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6731	Anti-ENO3 Rabbit pAb	WB/IHC-P
PTM-5280	Anti-EphA2 Rabbit mAb	WB/IHC-P
PTM-5017	Anti-ErbB2 Mouse mAb	WB
PTM-6304	Anti-ESRRB Mouse mAb	WB/ICC/IF
PTM-6967	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6838	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6792	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6816	Anti-Estrogen Related Receptor alpha Rabbit mAb	WB
PTM-6290	Anti-EXOC3L2 Mouse mAb	WB
PTM-6310	Anti-EZH2 Mouse mAb	WB/IHC-P
PTM-6272	Anti-EZH2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5227	Anti-EZH2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6364	Anti-Ezrin Rabbit mAb	WB/IHC-P/IP
PTM-5732	Anti-F4/80 Rabbit mAb	IHC-P
PTM-6416	Anti-FABP4 Rabbit mAb	WB/IHC-P
PTM-5574	Anti-Factor XIIIa Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6367	Anti-FADD (phospho S194) Rabbit mAb	WB
PTM-5717	Anti-FAK Mouse mAb	WB
PTM-5358	Anti-FAP Rabbit mAb	IHC-P
PTM-5444	Anti-FASN Rabbit mAb	WB/IHC-P/FC/IP
PTM-5269	Anti-FBXW7 Rabbit mAb	IHC-P/ICC/IF
PTM-5276	Anti-FGFR2 Rabbit mAb	IHC-P
PTM-6408	Anti-FGFR3 Rabbit mAb	WB/IP
PTM-6790	Anti-Fibronectin Rabbit pAb	WB/IHC-P
PTM-5942	Anti-FIS1 Mouse mAb	WB/ICC/IF
PTM-5582	Anti-FLI1 Rabbit mAb	WB
PTM-5743	Anti-FOXA1 Rabbit mAb	WB/IHC-P/FC
PTM-5666	Anti-FOXP1 Rabbit mAb	WB/IHC-P/IP
PTM-6229	Anti-FSH beta Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5688	Anti-FYN Mouse mAb	WB
PTM-6818	Anti-Galectin-9 Rabbit pAb	WB/IHC-P
PTM-6684	Anti-GAPDH Mouse mAb	WB/IHC-P/ICC/IF/FC
PTM-6420	Anti-Gastric Mucin-6/MUC6 Rabbit mAb	IHC-P/mIHC-P
PTM-7059	Anti-GBA Rabbit mAb	WB/IHC-P
PTM-5133	Anti-GCDFP-15 Rabbit mAb	IHC-P
PTM-5817	Anti-GCK Mouse mAb	WB/IHC-P/ICC/IF
PTM-7112	Anti-Glucagon Rabbit mAb	IHC-P
PTM-7446	Anti-Glucocorticoid receptor Rabbit mAb	WB/ICC/IF
PTM-6209	Anti-GLUT2 Rabbit mAb	WB/FC
PTM-6992	Anti-GLUT3+GLUT14 Rabbit mAb	WB
PTM-6775	Anti-Glutathione Peroxidase 4 Rabbit mAb	WB
PTM-6297	Anti-Glycogen synthase 1/GYS1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6920	Anti-Glypican 1 Rabbit mAb	WB
PTM-5669	Anti-Glypican 3 Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6751	Anti-Glypican-3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5959	Anti-GNB2 Rabbit mAb	WB/ICC/IF
PTM-5649	Anti-GNB2 Rabbit mAb	WB
PTM-5678	Anti-GOT1 Rabbit mAb	WB/IHC-P
PTM-6037	Anti-GOT2 Rabbit mAb	WB/IHC-P
PTM-6186	Anti-GPA33 Rabbit mAb	WB/IHC-P
PTM-6767	Anti-GPI Rabbit pAb	WB/IHC-P
PTM-5506	Anti-GPNMB Mouse mAb	WB/IHC-P/ICC/IF
PTM-5647	Anti-GPNMB Rabbit mAb	WB/IHC-P/FC
PTM-7326	Anti-Granulin Mouse mAb	WB/IHC-P
PTM-5515	Anti-GRP75 Mouse mAb	WB/IHC-P
PTM-6180	Anti-GRP75 Rabbit mAb	WB/ICC/IF
PTM-5727	Anti-GRP94 Mouse mAb	WB
PTM-5800	Anti-GRP94 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7313	Anti-GSK3 beta Mouse mAb	WB/IHC-P
PTM-7212	Anti-HACE1 Rabbit mAb	WB
PTM-7193	Anti-Hamartin Rabbit mAb	WB
PTM-6196	Anti-hCG beta Rabbit mAb	IHC-P/FC
PTM-6380	Anti-HEC1 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5225	Anti-hENT1 Rabbit mAb	IHC-P
PTM-6789	Anti-HER4/ErbB4 Rabbit mAb	IHC-P/ICC/IF
PTM-5851	Anti-HIF-1 alpha Mouse mAb	WB/IHC-P
PTM-5623	Anti-HIF-1 beta / ARNT Mouse mAb	WB
PTM-5371	Anti-HK2 Mouse mAb	WB
PTM-6303	Anti-HMGA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6018	Anti-HMGCR Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6937	Anti-HOXB13 Rabbit mAb	WB
PTM-6188	Anti-HSP27 Rabbit mAb	WB/ICC/IF/FC
PTM-5692	Anti-HSP70 Mouse mAb	WB
PTM-7134	Anti-Hsp70 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5949	Anti-IDH1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6878	Anti-IDH2 Rabbit pAb	WB/IHC-P
PTM-5083	Anti-IDO-1 Rabbit mAb	IHC-P
PTM-5082	Anti-IDO1 Rabbit mAb	IHC-P
PTM-6698	Anti-IFITM1 Rabbit pAb	WB
PTM-7150	Anti-IGF1 Receptor Rabbit mAb	WB/FC
PTM-6682	Anti-IGF2BP3 Mouse mAb	WB
PTM-6475	Anti-IGFBP3 Rabbit mAb	WB
PTM-6401	Anti-IkB alpha (phospho S32/36) Rabbit mAb	WB/IHC-P
PTM-6005	Anti-IkB alpha Mouse mAb	WB/IHC-P/ICC/IF
PTM-6274	Anti-IkB alpha Rabbit mAb	WB
PTM-7107	Anti-IKK alpha Rabbit mAb	WB
PTM-6546	Anti-IKK gamma Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6117	Anti-IL-12 Rabbit mAb	ICC/IF/FC
PTM-6235	Anti-IL-18 Mouse mAb	WB
PTM-6712	Anti-IL-4 Rabbit mAb	WB/ICC/IF/FC
PTM-6138	Anti-IL-9 Mouse mAb	FC
PTM-5898	Anti-ILK Mouse mAb	WB/IHC-P/ICC/IF
PTM-7379	Anti-IMPDPH2 Mouse mAb	WB/IHC-P
PTM-6170	Anti-IMPDPH2 Rabbit mAb	WB/ICC/IF/IP
PTM-6460	Anti-ING1 Rabbit mAb	WB
PTM-5993	Anti-Inhibin alpha Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6217	Anti-iNOS Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5226	Anti-iNOS Rabbit mAb	IHC-P
PTM-5510	Anti-Integrin alpha 3 / CD49c Mouse mAb	FC
PTM-5429	Anti-Integrin alpha 3 Mouse mAb	WB/IHC-P
PTM-5529	Anti-Integrin beta 3 / CD61 Mouse mAb	FC
PTM-6563	Anti-Interferon alpha-2 Rabbit mAb	WB
PTM-7135	Anti-JAK2 Rabbit mAb	WB
PTM-6456	Anti-KAP1 Mouse mAb	WB/IHC-P
PTM-6647	Anti-Kappa light chain Rabbit mAb	WB/IHC-P
PTM-5209	Anti-Kappa Light Chain Rabbit mAb	WB/IHC-P
PTM-6678	Anti-Keratin 10 Rabbit mAb	WB/IHC-P
PTM-5048	Anti-Keratin 10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5943	Anti-Keratin 14 Mouse mAb	IHC-P/ICC/IF
PTM-5391	Anti-Keratin 14 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5049	Anti-Keratin 14 Rabbit mAb	WB/IHC-P
PTM-5986	Anti-Keratin 15 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5067	Anti-Keratin 17 Rabbit mAb	WB/IHC-P/IP
PTM-5062	Anti-Keratin 18 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5645	Anti-Keratin 19 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5611	Anti-Keratin 20 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5454	Anti-Keratin 20 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5050	Anti-Keratin 20 Rabbit mAb	IHC-P
PTM-5012	Anti-Keratin 4 Rabbit mAb	IHC-P/mIHC-P
PTM-5004	Anti-Keratin 5 Mouse mAb	IHC-P
PTM-5875	Anti-Keratin 5 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5080	Anti-Keratin 5 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5040	Anti-Keratin 5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5492	Anti-Keratin 5/6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5079	Anti-Keratin 6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5298	Anti-Keratin 7 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5011	Anti-Keratin 7 Rabbit mAb	WB/IHC-P/IP
PTM-5107	Anti-Keratin 8 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5069	Anti-Keratin 8 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6058	Anti-Keratin 8 Rabbit mAb	WB/ICC/IF/FC
PTM-5491	Anti-Keratin HMW Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5490	Anti-Keratin LMW Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5032	Anti-Ki67 Rabbit mAb	IHC-P/mIHC-P/ICC/IF
PTM-5725	Anti-KLK11 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5425	Anti-KRAS Mouse mAb	WB/ICC/IF
PTM-6184	Anti-KRAS Rabbit mAb	WB/ICC/IF/IP
PTM-5287	Anti-KRT15 Mouse mAb	WB/ICC/IF
PTM-6244	Anti-Lactoferrin Mouse mAb	WB/ICC/IF
PTM-5190	Anti-Lambda Light Chain Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7332	Anti-Lamin B Receptor/LBR Mouse mAb	WB/IHC-P
PTM-5495	Anti-Lamin B1 Rabbit mAb	WB/IHC-P/mIHC-P
PTM-6357	Anti-Lamin B2 Rabbit mAb	WB/IHC-P/IP
PTM-6660	Anti-Laminin beta 1 Rabbit mAb	WB
PTM-5680	Anti-Laminin beta 1 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-6777	Anti-LATS1 Rabbit mAb	WB
PTM-6384	Anti-LC3B Rabbit mAb	WB/IHC-P
PTM-5869	Anti-LDHB Mouse mAb	WB/ICC/IF
PTM-6302	Anti-LDHC Rabbit mAb	WB
PTM-5861	Anti-LDL Receptor Mouse mAb	WB/ICC/IF
PTM-6258	Anti-LDL Receptor Rabbit mAb	WB/IHC-P/FC
PTM-6596	Anti-Legumain Mouse mAb	WB/IHC-P
PTM-7359	Anti-LIF Rat mAb	WB/IHC-P
PTM-7340	Anti-LIMK1 Rabbit mAb	WB
PTM-6651	Anti-LMO4 Rabbit mAb	WB
PTM-5603	Anti-LOX Rabbit mAb	WB/ICC/IF
PTM-6610	Anti-LRP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6024	Anti-LTA4H Rabbit mAb	WB
PTM-6619	Anti-Lumican Rabbit mAb	WB
PTM-5290	Anti-LYPD3 Rabbit mAb	IHC-P
PTM-5255	Anti-LYPD3 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7074	Anti-LYRIC/AEG1 Rabbit mAb	WB/IHC-P
PTM-6179	Anti-MAD2L1 Mouse mAb	WB/IHC-P
PTM-5516	Anti-MAD2L2 Mouse mAb	WB
PTM-5235	Anti-MAGEA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5234	Anti-MAGEA1 Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5896	Anti-MAP1LC3A Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6488	Anti-Maspin Mouse mAb	WB
PTM-5385	Anti-MDH1 Mouse mAb	WB
PTM-7310	Anti-MDM2 Mouse mAb	WB
PTM-7302	Anti-MED1 Rabbit mAb	WB
PTM-5983	Anti-MEK2 Rabbit mAb	WB/IHC-P
PTM-6890	Anti-MEK5 Rabbit mAb	WB
PTM-6153	Anti-MEK6/MKK6 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6537	Anti-Mesothelin Rabbit mAb	WB
PTM-5219	Anti-Mesothelin Rabbit mAb	WB/IP
PTM-7177	Anti-Met Mouse mAb	WB
PTM-5448	Anti-Met Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7185	Anti-MFF Rabbit mAb	WB
PTM-7151	Anti-Mitofusin2+Mitofusin1 Mouse mAb	WB/IHC-P
PTM-6512	Anti-MLH1 Rabbit mAb	WB
PTM-7144	Anti-MMP13 Mouse mAb	WB
PTM-7375	Anti-MMP7 Rabbit mAb	WB/IHC-P
PTM-6764	Anti-MMP8 Rabbit mAb	IHC-P
PTM-7259	Anti-MMP9 Rabbit mAb	WB/IHC-P
PTM-7404	Anti-MORF4L2 Rabbit mAb	WB/IHC-P
PTM-5342	Anti-MSH2 Rabbit mAb	IHC-P
PTM-5034	Anti-MSH2 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5196	Anti-MSH6 Mouse mAb	WB/IHC-P/FC
PTM-5014	Anti-MSH6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6594	Anti-mTOR Mouse mAb	WB/IHC-P/ICC/IF
PTM-5373	Anti-MUC1 Rabbit mAb	WB/IHC-P
PTM-6386	Anti-MUC13 Rabbit mAb	WB/ICC/IF
PTM-5995	Anti-MUC2 Rabbit mAb	WB/IHC-P/mIHC-P/FC
PTM-5372	Anti-MUC2 Rabbit pAb	IHC-P/ICC/IF
PTM-6334	Anti-Mucin5AC Mouse mAb	IHC-P/mIHC-P/ICC/IF
PTM-6905	Anti-MVP Rabbit pAb	WB
PTM-7079	Anti-Myoglobin Rabbit mAb	WB/IHC-P
PTM-5221	Anti-N-Cadherin Rabbit mAb	WB/IHC-P/IP
PTM-7356	Anti-N-Myc Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-5348	Anti-NaPi2b Rabbit mAb	IHC-P
PTM-7299	Anti-NCK1 Rabbit mAb	WB
PTM-7328	Anti-NCoR1 Rabbit mAb	WB
PTM-5402	Anti-NDRG2 Mouse mAb	WB/ICC/IF
PTM-7189	Anti-NDUFS3 Rabbit mAb	WB/IHC-P
PTM-5070	Anti-Nestin Rabbit mAb	IHC-P
PTM-6181	Anti-NEU3 Mouse mAb	WB/IHC-P
PTM-6787	Anti-Neutrophil Elastase Rabbit mAb	WB
PTM-5706	Anti-NF-κB1 p105/p50 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5755	Anti-NF-κB1 p105/p50 Rabbit mAb	WB/IHC-P
PTM-6593	Anti-NF-κB2 p100/p52 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5250	Anti-NKG2D Mouse mAb	FC
PTM-5330	Anti-NO66 Mouse mAb	WB/IHC-P
PTM-6996	Anti-NOX2/gp91phox Rabbit mAb	WB
PTM-5352	Anti-NPM1 Mouse mAb	WB
PTM-5511	Anti-NRF1 Mouse mAb	WB/ICC/IF
PTM-6449	Anti-NSE Mouse mAb	WB
PTM-5519	Anti-NSE Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6645	Anti-NSE Rabbit pAb	WB/FC
PTM-5716	Anti-NUP50 Mouse mAb	WB/IHC-P
PTM-5547	Anti-NUP62 Mouse mAb	WB/IHC-P
PTM-6051	Anti-OPCML Rabbit mAb	WB
PTM-5571	Anti-p130 Cas Rabbit mAb	WB/FC
PTM-7125	Anti-p21 Rabbit mAb	WB
PTM-6510	Anti-p21 Rabbit mAb	WB/IHC-P/IP
PTM-6471	Anti-p27 Kip1 Rabbit mAb	WB/IHC-P
PTM-6324	Anti-p44/42 MAPK (Erk1/2) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5850	Anti-p44/42 MAPK (Erk1/2) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6466	Anti-p53 Mouse mAb	WB
PTM-6517	Anti-p53 Rabbit mAb	WB
PTM-6319	Anti-p53 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5084	Anti-p53 Rabbit mAb	IHC-P
PTM-7260	Anti-p53R2 Rabbit mAb	WB
PTM-6257	Anti-p63 Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7099	Anti-Pan Trk Rabbit mAb	WB/IHC-P
PTM-6249	Anti-PAP Rabbit mAb	IHC-P
PTM-6015	Anti-PAP Rabbit mAb	WB/IHC-P/FC
PTM-5989	Anti-PARP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7431	Anti-PARP4 Rabbit mAb	WB
PTM-7249	Anti-Paxillin (Phospho Tyr118) Rabbit mAb	WB
PTM-6017	Anti-Paxillin Rabbit mAb	WB/IP
PTM-5650	Anti-PCK2 Rabbit mAb	WB/ICC/IF
PTM-6103	Anti-PD-L1 Rabbit mAb	WB/IHC-P
PTM-5075	Anti-PD-L1 Rabbit mAb	WB/IHC-P
PTM-6716	Anti-PD1 Mouse mAb	IHC-P
PTM-5445	Anti-PD1 Rabbit mAb	IHC-P/ICC/IF
PTM-5077	Anti-PD1 Rabbit mAb	IHC-P
PTM-6410	Anti-PDCD2 Rabbit mAb	WB
PTM-6482	Anti-PDCD4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6564	Anti-PDGFR beta Rabbit pAb	WB/IHC-P
PTM-6541	Anti-PDHA1 Rabbit mAb	WB
PTM-6321	Anti-PDK2 Rabbit mAb	WB
PTM-5967	Anti-PDP1 Rabbit mAb	WB/ICC/IF
PTM-6907	Anti-PDX1 Rabbit pAb	WB/IHC-P
PTM-7163	Anti-PEDF Rabbit mAb	WB
PTM-7066	Anti-PEG10 Rabbit mAb	WB
PTM-6723	Anti-Pepsinogen II/PGC Rabbit pAb	WB/IHC-P
PTM-7174	Anti-Perilipin-1 Rabbit mAb	WB
PTM-7027	Anti-Perilipin-2 Rabbit mAb	WB/ICC/IF
PTM-6699	Anti-periostin Mouse mAb	WB/IHC-P
PTM-6275	Anti-PGAM1 Mouse mAb	WB/ICC/IF
PTM-5888	Anti-PGAM1 Rabbit mAb	WB/ICC/IF
PTM-6861	Anti-PGAM1 Rabbit pAb	WB/IHC-P
PTM-6177	Anti-PGRMC1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6615	Anti-Phospho-(Ser/Thr) PDK1 Docking Motif Mouse mAb	WB
PTM-7428	Anti-Phospho-AMPK alpha 2 (Ser345) Rabbit mAb	WB
PTM-6785	Anti-Phospho-AMPK beta 1 (Ser182) Rabbit pAb	WB
PTM-7083	Anti-Phospho-Bcl-2 (Thr56) Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-7423	Anti-Phospho-beta Catenin (Thr41/Ser45) Rabbit mAb	WB
PTM-7209	Anti-Phospho-c-Myb (Ser11) Rabbit mAb	WB
PTM-6921	Anti-Phospho-c-Myc (Ser62) Rabbit mAb	WB
PTM-6887	Anti-Phospho-CAD (Ser1859) Rabbit pAb	WB
PTM-6855	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit mAb	WB
PTM-6801	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit pAb	WB/ICC/IF
PTM-6931	Anti-Phospho-CD31 (Tyr713) Rabbit mAb	WB
PTM-6955	Anti-Phospho-CDCP1 (Tyr743) Rabbit pAb	WB
PTM-7117	Anti-Phospho-Chk2 (Ser19) Rabbit pAb	WB/ICC/IF
PTM-7136	Anti-Phospho-Chk2 (Thr68) Rabbit mAb	WB
PTM-6185	Anti-Phospho-Cyclin D1 (Thr286) Rabbit mAb	WB/IHC-P/FC
PTM-6870	Anti-Phospho-EGFR (Tyr1092) Rabbit pAb	WB/ICC/IF
PTM-7373	Anti-Phospho-EGFR (Tyr1173) Rabbit mAb	WB/ICC/IF
PTM-6065	Anti-Phospho-FGFR1 (Tyr653/654) Mouse mAb	WB
PTM-6652	Anti-Phospho-HSP27 (Ser78) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7236	Anti-Phospho-IkB alpha (Tyr42) Rabbit mAb	WB
PTM-7130	Anti-Phospho-p38 MAPK (Thr180/Tyr182) Mouse mAb	WB/ICC/IF
PTM-6047	Anti-Phospho-p53 (Ser15) Mouse mAb	WB/ICC/IF
PTM-7425	Anti-Phospho-p53 (Thr377) Rabbit mAb	WB
PTM-7126	Anti-Phospho-PAK1 (Thr212) Rabbit mAb	WB
PTM-7381	Anti-Phospho-Rb(Thr826) Rabbit mAb	WB
PTM-6224	Anti-Phospho-RIPK1(Ser161) Mouse mAb	WB
PTM-7419	Anti-Phospho-STAT1 (Ser727) Rabbit mAb	WB/IHC-P
PTM-7024	Anti-Phospho-Stat3 (Ser727) Rabbit mAb	WB
PTM-5644	Anti-Phospho-Tau (Ser199) Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6714	Anti-Phospho-TBK1/NAK (Ser172) Rabbit mAb	WB/IP
PTM-7432	Anti-Phospho-Vimentin (Ser39) Rabbit mAb	WB/ICC/IF
PTM-7041	Anti-Phospho-Vimentin (Ser39) Rabbit mAb	WB
PTM-6211	Anti-Phospho-YAP1 (Ser127) Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5969	Anti-PI3 kinase p110 alpha Mouse mAb	WB
PTM-5228	Anti-PI3 Kinase p110 alpha Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5670	Anti-PI3 Kinase p85 alpha Rabbit mAb	IHC-P/ICC/IF
PTM-5426	Anti-PI3 Kinase p85 Beta Mouse mAb	WB/IHC-P
PTM-7293	Anti-Pim-1 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6291	Anti-PKC alpha Rabbit mAb	WB/ICC/IF
PTM-6206	Anti-PKC delta Rabbit mAb	WB/ICC/IF
PTM-7161	Anti-PKC mu Rabbit mAb	WB
PTM-5087	Anti-PKM2 Mouse mAb	WB/ICC/IF/FC
PTM-7010	Anti-PLAG1 Rabbit mAb	WB/IHC-P
PTM-7180	Anti-PLCG1 Rabbit mAb	WB
PTM-5952	Anti-PLK1 Mouse mAb	WB/ICC/IF
PTM-7416	Anti-PLK1 Rabbit mAb	WB
PTM-5659	Anti-PLK1 Rabbit mAb	WB
PTM-6458	Anti-PLOD3 Mouse mAb	WB
PTM-6262	Anti-PPP1CA Mouse mAb	WB/ICC/IF
PTM-6486	Anti-PPP1CB Rabbit mAb	WB/IHC-P
PTM-6128	Anti-PPP1CB Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5317	Anti-PPP4R1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6523	Anti-PRAME Rabbit mAb	WB
PTM-6525	Anti-PRDX3 Mouse mAb	WB/ICC/IF
PTM-6538	Anti-PRK2 Rabbit mAb	WB
PTM-6120	Anti-Prohibitin Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5846	Anti-PSA Rabbit mAb	WB/IHC-P
PTM-5747	Anti-PSA Rabbit pAb	WB/IHC-P
PTM-6242	Anti-PSA/KLK3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5757	Anti-PSMA/FOLH1 Rabbit mAb	WB/IHC-P/FC
PTM-5940	Anti-PTBP1 Rabbit mAb	WB
PTM-5748	Anti-PTEN Rabbit mAb	WB/IHC-P/FC
PTM-5119	Anti-PTEN Rabbit mAb	WB/IHC-P
PTM-7045	Anti-PYK2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6197	Anti-Rab11A Rabbit mAb	WB/ICC/IF
PTM-6336	Anti-Rac1 Mouse mAb	WB
PTM-7443	Anti-Rad21 Rabbit mAb	WB/IHC-P/ICC/IF/FC/ChIP
PTM-5956	Anti-RALB Mouse mAb	WB/ICC/IF
PTM-5624	Anti-RAP1A+RAP1B Rabbit mAb	WB/ICC/IF
PTM-6773	Anti-Ras Mouse mAb	WB/IHC-P
PTM-6246	Anti-REDD1 Mouse mAb	WB/ICC/IF
PTM-6657	Anti-Ret Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-7009	Anti-RING2/RING1B/RNF2 Rabbit pAb	WB
PTM-5259	Anti-RMP/C19orf2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6989	Anti-Robo1 Rabbit mAb	WB
PTM-6454	Anti-ROCK1 Mouse mAb	WB/IHC-P
PTM-6169	Anti-ROCK2 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6926	Anti-ROCK2+ROCK1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7303	Anti-RUNX1/AML1 Rabbit mAb	WB
PTM-7025	Anti-RXRA Rabbit mAb	WB/ICC/IF
PTM-7164	Anti-Ryanodine Receptor Rabbit mAb	IHC-P
PTM-6175	Anti-S100A4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6283	Anti-S100A8 Rabbit mAb	WB/IHC-P/FC
PTM-6343	Anti-S100A9 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5451	Anti-S100P Rabbit mAb	IHC-P/FC
PTM-7308	Anti-SAPK/JNK (phospho T183/Y185) Mouse mAb	WB
PTM-5293	Anti-SATB1 Rabbit mAb	IHC-P
PTM-5295	Anti-SATB2 Rabbit mAb	IHC-P
PTM-5132	Anti-SATB2 Rabbit mAb	IHC-P
PTM-6988	Anti-Scavenging Receptor SR-BI Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5909	Anti-SDC2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7234	Anti-SET/TAF-I Rabbit mAb	WB/IHC-P
PTM-5365	Anti-SETDB1 Mouse mAb	WB
PTM-6462	Anti-SIRT3 Rabbit mAb	WB
PTM-5683	Anti-SIRT3 Rabbit mAb	WB
PTM-6658	Anti-SIRT6 Rabbit mAb	WB
PTM-7231	Anti-SIRT7 Rabbit mAb	WB
PTM-6813	Anti-SKP2 Mouse mAb	WB
PTM-6973	Anti-SLC16A1 Rabbit mAb	WB/IHC-P
PTM-6971	Anti-SLC45A3 Rabbit mAb	IHC-P
PTM-6604	Anti-SLC45A3 Rabbit mAb	IHC-P
PTM-6127	Anti-Smac/Diablo Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5977	Anti-Smac/Diablo Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6929	Anti-Smad2 (phospho S255) Rabbit mAb	WB
PTM-7336	Anti-SMARCB1/BAF47 Rabbit mAb	WB
PTM-6688	Anti-Snail Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-7288	Anti-SOCS3 Rabbit mAb	WB
PTM-5681	Anti-SOD1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5876	Anti-SOX11 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-6736	Anti-SPARC Rat mAb	WB
PTM-6434	Anti-SQSTM1/p62 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5483	Anti-SQSTM1/p62 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5457	Anti-ST13 Rabbit mAb	WB/IHC-P/FC
PTM-7211	Anti-StAR Rabbit mAb	WB/IHC-P
PTM-6599	Anti-STAT1 Mouse mAb	WB
PTM-5754	Anti-STAT1 Rabbit mAb	WB/IHC-P/IP
PTM-5318	Anti-STAT5b Mouse mAb	WB/ICC/IF
PTM-5828	Anti-STAT5b Rabbit mAb	WB
PTM-6057	Anti-STAT6 Rabbit mAb	WB/ICC/IF/FC
PTM-6805	Anti-STAT6 Rabbit pAb	WB
PTM-5398	Anti-STK11 Mouse mAb	WB
PTM-5521	Anti-STK4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5838	Anti-STOML2 Mouse mAb	WB/ICC/IF
PTM-6742	Anti-Syntaxin 1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7214	Anti-Syntaxin 3 Rabbit mAb	WB/IHC-P
PTM-5764	Anti-TAK1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6946	Anti-TBLR1/TBL1XR1 Rabbit mAb	WB/IHC-P
PTM-5664	Anti-TCL1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6897	Anti-TGF beta 1 Rabbit mAb	WB/IHC-P
PTM-7029	Anti-TGF beta 2 Mouse mAb	WB/IHC-P
PTM-5773	Anti-TGF beta Receptor II Mouse mAb	WB/IHC-P/ICC/IF
PTM-5265	Anti-Thymidylate Synthase Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5667	Anti-Thyroglobulin Rabbit mAb	IHC-P/FC
PTM-6759	Anti-TIE2 Rabbit pAb	WB/IHC-P
PTM-7285	Anti-Tissue Plasminogen Activator Rabbit mAb	WB/IHC-P
PTM-6776	Anti-TM4SF1 Mouse mAb	FC
PTM-6835	Anti-TMEM192 Rabbit mAb	WB
PTM-5994	Anti-TMPRSS2 Rabbit mAb	WB/IHC-P/FC
PTM-6894	Anti-TMS1/ASC Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6834	Anti-TNF receptor II Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6388	Anti-Topoisomerase I Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7407	Anti-Topoisomerase II alpha Mouse mAb	WB/ICC/IF
PTM-5901	Anti-TPI1 Mouse mAb	WB/ICC/IF
PTM-5414	Anti-TPT1 Mouse mAb	WB/IHC-P
PTM-6056	Anti-TPX2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7354	Anti-TRADD Mouse mAb	WB/IHC-P
PTM-5534	Anti-TRAF2 Mouse mAb	WB/ICC/IF
PTM-6081	Anti-Transferrin Mouse mAb	WB/ICC/IF
PTM-6311	Anti-Transketolase Mouse mAb	WB/IHC-P/FC
PTM-6325	Anti-TrkA (phospho Y791) Rabbit mAb	WB
PTM-6852	Anti-TrkC Mouse mAb	WB
PTM-5344	Anti-TROP2 Rabbit mAb	WB/IHC-P
PTM-6238	Anti-TWEAKR/FN14 Mouse mAb	FC
PTM-5689	Anti-TXN Mouse mAb	WB/IHC-P/ICC/IF
PTM-6053	Anti-UCP1 Rabbit mAb	WB/IHC-P
PTM-6847	Anti-UQCRFS1/RISP Rabbit mAb	WB/IHC-P/IP
PTM-6689	Anti-VAV3 Rabbit mAb	WB
PTM-5086	Anti-VCAM1 Mouse mAb	IHC-P
PTM-7105	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-6623	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-6527	Anti-VDAC1/Porin Mouse mAb	WB/IHC-P
PTM-6157	Anti-VDAC1/Porin Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6935	Anti-VDAC1/Porin Rabbit pAb	WB/IHC-P
PTM-5278	Anti-VEGF Receptor2 Rabbit mAb	WB
PTM-5392	Anti-Villin Rabbit mAb	WB/IHC-P
PTM-5376	Anti-Vimentin Rabbit mAb	WB/IHC-P/ICC/IF/C/IP/mlIHC-P
PTM-5168	Anti-Vinculin Mouse mAb	WB/IHC-P
PTM-5224	Anti-Vinculin Rabbit mAb	WB/IHC-P
PTM-5361	Anti-VISTA Rabbit mAb	IHC-P
PTM-7086	Anti-Von Willebrand Factor Rabbit mAb	WB/IHC-P
PTM-5267	Anti-VPAC2 Rabbit mAb	IHC-P/ICC/IF
PTM-5561	Anti-WBP2 Mouse mAb	WB/IHC-P
PTM-6355	Anti-WDR19 Mouse mAb	WB
PTM-5814	Anti-WNT2 Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7298	Anti-WT1 Rabbit mAb	WB/IHC-P
PTM-7035	Anti-WT1 Rabbit pAb	WB
PTM-6039	Anti-WWOX Mouse mAb	WB/ICC/IF
PTM-7202	Anti-YAP1 Rabbit mAb	WB/IHC-P

Immunology Research Antibodies

PTM-5213	Anti-4-1BBL/TNFSF9 Mouse mAb	WB
PTM-5622	Anti-ABCE1 Rabbit mAb	WB/ICC/IF
PTM-7208	Anti-ADRA1A Rabbit mAb	WB/FC/ICC/IF
PTM-7200	Anti-AID Rabbit mAb	WB
PTM-7179	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-6831	Anti-ALPP Mouse mAb	WB/IHC-P
PTM-5325	Anti-Amy2 Rabbit mAb	WB
PTM-5176	Anti-B7-H3 Mouse mAb	WB
PTM-5218	Anti-B7-H3/CD276 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5275	Anti-B7-H3/CD276 Rabbit mAb	WB
PTM-6797	Anti-B7-H4 Rabbit mAb	WB
PTM-6094	Anti-BAFF/TNFSF13B Mouse mAb	ICC/IF/FC
PTM-7077	Anti-BANK1 Rabbit mAb	WB/IHC-P
PTM-6266	Anti-BAP29 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5332	Anti-Bcl-2 Mouse mAb	WB/ICC/IF
PTM-5202	Anti-Bcl-2 Rabbit mAb	WB/IHC-P/IP
PTM-5587	Anti-Bcl-2 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7145	Anti-Bim Rabbit mAb	WB
PTM-5205	Anti-BOB-1 Rabbit mAb	WB/IHC-P
PTM-6282	Anti-BOB-1/OBF-1 Rabbit mAb	WB/ICC/IF/IP
PTM-5419	Anti-BRD8 Mouse mAb	WB/IHC-P
PTM-5214	Anti-BTK Mouse mAb	WB
PTM-6936	Anti-BTN2A2/BTF2 Rabbit mAb	WB
PTM-6677	Anti-c-Kit Rabbit mAb	IHC-P
PTM-5675	Anti-c-Kit Rabbit mAb	IHC-P
PTM-6360	Anti-c-Maf Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6524	Anti-C1q Rabbit mAb	WB
PTM-7277	Anti-C5 Rabbit mAb	WB/IHC-P
PTM-6690	Anti-Calcitonin Rabbit pAb	WB
PTM-5635	Anti-CAPG Mouse mAb	WB/ICC/IF
PTM-5056	Anti-Caveolin-1 Rabbit mAb	WB/IHC-P
PTM-6398	Anti-Caveolin-1 Rabbit mAb	WB
PTM-7203	Anti-CCR2 Rabbit mAb	IHC-P
PTM-5736	Anti-CCR6 Rabbit mAb	WB/IHC-P
PTM-6301	Anti-CCR7 Rabbit mAb	WB/ICC/IF/FC
PTM-6076	Anti-CCT7 Mouse mAb	WB
PTM-5139	Anti-CD10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5078	Anti-CD10 Rabbit mAb	WB/IHC-P/FC/IP
PTM-5036	Anti-CD10 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5569	Anti-CD102/ICAM-2 Mouse mAb	FC
PTM-5052	Anti-CD103 Rabbit mAb	IHC-P
PTM-5331	Anti-CD115 Mouse mAb	WB
PTM-5536	Anti-CD11a Mouse mAb	ICC/IF/FC
PTM-6625	Anti-CD11b Mouse mAb	FC
PTM-5091	Anti-CD11b Rabbit mAb	IHC-P
PTM-5104	Anti-CD11b Rabbit mAb	IHC-P
PTM-6088	Anti-CD11b Rabbit mAb	IHC-P/ICC/IF/FC
PTM-6232	Anti-CD11b Rat mAb	FC
PTM-6875	Anti-CD11c Mouse mAb	WB/IHC-P
PTM-5126	Anti-CD13 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5456	Anti-CD13 Rabbit mAb	WB/IHC-P
PTM-5277	Anti-CD134/OX40 Rabbit mAb	IHC-P
PTM-7284	Anti-CD137 Rabbit mAb	WB
PTM-5101	Anti-CD138/Syndecan-1 Rabbit mAb	IHC-P
PTM-6644	Anti-CD14 Rabbit mAb	WB/IHC-P
PTM-5129	Anti-CD14 Rabbit mAb	WB/IHC-P
PTM-5114	Anti-CD140a Mouse mAb	FC
PTM-6393	Anti-CD146 Rabbit mAb	WB/IP
PTM-5175	Anti-CD147 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5615	Anti-CD147 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-20129	Anti-CD147 Rabbit mAb	WB
PTM-5526	Anti-CD148 Mouse mAb	FC
PTM-6467	Anti-CD15/SSEA1 Mouse mAb	IHC-P/ICC/IF/FC
PTM-6923	Anti-CD15/SSEA1 Mouse mAb	FC
PTM-5100	Anti-CD15/SSEA1 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5499	Anti-CD154 Mouse mAb	FC
PTM-5523	Anti-CD156c Mouse mAb	FC
PTM-5546	Anti-CD162/PSGL-1 Mouse mAb	ICC/IF/FC
PTM-6413	Anti-CD163 Rabbit mAb	WB/IHC-P
PTM-6747	Anti-CD163L1 Rabbit pAb	IHC-P
PTM-5607	Anti-CD168 Rabbit mAb	WB/ICC/IF/FC
PTM-5362	Anti-CD169 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5118	Anti-CD169 Rabbit mAb	IHC-P
PTM-5127	Anti-CD16a Rabbit mAb	WB/IHC-P
PTM-5103	Anti-CD16a Rabbit mAb	IHC-P
PTM-5567	Anti-CD18 Mouse mAb	FC
PTM-6626	Anti-CD180 Mouse mAb	FC
PTM-6402	Anti-CD19 (FITC) Mouse mAb	ICC/IF/FC
PTM-5932	Anti-CD19 (PE) Mouse mAb	FC
PTM-5479	Anti-CD19 Mouse mAb	FC
PTM-5240	Anti-CD19 Mouse mAb	FC
PTM-5731	Anti-CD19 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5322	Anti-CD19 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5441	Anti-CD19 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7091	Anti-CD1a (PE) Mouse mAb	FC
PTM-5201	Anti-CD1a Mouse mAb	FC
PTM-5015	Anti-CD1a Rabbit mAb	IHC-P/mIHC-P/FC
PTM-5503	Anti-CD1b Mouse mAb	FC
PTM-6711	Anti-CD1c Rabbit mAb	WB/IHC-P/IP
PTM-5906	Anti-CD1d Mouse mAb	ICC/IF/FC
PTM-5528	Anti-CD2 Mouse mAb	FC
PTM-5319	Anti-CD2 Rabbit mAb	IHC-P
PTM-6953	Anti-CD20 (APC) Mouse mAb	FC
PTM-6819	Anti-CD20 (FITC) Mouse mAb	FC

Cat.#	Product	Recommended Applications
PTM-6725	Anti-CD20 (PE) Mouse mAb	FC
PTM-6096	Anti-CD20 Mouse mAb	ICC/IF/FC
PTM-5437	Anti-CD20 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5037	Anti-CD20 Rabbit mAb	WB/IHC-P/FC/IP
PTM-6808	Anti-CD200 Mouse mAb	FC
PTM-6020	Anti-CD206 Mouse mAb	WB/IHC-P
PTM-5343	Anti-CD206 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5125	Anti-CD21 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5124	Anti-CD21 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6287	Anti-CD226 Mouse mAb	FC
PTM-5634	Anti-CD228 Mouse mAb	FC
PTM-5117	Anti-CD23 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5640	Anti-CD23 Rabbit mAb	IHC-P/ICC/IF
PTM-6255	Anti-CD239/BCAM Rabbit mAb	WB/IHC-P/FC
PTM-7120	Anti-CD24 (FITC) Rat mAb	FC
PTM-6554	Anti-CD24 Rat mAb	FC
PTM-5756	Anti-CD25 Rabbit mAb	ICC/IF/FC
PTM-5556	Anti-CD27 Mouse mAb	ICC/IF/FC
PTM-5538	Anti-CD270 Mouse mAb	FC
PTM-5772	Anti-CD272/BTLA Mouse mAb	FC
PTM-7116	Anti-CD272/BTLA Rabbit mAb	WB/ICC/IF
PTM-5055	Anti-CD278/ICOS Rabbit mAb	IHC-P
PTM-5540	Anti-CD28 Mouse mAb	ICC/IF/FC
PTM-6732	Anti-CD3 (PE) Rabbit mAb	ICC/IF/FC
PTM-7018	Anti-CD3 epsilon (APC) Mouse mAb	FC
PTM-5248	Anti-CD3 epsilon Mouse mAb	ICC/IF/FC
PTM-5013	Anti-CD3 epsilon Rabbit mAb	WB/IHC-P/mIHC-P/FC/IP
PTM-5643	Anti-CD31 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5031	Anti-CD31 Rabbit mAb	IHC-P/mIHC-P/FC
PTM-6225	Anti-CD31 Rat mAb	FC
PTM-7080	Anti-CD31(FITC) Rat mAb	FC
PTM-5827	Anti-CD32A Mouse mAb	FC
PTM-6108	Anti-CD32B Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6633	Anti-CD32B Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5207	Anti-CD33 Mouse mAb	WB
PTM-5189	Anti-CD33 Rabbit mAb	IHC-P
PTM-5336	Anti-CD34 Mouse mAb	IHC-P/FC
PTM-5811	Anti-CD34 Rabbit mAb	WB/IHC-P/mlIHC-P
PTM-5288	Anti-CD35 Rabbit mAb	IHC-P/ICC/IF
PTM-5646	Anti-CD35 Rabbit mAb	IHC-P
PTM-5128	Anti-CD35 Rabbit mAb	IHC-P
PTM-6898	Anti-CD36 Rabbit mAb	WB/IHC-P
PTM-5532	Anti-CD37 Mouse mAb	FC
PTM-7121	Anti-CD38 (APC) Mouse mAb	ICC/IF/FC
PTM-7019	Anti-CD38 (APC) Mouse mAb	FC
PTM-5931	Anti-CD38 (PE) Mouse mAb	ICC/IF/FC
PTM-5286	Anti-CD38 Mouse mAb	ICC/IF/FC
PTM-5239	Anti-CD38 Mouse mAb	FC
PTM-5662	Anti-CD38 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5102	Anti-CD38 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7232	Anti-CD39 Rabbit mAb	WB/IHC-P
PTM-7118	Anti-CD39 Rabbit mAb	WB
PTM-6559	Anti-CD3D Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7152	Anti-CD3E Rabbit mAb	WB
PTM-5488	Anti-CD3G Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5249	Anti-CD4 Mouse mAb	FC
PTM-7143	Anti-CD4 Mouse mAb	FC
PTM-5053	Anti-CD4 Rabbit mAb	WB/IHC-P/FC
PTM-6916	Anti-CD4(FITC) Mouse mAb	FC
PTM-7081	Anti-CD40 (FITC) Mouse mAb	FC
PTM-5537	Anti-CD40 Mouse mAb	ICC/IF/FC
PTM-6285	Anti-CD40 Rabbit mAb	WB/IHC-P
PTM-7078	Anti-CD40 Rabbit mAb	WB/ICC/IF/FC
PTM-6624	Anti-CD41 Mouse mAb	FC
PTM-6019	Anti-CD41 Rabbit mAb	WB
PTM-5524	Anti-CD42b Mouse mAb	FC
PTM-6674	Anti-CD42b Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7037	Anti-CD42b Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5237	Anti-CD42b Rabbit mAb	WB/IHC-P
PTM-5289	Anti-CD42b Rabbit mAb	WB/IHC-P
PTM-5541	Anti-CD43 Mouse mAb	ICC/IF/FC
PTM-5093	Anti-CD43 Rabbit mAb	WB/FC
PTM-6910	Anti-CD44 (PE) Rabbit mAb	ICC/IF/FC
PTM-5051	Anti-CD44 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6106	Anti-CD44 Rabbit mAb	WB/IHC-P
PTM-7071	Anti-CD45 (FITC) Mouse mAb	FC
PTM-6733	Anti-CD45 (PE) Rabbit mAb	ICC/IF/FC
PTM-5184	Anti-CD45 Mouse mAb	WB
PTM-5194	Anti-CD45 Mouse mAb	IHC-P/ICC/IF/FC
PTM-6584	Anti-CD45 Mouse mAb	IHC-P/ICC/IF/FC
PTM-6778	Anti-CD45 Mouse mAb	FC
PTM-5320	Anti-CD45 Rabbit mAb	WB/IHC-P/mlHC-P/ICC/IF/FC
PTM-5749	Anti-CD45 Rabbit mAb	WB/IHC-P
PTM-6974	Anti-CD47 Rabbit mAb	WB
PTM-5120	Anti-CD47 Rabbit mAb	IHC-P/FC
PTM-6230	Anti-CD48 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5507	Anti-CD49a Mouse mAb	FC
PTM-5544	Anti-CD5 Mouse mAb	FC
PTM-5568	Anti-CD5 Mouse mAb	ICC/IF/FC
PTM-5549	Anti-CD50/ICAM-3 Mouse mAb	FC
PTM-5807	Anti-CD50/ICAM3 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5481	Anti-CD55 Mouse mAb	FC
PTM-6144	Anti-CD57 Mouse mAb	IHC-P
PTM-6183	Anti-CD58 Rabbit mAb	WB/ICC/IF/FC
PTM-6803	Anti-CD58 Rabbit pAb	WB/ICC/IF
PTM-5686	Anti-CD59 Rabbit mAb	WB/IHC-P
PTM-6418	Anti-CD6/T12 Rabbit mAb	WB
PTM-5925	Anti-CD62L Mouse mAb	FC
PTM-5982	Anti-CD62L Rabbit mAb	WB/ICC/IF/FC
PTM-6908	Anti-CD62L Rabbit mAb	WB
PTM-6715	Anti-CD62P Mouse mAb	FC
PTM-5633	Anti-CD62P Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6798	Anti-CD63 Mouse mAb	WB
PTM-6863	Anti-CD63 Rabbit mAb	WB
PTM-6145	Anti-CD66b Mouse mAb	IHC-P
PTM-5185	Anti-CD68 Mouse mAb	IHC-P/ICC/IF
PTM-6382	Anti-CD68 Mouse mAb	IHC-P/ICC/IF
PTM-5130	Anti-CD68 Rabbit mAb	IHC-P/mIHC-P/FC
PTM-6495	Anti-CD69 Mouse mAb	ICC/IF/FC
PTM-5242	Anti-CD7 Mouse mAb	WB/IHC-P
PTM-6361	Anti-CD7 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6710	Anti-CD7 Rabbit mAb	WB
PTM-5066	Anti-CD7 Rabbit mAb	IHC-P/FC
PTM-6642	Anti-CD70 Rabbit mAb	WB
PTM-6455	Anti-CD71 Mouse mAb	WB/ICC/IF/FC
PTM-6294	Anti-CD71 Rabbit mAb	WB/IHC-P/FC
PTM-6940	Anti-CD72 Rabbit mAb	WB
PTM-5461	Anti-CD73 Mouse mAb	FC
PTM-6631	Anti-CD74 Rabbit mAb	WB/IHC-P/FC
PTM-5674	Anti-CD79a Rabbit mAb	WB/IHC-P
PTM-5616	Anti-CD79b Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5186	Anti-CD8 alpha Mouse mAb	IHC-P
PTM-5081	Anti-CD8 alpha Rabbit mAb	IHC-P/mIHC-P
PTM-6399	Anti-CD8 alpha Rabbit mAb	WB/ICC/IF
PTM-6268	Anti-CD8 alpha Rabbit mAb	IHC-P/FC
PTM-7132	Anti-CD8 beta Rabbit mAb	WB/IHC-P
PTM-7016	Anti-CD80 (APC) Mouse mAb	FC
PTM-7023	Anti-CD80 (FITC) Mouse mAb	FC
PTM-5247	Anti-CD80 Mouse mAb	IHC-P/FC
PTM-7021	Anti-CD81 (FITC) Mouse mAb	ICC/IF
PTM-6586	Anti-CD81 Mouse mAb	ICC/IF/FC
PTM-6533	Anti-CD81 Mouse mAb	ICC/IF/FC
PTM-6220	Anti-CD81 Rabbit mAb	WB
PTM-7076	Anti-CD81 Rabbit pAb	IHC-P/FC
PTM-6627	Anti-CD82 Mouse mAb	FC
PTM-7062	Anti-CD86 (APC) Mouse mAb	ICC/IF

Cat.#	Product	Recommended Applications
PTM-5570	Anti-CD86 Mouse mAb	ICC/IF/FC
PTM-5334	Anti-CD86 Mouse mAb	ICC/IF/FC
PTM-5553	Anti-CD89 Mouse mAb	ICC/IF/FC
PTM-6251	Anti-CD9 Rabbit mAb	WB/IHC-P/IP
PTM-6820	Anti-CD90 (FITC) Mouse mAb	FC
PTM-5710	Anti-CD90/Thy1 Mouse mAb	FC
PTM-6368	Anti-CD90/Thy1 Rabbit mAb	WB/IHC-P
PTM-7075	Anti-CD95 (FITC) Mouse mAb	FC
PTM-5484	Anti-CD95 Mouse mAb	ICC/IF/FC
PTM-5485	Anti-CD98 Mouse mAb	FC
PTM-5558	Anti-CD99 Mouse mAb	ICC/IF/FC
PTM-5349	Anti-CD99 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5991	Anti-CD99 Rabbit mAb	IHC-P/ICC/IF
PTM-5741	Anti-CDK8 Rabbit mAb	WB/IHC-P
PTM-5739	Anti-CDK9 Rabbit mAb	WB/IHC-P
PTM-5446	Anti-CEACAM1 Rabbit mAb	WB/IHC-P/FC
PTM-6377	Anti-CHD3 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6982	Anti-CHD4 Rabbit mAb	WB/IHC-P
PTM-6762	Anti-CHD4 Rabbit pAb	WB/IHC-P
PTM-7131	Anti-Claudin-11 Rabbit mAb	WB
PTM-7239	Anti-Claudin-11 Rabbit mAb	WB
PTM-7095	Anti-CLEC2D Rabbit mAb	WB
PTM-6873	Anti-Clusterin beta Chain Rabbit mAb	WB
PTM-7261	Anti-Clusterin Rabbit mAb	WB/IHC-P
PTM-6841	Anti-Complement C3 Rabbit pAb	WB
PTM-5217	Anti-Complement C4-B Rabbit mAb	WB/IHC-P
PTM-7330	Anti-Complement factor 8 beta/C8B Rabbit mAb	WB/IHC-P
PTM-6951	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-7147	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-6139	Anti-CRACC/SLAMF7/CD319 Mouse mAb	ICC/IF/FC
PTM-6115	Anti-CRACR2A Mouse mAb	WB
PTM-5810	Anti-CREB Rabbit mAb	WB/IHC-P
PTM-6411	Anti-CSF-1R Rabbit mAb	WB/ICC/IF
PTM-5257	Anti-CSF-1R Rabbit mAb	IHC-P

Cat.#	Product	Recommended Applications
PTM-7157	Anti-CTGF Rabbit mAb	WB
PTM-6333	Anti-CX3CR1 Rabbit mAb	WB
PTM-7153	Anti-CXCL7/PBP Rabbit mAb	WB
PTM-5771	Anti-CXCR4 Mouse mAb	FC
PTM-6901	Anti-CXCR5 Rabbit pAb	WB
PTM-6444	Anti-CYR61 Mouse mAb	WB/IHC-P
PTM-7183	Anti-DcR1 Rabbit mAb	WB/ICC/IF
PTM-7057	Anti-DPP4/CD26 Rabbit mAb	WB/IHC-P
PTM-7394	Anti-EBI3 Rabbit mAb	WB/IHC-P
PTM-5464	Anti-eIF3M/PCID1 Mouse mAb	WB/ICC/IF
PTM-5971	Anti-ETEA/U8XD8 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6027	Anti-EXOSC9 Mouse mAb	WB
PTM-6269	Anti-FER Rabbit mAb	WB
PTM-5387	Anti-FKBP52 Mouse mAb	WB/IHC-P
PTM-5016	Anti-FOXP3 Rabbit mAb	IHC-P
PTM-5700	Anti-G3BP1 Mouse mAb	WB/ICC/IF
PTM-6077	Anti-GATA3 Mouse mAb	WB/ICC/IF
PTM-6767	Anti-GPI Rabbit pAb	WB/IHC-P
PTM-5506	Anti-GPNMB Mouse mAb	WB/IHC-P/ICC/IF
PTM-5647	Anti-GPNMB Rabbit mAb	WB/IHC-P/FC
PTM-7326	Anti-Granulin Mouse mAb	WB/IHC-P
PTM-5246	Anti-Granzyme B Mouse mAb	WB
PTM-6146	Anti-HARS Mouse mAb	WB/ICC/IF
PTM-5262	Anti-HLA-DPA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5263	Anti-HLA-DPB1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-7070	Anti-HLA-DPB1 Rabbit pAb	WB
PTM-5578	Anti-HLA-DRA Rabbit mAb	IHC-P/FC
PTM-5517	Anti-HLA-DRB1 Rabbit mAb	WB/ICC/IF/FC
PTM-5841	Anti-HLA-E Mouse mAb	WB/IHC-P/ICC/IF
PTM-5208	Anti-Human alpha 1-Acid Glycoprotein Mouse mAb	WB
PTM-7378	Anti-Human alpha 1-Acid Glycoprotein Mouse mAb	WB/IHC-P
PTM-6187	Anti-Human IgD Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7108	Anti-Iba1 Rabbit mAb	IHC-P/mIHC-P
PTM-6749	Anti-Iba1/AIF-1 Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6415	Anti-ICAM-2 Rabbit mAb	WB/IHC-P/FC
PTM-5083	Anti-IDO-1 Rabbit mAb	IHC-P
PTM-5082	Anti-IDO1 Rabbit mAb	IHC-P
PTM-5792	Anti-IFIT3 Mouse mAb	WB
PTM-6698	Anti-IFITM1 Rabbit pAb	WB
PTM-5555	Anti-IgA Mouse mAb	WB/ICC/IF
PTM-6717	Anti-IGF2R/CI-M6PR Rabbit mAb	WB/IHC-P
PTM-6513	Anti-IGJ Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5222	Anti-IGJ Rabbit mAb	IHC-P
PTM-6095	Anti-IgM Mouse mAb	WB
PTM-6648	Anti-IgM Rabbit pAb	WB/ICC/IF/FC
PTM-6401	Anti-IkB alpha (phospho S32/36) Rabbit mAb	WB/IHC-P
PTM-6005	Anti-IkB alpha Mouse mAb	WB/IHC-P/ICC/IF
PTM-5020	Anti-IkB alpha Mouse mAb	WB/IHC-P
PTM-6274	Anti-IkB alpha Rabbit mAb	WB
PTM-6591	Anti-IkB beta Rabbit mAb	WB/ICC/IF/FC
PTM-5829	Anti-IKK alpha Rabbit mAb	WB/ICC/IF
PTM-7107	Anti-IKK alpha Rabbit mAb	WB
PTM-5701	Anti-IKK gamma Mouse mAb	WB/IHC-P
PTM-6546	Anti-IKK gamma Rabbit mAb	WB/ICC/IF/FC/IP
PTM-5116	Anti-IL-1 beta Mouse mAb	WB
PTM-6799	Anti-IL-1 beta Rabbit mAb	WB
PTM-5750	Anti-IL-10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5199	Anti-IL-10RB Mouse mAb	FC
PTM-6117	Anti-IL-12 Rabbit mAb	ICC/IF/FC
PTM-6040	Anti-IL-12RB1 Mouse mAb	WB
PTM-5596	Anti-IL-13RA2/CD123a2 Mouse mAb	ICC/IF/FC
PTM-6235	Anti-IL-18 Mouse mAb	WB
PTM-6597	Anti-IL-1RA Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6869	Anti-IL-22 Rabbit pAb	WB
PTM-7256	Anti-IL-26 Rabbit mAb	WB
PTM-6712	Anti-IL-4 Rabbit mAb	WB/ICC/IF/FC
PTM-5241	Anti-IL-6 Mouse mAb	WB
PTM-5212	Anti-IL-7R Mouse mAb	FC

Cat.#	Product	Recommended Applications
PTM-6359	Anti-IL-7R/CD127 Rabbit mAb	ICC/IF/FC
PTM-6138	Anti-IL-9 Mouse mAb	FC
PTM-6849	Anti-IL3RA/CD123 Rabbit pAb	WB/IHC-P
PTM-5226	Anti-iNOS Rabbit mAb	IHC-P
PTM-5527	Anti-Integrin alpha 1 / CD49a Mouse mAb	FC
PTM-5550	Anti-Integrin alpha 2 / CD49b Mouse mAb	FC
PTM-5510	Anti-Integrin alpha 3 / CD49c Mouse mAb	FC
PTM-5548	Anti-Integrin alpha 4 / CD49d Mouse mAb	FC
PTM-6563	Anti-Interferon alpha-2 Rabbit mAb	WB
PTM-5910	Anti-IPO7 Mouse mAb	WB/IHC-P
PTM-5422	Anti-IRAK1 Mouse mAb	WB
PTM-6696	Anti-IRAK1 Rabbit mAb	WB/IHC-P
PTM-5562	Anti-IRF3 Mouse mAb	WB
PTM-6237	Anti-IRF6 Rabbit mAb	WB/ICC/IF
PTM-6442	Anti-IRF8 Mouse mAb	WB/ICC/IF
PTM-6769	Anti-IRG1 Rabbit pAb	WB
PTM-6927	Anti-ISG20 Rabbit mAb	WB/IHC-P
PTM-7290	Anti-JMJD6 Mouse mAb	WB/IHC-P
PTM-6948	Anti-JNK2 Rabbit mAb	WB/IHC-P
PTM-6647	Anti-Kappa light chain Rabbit mAb	WB/IHC-P
PTM-5209	Anti-Kappa Light Chain Rabbit mAb	WB/IHC-P
PTM-5341	Anti-Keratin 16 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5340	Anti-Keratin 16 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6867	Anti-KIR3DL1 Rabbit mAb	WB/IHC-P
PTM-5279	Anti-Ku80 Rabbit mAb	IHC-P
PTM-5190	Anti-Lambda Light Chain Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6660	Anti-Laminin beta 1 Rabbit mAb	WB
PTM-6596	Anti-Legumain Mouse mAb	WB/IHC-P
PTM-7359	Anti-LIF Rat mAb	WB/IHC-P
PTM-6431	Anti-LILRB1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6024	Anti-LTA4H Rabbit mAb	WB
PTM-6264	Anti-Ly-6K Rabbit mAb	WB/FC
PTM-5658	Anti-LY75/DEC-205 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-6763	Anti-Lysozyme Rabbit mAb	WB/IHC-P/IP

Cat.#	Product	Recommended Applications
PTM-5235	Anti-MAGEA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5234	Anti-MAGEA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7286	Anti-MARCO Rabbit mAb	WB/IHC-P
PTM-7102	Anti-MAVS Rabbit pAb	WB
PTM-7115	Anti-MBL2 Rabbit mAb	WB/IHC-P
PTM-6153	Anti-MEK6/MKK6 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6172	Anti-MHC Class I Rabbit mAb	WB/IHC-P
PTM-6400	Anti-MHC class I Rabbit mAb	WB
PTM-5271	Anti-MUC4 Rabbit mAb	IHC-P
PTM-5487	Anti-Myeloperoxidase Rabbit mAb	WB/IP
PTM-5065	Anti-Myeloperoxidase Rabbit mAb	IHC-P
PTM-5305	Anti-NCAM1/CD56 Mouse mAb	WB/ICC/IF/FC
PTM-6532	Anti-Nectin-2/CD112 Mouse mAb	FC
PTM-6787	Anti-Neutrophil Elastase Rabbit mAb	WB
PTM-5706	Anti-NF-κB1 p105/p50 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5755	Anti-NF-κB1 p105/p50 Rabbit mAb	WB/IHC-P
PTM-5254	Anti-NFκB p65 Mouse mAb	WB
PTM-5591	Anti-NFκB p65 Mouse mAb	WB/ICC/IF
PTM-5250	Anti-NKG2D Mouse mAb	FC
PTM-6996	Anti-NOX2/gp91phox Rabbit mAb	WB
PTM-7392	Anti-OAS2 Mouse mAb	WB/ICC/IF
PTM-6902	Anti-Optineurin Rabbit mAb	WB
PTM-5920	Anti-ORM1/2 Mouse mAb	WB/ICC/IF
PTM-5033	Anti-PAX5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6663	Anti-PAX5 Rabbit mAb	WB/IHC-P/IP
PTM-6103	Anti-PD-L1 Rabbit mAb	WB/IHC-P
PTM-5075	Anti-PD-L1 Rabbit mAb	WB/IHC-P
PTM-6716	Anti-PD1 Mouse mAb	IHC-P
PTM-5445	Anti-PD1 Rabbit mAb	IHC-P/ICC/IF
PTM-5077	Anti-PD1 Rabbit mAb	IHC-P
PTM-6723	Anti-Pepsinogen II/PGC Rabbit pAb	WB/IHC-P
PTM-6387	Anti-PF4 Rabbit mAb	WB/IHC-P
PTM-7376	Anti-Phospho-BTK (Tyr551) Rabbit mAb	FC
PTM-7100	Anti-Phospho-CD19 (Tyr531) Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-6359	Anti-IL-7R/CD127 Rabbit mAb	ICC/IF/FC
PTM-6138	Anti-IL-9 Mouse mAb	FC
PTM-6849	Anti-IL3RA/CD123 Rabbit pAb	WB/IHC-P
PTM-5226	Anti-iNOS Rabbit mAb	IHC-P
PTM-5527	Anti-Integrin alpha 1 / CD49a Mouse mAb	FC
PTM-5550	Anti-Integrin alpha 2 / CD49b Mouse mAb	FC
PTM-5510	Anti-Integrin alpha 3 / CD49c Mouse mAb	FC
PTM-5548	Anti-Integrin alpha 4 / CD49d Mouse mAb	FC
PTM-6563	Anti-Interferon alpha-2 Rabbit mAb	WB
PTM-5910	Anti-IPO7 Mouse mAb	WB/IHC-P
PTM-5422	Anti-IRAK1 Mouse mAb	WB
PTM-6696	Anti-IRAK1 Rabbit mAb	WB/IHC-P
PTM-5562	Anti-IRF3 Mouse mAb	WB
PTM-6237	Anti-IRF6 Rabbit mAb	WB/ICC/IF
PTM-6442	Anti-IRF8 Mouse mAb	WB/ICC/IF
PTM-6769	Anti-IRG1 Rabbit pAb	WB
PTM-6927	Anti-ISG20 Rabbit mAb	WB/IHC-P
PTM-7290	Anti-JMJD6 Mouse mAb	WB/IHC-P
PTM-6948	Anti-JNK2 Rabbit mAb	WB/IHC-P
PTM-6647	Anti-Kappa light chain Rabbit mAb	WB/IHC-P
PTM-5209	Anti-Kappa Light Chain Rabbit mAb	WB/IHC-P
PTM-5341	Anti-Keratin 16 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5340	Anti-Keratin 16 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6867	Anti-KIR3DL1 Rabbit mAb	WB/IHC-P
PTM-5279	Anti-Ku80 Rabbit mAb	IHC-P
PTM-5190	Anti-Lambda Light Chain Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6660	Anti-Laminin beta 1 Rabbit mAb	WB
PTM-6596	Anti-Legumain Mouse mAb	WB/IHC-P
PTM-7359	Anti-LIF Rat mAb	WB/IHC-P
PTM-6431	Anti-LILRB1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6024	Anti-LTA4H Rabbit mAb	WB
PTM-6264	Anti-Ly-6K Rabbit mAb	WB/FC
PTM-5658	Anti-LY75/DEC-205 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-6763	Anti-Lysozyme Rabbit mAb	WB/IHC-P/IP

Cat.#	Product	Recommended Applications
PTM-5235	Anti-MAGEA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5234	Anti-MAGEA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7286	Anti-MARCO Rabbit mAb	WB/IHC-P
PTM-7102	Anti-MAVS Rabbit pAb	WB
PTM-7115	Anti-MBL2 Rabbit mAb	WB/IHC-P
PTM-6153	Anti-MEK6/MKK6 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6172	Anti-MHC Class I Rabbit mAb	WB/IHC-P
PTM-6400	Anti-MHC class I Rabbit mAb	WB
PTM-5271	Anti-MUC4 Rabbit mAb	IHC-P
PTM-5487	Anti-Myeloperoxidase Rabbit mAb	WB/IP
PTM-5065	Anti-Myeloperoxidase Rabbit mAb	IHC-P
PTM-5305	Anti-NCAM1/CD56 Mouse mAb	WB/ICC/IF/FC
PTM-6532	Anti-Nectin-2/CD112 Mouse mAb	FC
PTM-6787	Anti-Neutrophil Elastase Rabbit mAb	WB
PTM-5706	Anti-NF-κB1 p105/p50 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5755	Anti-NF-κB1 p105/p50 Rabbit mAb	WB/IHC-P
PTM-5254	Anti-NFκB p65 Mouse mAb	WB
PTM-5591	Anti-NFκB p65 Mouse mAb	WB/ICC/IF
PTM-5250	Anti-NKG2D Mouse mAb	FC
PTM-6996	Anti-NOX2/gp91phox Rabbit mAb	WB
PTM-7392	Anti-OAS2 Mouse mAb	WB/ICC/IF
PTM-6902	Anti-Optineurin Rabbit mAb	WB
PTM-5920	Anti-ORM1/2 Mouse mAb	WB/ICC/IF
PTM-5033	Anti-PAX5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6663	Anti-PAX5 Rabbit mAb	WB/IHC-P/IP
PTM-6103	Anti-PD-L1 Rabbit mAb	WB/IHC-P
PTM-5075	Anti-PD-L1 Rabbit mAb	WB/IHC-P
PTM-6716	Anti-PD1 Mouse mAb	IHC-P
PTM-5445	Anti-PD1 Rabbit mAb	IHC-P/ICC/IF
PTM-5077	Anti-PD1 Rabbit mAb	IHC-P
PTM-6723	Anti-Pepsinogen II/PGC Rabbit pAb	WB/IHC-P
PTM-6387	Anti-PF4 Rabbit mAb	WB/IHC-P
PTM-7376	Anti-Phospho-BTK (Tyr551) Rabbit mAb	FC
PTM-7100	Anti-Phospho-CD19 (Tyr531) Rabbit pAb	WB

Cat.#	Product	Recommended Applications
PTM-5751	Anti-Phospho-CREB (Ser133) Rabbit mAb	WB/IHC-P/FC
PTM-7236	Anti-Phospho-IkB alpha (Tyr42) Rabbit mAb	WB
PTM-7130	Anti-Phospho-p38 MAPK (Thr180/Tyr182) Mouse mAb	WB/ICC/IF
PTM-6714	Anti-Phospho-TBK1/NAK (Ser172) Rabbit mAb	WB/IP
PTM-6639	Anti-PI3 Kinase Class III Rabbit mAb	WB/ICC/IF
PTM-5969	Anti-PI3 kinase p110 alpha Mouse mAb	WB
PTM-6358	Anti-PI3 Kinase p110 delta Rabbit mAb	WB
PTM-6943	Anti-Placental Alkaline Phosphatase (PLAP) Rabbit pAb	WB/IHC-P
PTM-6337	Anti-PRDM1/Blimp1 Rabbit mAb	WB
PTM-5270	Anti-PU.1 Rabbit mAb	WB/IHC-P
PTM-5722	Anti-PVR/CD155 Mouse mAb	FC
PTM-5407	Anti-RAB27A Mouse mAb	WB
PTM-6336	Anti-Rac1 Mouse mAb	WB
PTM-6756	Anti-RAET1E Rabbit mAb	WB/ICC/IF/IP
PTM-6092	Anti-RANKL Mouse mAb	ICC/IF/FC
PTM-6061	Anti-RIG-I/DDX58 Mouse mAb	WB/ICC/IF
PTM-5394	Anti-RPL3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6343	Anti-S100A9 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7308	Anti-SAPK/JNK (phospho T183/Y185) Mouse mAb	WB
PTM-5132	Anti-SATB2 Rabbit mAb	IHC-P
PTM-6740	Anti-SERPINB9 Mouse mAb	WB
PTM-7234	Anti-SET/TAF-I Rabbit mAb	WB/IHC-P
PTM-6837	Anti-SHIP-1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7371	Anti-SHIP-2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5801	Anti-SIGIRR Rabbit mAb	WB
PTM-5559	Anti-Siglec-9 Mouse mAb	ICC/IF/FC
PTM-5554	Anti-SLAMF6/CD352 Mouse mAb	ICC/IF/FC
PTM-5854	Anti-SND1 Mouse mAb	WB/ICC/IF
PTM-7288	Anti-SOCS3 Rabbit mAb	WB
PTM-6483	Anti-ST2 Mouse mAb	WB
PTM-5090	Anti-STAT3 Mouse mAb	WB/IHC-P
PTM-5360	Anti-STING Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5359	Anti-STING Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5244	Anti-T-bet Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5912	Anti-TAB1 Mouse mAb	WB
PTM-5764	Anti-TAK1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7186	Anti-Tau(phospho s396) Rabbit mAb	WB/IHC-P
PTM-5877	Anti-TBX21 Rabbit mAb	WB
PTM-5357	Anti-Thrombomodulin Rabbit mAb	IHC-P/ICC/IF
PTM-5518	Anti-TIA1 Rabbit mAb	WB/FC/IP
PTM-5192	Anti-TLR4 Mouse mAb	WB
PTM-6834	Anti-TNF receptor II Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5617	Anti-TPSAB1 Rabbit mAb	IHC-P
PTM-5047	Anti-TRAF4 Mouse mAb	WB
PTM-6928	Anti-TRAF6 Rabbit mAb	WB/IHC-P
PTM-6081	Anti-Transferrin Mouse mAb	WB/ICC/IF
PTM-5346	Anti-TRBC1 Rabbit mAb	WB
PTM-6995	Anti-TREX1 Rabbit mAb	WB
PTM-5637	Anti-TRIM21 Mouse mAb	WB/IHC-P
PTM-6550	Anti-TRIM5 Mouse mAb	WB
PTM-5857	Anti-TRK Fused Gene Mouse mAb	WB/ICC/IF
PTM-5086	Anti-VCAM1 Mouse mAb	IHC-P
PTM-5361	Anti-VISTA Rabbit mAb	IHC-P
PTM-6189	Anti-Vitronectin/S-Protein Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5267	Anti-VPAC2 Rabbit mAb	IHC-P/ICC/IF
PTM-5836	Anti-Zap-70 Mouse mAb	WB/ICC/IF

Cardiovascular Research Antibodies

PTM-6942	Anti-AACT Rabbit mAb	WB/IHC-P
PTM-6906	Anti-ACAA2 Rabbit pAb	WB/IHC-P
PTM-6909	Anti-ACADS Rabbit pAb	WB/IHC-P
PTM-7060	Anti-ACADVL/VLCAD Rabbit pAb	WB/IHC-P
PTM-7158	Anti-ACE2 Rabbit mAb	WB/IHC-P
PTM-6542	Anti-ACE2 Rabbit pAb	WB/IHC-P
PTM-6540	Anti-Acetyl-CoA Carboxylase 1 Mouse mAb	WB/IHC-P
PTM-6794	Anti-Acetyl-SOD2 (Lys68) Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6766	Anti-Aconitase 1/ACO1 Rabbit mAb	WB
PTM-6004	Anti-ACTC1 Mouse mAb	WB/IHC-P
PTM-5417	Anti-ADH1B Mouse mAb	WB
PTM-6848	Anti-Adiponectin Rabbit mAb	WB/IHC-P
PTM-7208	Anti-ADRA1A Rabbit mAb	WB/FC/ICC/IF
PTM-6034	Anti-AFM Mouse mAb	WB/IHC-P/ICC/IF
PTM-6243	Anti-AFP Mouse mAb	WB/ICC/IF
PTM-6397	Anti-AFP Rabbit mAb	WB/IHC-P/FC
PTM-6322	Anti-AFP Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6735	Anti-AKT (phospho T308) Rabbit mAb	WB
PTM-6925	Anti-AKT2 Mouse mAb	WB
PTM-6072	Anti-Albumin Mouse mAb	WB/ICC/IF
PTM-7179	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-7397	Anti-alpha 1 Antitrypsin Mouse mAb	WB/IHC-P
PTM-5785	Anti-alpha Actinin 2 Mouse mAb	WB/IHC-P
PTM-6373	Anti-alpha Actinin 2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6721	Anti-Alpha-1-Antitrypsin Mouse mAb	WB/IHC-P/ICC/IF
PTM-6739	Anti-AMPK alpha 1+AMPK alpha 2 Rabbit mAb	WB
PTM-6526	Anti-AMPK beta 2 Mouse mAb	WB/ICC/IF
PTM-6149	Anti-AMPK gamma 1 Mouse mAb	WB
PTM-7271	Anti-APE1 Rabbit mAb	WB/IHC-P
PTM-6082	Anti-ApoA1 Mouse mAb	WB
PTM-7263	Anti-ApoH Mouse mAb	WB/IHC-P
PTM-7043	Anti-ATG16L1 Mouse mAb	WB
PTM-6296	Anti-ATG4B Rabbit mAb	WB/IHC-P/FC
PTM-6142	Anti-ATG5 Mouse mAb	WB/ICC/IF
PTM-5868	Anti-ATG5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6267	Anti-ATG7 Rabbit mAb	WB
PTM-6514	Anti-ATG9A Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6588	Anti-beta Catenin Rabbit pAb	WB/IHC-P
PTM-6016	Anti-beta-2 Adrenergic Receptor Rabbit mAb	WB/ICC/IF/FC
PTM-6608	Anti-beta-2-Microglobulin Rabbit mAb	WB/IHC-P
PTM-6709	Anti-BNP Rabbit pAb	WB
PTM-5900	Anti-BTN3A3 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-7314	Anti-C/EBP beta Mouse mAb	WB/IHC-P
PTM-7216	Anti-C/EBP delta Mouse mAb	WB
PTM-6208	Anti-Caldesmon/CDM Rabbit mAb	WB/IHC-P
PTM-6314	Anti-Calmodulin 1/2/3 Rabbit mAb	WB
PTM-6629	Anti-Calponin 1 Mouse mAb	WB
PTM-6126	Anti-Calponin 1 Rabbit mAb	WB/IHC-P/FC
PTM-6045	Anti-CAMKII alpha Mouse mAb	WB/ICC/IF
PTM-6299	Anti-CAMKII alpha Rabbit mAb	WB/IHC-P
PTM-6202	Anti-Carbonic Anhydrase 2/CA2 Rat mAb	WB
PTM-6807	Anti-Carbonic Anhydrase 9 / CA9 Rabbit pAb	IHC-P
PTM-6938	Anti-Caveolin-3 Rabbit mAb	WB/IHC-P
PTM-7203	Anti-CCR2 Rabbit mAb	IHC-P
PTM-5924	Anti-CEACAM1 Mouse mAb	IHC-P/ICC/IF
PTM-6873	Anti-Clusterin beta Chain Rabbit mAb	WB
PTM-7261	Anti-Clusterin Rabbit mAb	WB/IHC-P
PTM-7330	Anti-Complement factor 8 beta/C8B Rabbit mAb	WB/IHC-P
PTM-6543	Anti-Connexin 43 / GJA1 Rabbit mAb	WB/IHC-P
PTM-6951	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-6612	Anti-CPT1A Mouse mAb	WB/IHC-P
PTM-7069	Anti-CPT1B Rabbit mAb	WB/IHC-P
PTM-6411	Anti-CSF-1R Rabbit mAb	WB/ICC/IF
PTM-7157	Anti-CTGF Rabbit mAb	WB
PTM-6901	Anti-CXCR5 Rabbit pAb	WB
PTM-7192	Anti-CYP17A1 Rabbit mAb	WB
PTM-6444	Anti-CYR61 Mouse mAb	WB/IHC-P
PTM-6090	Anti-Cytoglobin Mouse mAb	WB/IHC-P/ICC/IF
PTM-6171	Anti-DHCR24 Mouse mAb	WB/ICC/IF
PTM-6695	Anti-DR5 Rabbit mAb	WB/IHC-P
PTM-5922	Anti-EPCR/CD201 Rat mAb	FC
PTM-6087	Anti-EPHB4 Rabbit mAb	WB/ICC/IF/FC
PTM-6168	Anti-FABP2 Mouse mAb	WB
PTM-6313	Anti-FABP3 Mouse mAb	WB/IHC-P
PTM-6416	Anti-FABP4 Rabbit mAb	WB/IHC-P
PTM-5865	Anti-FABP5 Mouse mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7247	Anti-Factor XII Rabbit mAb	WB
PTM-7439	Anti-Ferritin Heavy Chain Rabbit mAb	WB/IHC-P
PTM-7141	Anti-Ferritin Light Chain Rabbit mAb	WB
PTM-7266	Anti-ferritin light chain Rabbit mAb	WB
PTM-6746	Anti-Ferritin Light Chain Rabbit pAb	WB/IHC-P
PTM-6707	Anti-Fetuin B Mouse mAb	WB/IHC-P
PTM-6408	Anti-FGFR3 Rabbit mAb	WB/IP
PTM-6962	Anti-FGL1 Rabbit mAb	WB/IHC-P
PTM-6790	Anti-Fibronectin Rabbit pAb	WB/IHC-P
PTM-6697	Anti-Fibulin-5 Rabbit mAb	WB
PTM-6922	Anti-FTO Rabbit mAb	WB/IHC-P
PTM-7444	Anti-GATA1 Rabbit mAb	WB/IP/ChIP
PTM-6465	Anti-GATA2 Mouse mAb	WB
PTM-6297	Anti-Glycogen synthase 1/GYS1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5563	Anti-Glycophorin A Mouse mAb	IHC-P/FC
PTM-6107	Anti-GSK3 alpha Rabbit mAb	WB/ICC/IF/IP
PTM-6598	Anti-GST3/GSTP1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7140	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P
PTM-6691	Anti-GST3/GSTP1 Rabbit pAb	WB/IHC-P/ICC/IF
PTM-6478	Anti-Heme Oxygenase 1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6052	Anti-Hemoglobin alpha Rabbit mAb	WB/IHC-P
PTM-5258	Anti-Hemoglobin Alpha Rabbit mAb	IHC-P
PTM-6999	Anti-Hexokinase I Rabbit pAb	WB/IHC-P
PTM-7301	Anti-HSF1 Mouse mAb	WB
PTM-6188	Anti-HSP27 Rabbit mAb	WB/ICC/IF/FC
PTM-7378	Anti-Human alpha 1-Acid Glycoprotein Mouse mAb	WB/IHC-P
PTM-6475	Anti-IGFBP3 Rabbit mAb	WB
PTM-7107	Anti-IKK alpha Rabbit mAb	WB
PTM-6799	Anti-IL-1 beta Rabbit mAb	WB
PTM-5926	Anti-Integrin alpha V / CD51 Mouse mAb	ICC/IF/FC
PTM-6696	Anti-IRAK1 Rabbit mAb	WB/IHC-P
PTM-5723	Anti-JAM-A Mouse mAb	ICC/IF/FC
PTM-6244	Anti-Lactoferrin Mouse mAb	WB/ICC/IF
PTM-7382	Anti-LAP3 Mouse mAb	WB/IHC-P


Cat.#	Product	Recommended Applications
PTM-6384	Anti-LC3B Rabbit mAb	WB/IHC-P
PTM-5861	Anti-LDL Receptor Mouse mAb	WB/ICC/IF
PTM-6258	Anti-LDL Receptor Rabbit mAb	WB/IHC-P/FC
PTM-6610	Anti-LRP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6763	Anti-Lysozyme Rabbit mAb	WB/IHC-P/IP
PTM-5896	Anti-MAP1LC3A Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6977	Anti-MAPKAPK-2 Rabbit mAb	WB
PTM-7341	Anti-MEF2C Rabbit mAb	WB/IHC-P
PTM-6890	Anti-MEK5 Rabbit mAb	WB
PTM-7144	Anti-MMP13 Mouse mAb	WB
PTM-7375	Anti-MMP7 Rabbit mAb	WB/IHC-P
PTM-6764	Anti-MMP8 Rabbit mAb	IHC-P
PTM-7259	Anti-MMP9 Rabbit mAb	WB/IHC-P
PTM-6594	Anti-mTOR Mouse mAb	WB/IHC-P/ICC/IF
PTM-5718	Anti-Muscarinic Acetylcholine Receptor 2/CM2 Rabbit mAb	WB/IHC-P
PTM-5766	Anti-MYBPC3 Mouse mAb	WB/IHC-P
PTM-7079	Anti-Myoglobin Rabbit mAb	WB/IHC-P
PTM-7312	Anti-NeuroD1 Rabbit mAb	WB
PTM-6750	Anti-Neuropeptide Y Rabbit pAb	IHC-P
PTM-7324	Anti-Notch1 Rabbit mAb	WB/IHC-P
PTM-7225	Anti-Notch4 Rabbit mAb	WB
PTM-5974	Anti-OGDH Mouse mAb	WB/ICC/IF
PTM-6365	Anti-OLR1/LOX1 Rabbit mAb	WB/ICC/IF/FC
PTM-5920	Anti-ORM1/2 Mouse mAb	WB/ICC/IF
PTM-7067	Anti-PCSK9 Rabbit mAb	WB
PTM-6564	Anti-PDGFR beta Rabbit pAb	WB/IHC-P
PTM-7163	Anti-PEDF Rabbit mAb	WB
PTM-5780	Anti-PEPD Mouse mAb	WB/ICC/IF
PTM-7174	Anti-Perilipin-1 Rabbit mAb	WB
PTM-7027	Anti-Perilipin-2 Rabbit mAb	WB/ICC/IF
PTM-6699	Anti-periostin Mouse mAb	WB/IHC-P
PTM-6387	Anti-PF4 Rabbit mAb	WB/IHC-P
PTM-6177	Anti-PGRMC1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7428	Anti-Phospho-AMPK alpha 2 (Ser345) Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6785	Anti-Phospho-AMPK beta 1 (Ser182) Rabbit pAb	WB
PTM-7423	Anti-Phospho-beta Catenin (Thr41/Ser45) Rabbit mAb	WB
PTM-7429	Anti-Phospho-Cardiac Troponin I (Ser23/24) Rabbit mAb	WB
PTM-6855	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit mAb	WB
PTM-6801	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit pAb	WB/ICC/IF
PTM-6931	Anti-Phospho-CD31 (Tyr713) Rabbit mAb	WB
PTM-6065	Anti-Phospho-FGFR1 (Tyr653/654) Mouse mAb	WB
PTM-7123	Anti-Phospho-HSF1 (Ser326) Rabbit mAb	WB
PTM-6957	Anti-Phospho-HSF1 (Ser326) Rabbit pAb	WB
PTM-6652	Anti-Phospho-HSP27 (Ser78) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6896	Anti-Phospho-Myosin Light Chain 2 (Ser19) Mouse mAb	WB/ICC/IF
PTM-6193	Anti-Phospho-PTP1B (Ser352) Rabbit mAb	WB
PTM-6956	Anti-Phospho-PTP1B (Ser352) Rabbit pAb	WB
PTM-6223	Anti-Phospho-PTP1B (Ser378) Rabbit mAb	WB
PTM-6192	Anti-Phospho-PTP1B (Ser386) Rabbit mAb	WB/ICC/IF
PTM-7024	Anti-Phospho-Stat3 (Ser727) Rabbit mAb	WB
PTM-7178	Anti-phospholamban Mouse mAb	WB/IHC-P
PTM-6173	Anti-PITX2/RGS Rabbit mAb	IHC-P/FC
PTM-6291	Anti-PKC alpha Rabbit mAb	WB/ICC/IF
PTM-7180	Anti-PLCG1 Rabbit mAb	WB
PTM-6859	Anti-PLCG2 Rabbit mAb	WB
PTM-5705	Anti-PLEK Mouse mAb	WB/IHC-P/ICC/IF
PTM-6525	Anti-PRDX3 Mouse mAb	WB/ICC/IF
PTM-6522	Anti-PRDX6 Mouse mAb	WB
PTM-5552	Anti-Prealbumin/Transthyretin Mouse mAb	WB/IHC-P
PTM-7327	Anti-Prealbumin/Transthyretin Rabbit mAb	WB/IHC-P
PTM-6011	Anti-Prothrombin Mouse mAb	WB/IHC-P
PTM-6479	Anti-Prothrombin Rabbit mAb	WB/IHC-P
PTM-7297	Anti-PROX1 Mouse mAb	WB/IHC-P
PTM-5978	Anti-PSMA1 Rabbit mAb	WB/IP
PTM-6477	Anti-PTP1B Rabbit mAb	WB/IP
PTM-6344	Anti-PTP1B Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5626	Anti-PTPN2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6092	Anti-RANKL Mouse mAb	ICC/IF/FC


Cat.#	Product	Recommended Applications
PTM-6240	Anti-RBP4 Mouse mAb	WB/ICC/IF
PTM-6246	Anti-REDD1 Mouse mAb	WB/ICC/IF
PTM-7164	Anti-Ryanodine Receptor Rabbit mAb	IHC-P
PTM-6175	Anti-S100A4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5113	Anti-S100B Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5918	Anti-Serpin C1 / Antithrombin-III Mouse mAb	WB
PTM-6105	Anti-Serum Amyloid P/SAP Rabbit mAb	WB/IHC-P
PTM-7371	Anti-SHIP-2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6813	Anti-SKP2 Mouse mAb	WB
PTM-6688	Anti-Snail Mouse mAb	WB

Selected Antibodies for Cardiovascular Biomarkers


Anti-FABP3 Mouse mAb PTM-6313



WB sample: 1.Mouse heart tissue 2.U-87 MG cells




IHC sample: Human myocardium




IHC sample: Human skeletal


Anti-Myoglobin Rabbit mAb PTM-7079



WB sample: Human heart tissue



IHC sample: Human myocardium



Negative control tissue: Human colon

Cat.#	Product	Recommended Applications
PTM-6736	Anti-SPARC Rat mAb	WB
PTM-5744	Anti-SPTB Rabbit mAb	WB
PTM-6434	Anti-SQSTM1/p62 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6483	Anti-ST2 Mouse mAb	WB
PTM-7211	Anti-StAR Rabbit mAb	WB/IHC-P
PTM-3126	Anti-SUCLA2 Mouse mAb	WB
PTM-5764	Anti-TAK1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7186	Anti-Tau(phospho s396) Rabbit mAb	WB/IHC-P
PTM-6897	Anti-TGF beta 1 Rabbit mAb	WB/IHC-P
PTM-6295	Anti-Thrombomodulin Rabbit mAb	WB/ICC/IF
PTM-5542	Anti-Thrombospondin 1 Mouse mAb	WB
PTM-6759	Anti-TIE2 Rabbit pAb	WB/IHC-P
PTM-6779	Anti-Tissue Factor Rabbit pAb	WB
PTM-7285	Anti-Tissue Plasminogen Activator Rabbit mAb	WB/IHC-P
PTM-5880	Anti-TMEM16A Rabbit mAb	IHC-P/FC
PTM-6834	Anti-TNF receptor II Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6200	Anti-TNFAIP1 Mouse mAb	WB/ICC/IF
PTM-7354	Anti-TRADD Mouse mAb	WB/IHC-P
PTM-6238	Anti-TWEAKR/FN14 Mouse mAb	FC
PTM-7238	Anti-ULK1 (phospho S623) Rabbit mAb	WB
PTM-6998	Anti-USF1 Rabbit mAb	WB
PTM-6623	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-7105	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-6939	Anti-Vinculin Rabbit mAb	WB/IHC-P
PTM-6609	Anti-VISTA Rabbit mAb	WB
PTM-6189	Anti-Vitronectin/S-Protein Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7086	Anti-Von Willebrand Factor Rabbit mAb	WB/IHC-P

Metabolism Research Antibodies

PTM-5620	Anti-4E-BP1 Mouse mAb	WB/ICC/IF
PTM-6823	Anti-58K Golgi protein Rabbit pAb	WB/IHC-P
PTM-6166	Anti-ABAT Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5095	Anti-ABCA1 Mouse mAb	WB
PTM-5996	Anti-ABCB11/BSEP Rabbit mAb	IHC-P
PTM-6906	Anti-ACAA2 Rabbit pAb	WB/IHC-P
PTM-5778	Anti-ACADM Rabbit mAb	WB/ICC/IF
PTM-6909	Anti-ACADS Rabbit pAb	WB/IHC-P
PTM-7060	Anti-ACADVL/VLCAD Rabbit pAb	WB/IHC-P
PTM-6540	Anti-Acetyl-CoA Carboxylase 1 Mouse mAb	WB/IHC-P
PTM-6794	Anti-Acetyl-SOD2 (Lys68) Rabbit mAb	WB/IHC-P
PTM-6766	Anti-Aconitase 1/ACO1 Rabbit mAb	WB
PTM-6983	Anti-Aconitase 2 Rabbit pAb	WB
PTM-7056	Anti-ACTH Rabbit mAb	IHC-P
PTM-6520	Anti-ADH4 Mouse mAb	WB
PTM-6848	Anti-Adiponectin Rabbit mAb	WB/IHC-P
PTM-6692	Anti-ADK Rabbit mAb	WB/IHC-P
PTM-7208	Anti-ADRA1A Rabbit mAb	WB/FC/ICC/IF
PTM-6034	Anti-AFM Mouse mAb	WB/IHC-P/ICC/IF
PTM-6148	Anti-AKR7A2 Mouse mAb	WB/ICC/IF
PTM-6735	Anti-AKT (phospho T308) Rabbit mAb	WB
PTM-6071	Anti-AKT1 Rabbit mAb	WB/IHC-P
PTM-6470	Anti-AKT1 Rabbit mAb	WB
PTM-6925	Anti-AKT2 Mouse mAb	WB
PTM-6072	Anti-Albumin Mouse mAb	WB/ICC/IF
PTM-5794	Anti-Albumin Mouse mAb	WB
PTM-5345	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-7179	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-5337	Anti-ALDH3A1 Rabbit mAb	WB/IHC-P
PTM-6091	Anti-ALDH5A1/SSADH Mouse mAb	WB/IHC-P
PTM-6646	Anti-ALDH7A1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6811	Anti-ALDH9A1 Mouse mAb	WB/IHC-P
PTM-5790	Anti-ALDOA Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6729	Anti-Aldolase A Rabbit pAb	WB/IHC-P/ICC/IF
PTM-6730	Anti-Aldolase B Rabbit pAb	WB/IHC-P/ICC/IF
PTM-6078	Anti-Aldolase C Mouse mAb	WB/IHC-P/ICC/IF
PTM-5981	Anti-ALOX15 Rabbit mAb	IHC-P/ICC/IF/FC

Cat.#	Product	Recommended Applications
PTM-5566	Anti-alpha Actinin 4 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6721	Anti-Alpha-1-Antitrypsin Mouse mAb	WB/IHC-P/ICC/IF
PTM-5581	Anti-AMACR Rabbit mAb	WB/IHC-P
PTM-6739	Anti-AMPK alpha 1+AMPK alpha 2 Rabbit mAb	WB
PTM-6718	Anti-AMPK alpha 2 Mouse mAb	WB
PTM-5832	Anti-AMPK beta 1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6526	Anti-AMPK beta 2 Mouse mAb	WB/ICC/IF
PTM-6149	Anti-AMPK gamma 1 Mouse mAb	WB
PTM-5501	Anti-ANGPTL3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7149	Anti-APAF1 Rabbit mAb	WB/IHC-P
PTM-5307	Anti-ApoA1 Mouse mAb	WB/IHC-P
PTM-6082	Anti-ApoA1 Mouse mAb	WB
PTM-5915	Anti-ApoE Mouse mAb	WB/IHC-P/ICC/IF
PTM-5921	Anti-Arginase-1 Mouse mAb	WB/ICC/IF
PTM-6826	Anti-Arginase-1 Rabbit pAb	WB/IHC-P
PTM-6028	Anti-ASL Mouse mAb	WB
PTM-6812	Anti-Asparagine Synthetase Mouse mAb	WB/IHC-P
PTM-6650	Anti-ATG13 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7043	Anti-ATG16L1 Mouse mAb	WB
PTM-5606	Anti-ATG4B Rabbit mAb	WB/ICC/IF
PTM-6296	Anti-ATG4B Rabbit mAb	WB/IHC-P/FC
PTM-6142	Anti-ATG5 Mouse mAb	WB/ICC/IF
PTM-5868	Anti-ATG5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6267	Anti-ATG7 Rabbit mAb	WB
PTM-6514	Anti-ATG9A Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5769	Anti-ATP5A1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5163	Anti-ATP5A1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6535	Anti-ATP5D Mouse mAb	WB/IHC-P
PTM-5468	Anti-ATP5O Mouse mAb	WB/ICC/IF
PTM-5600	Anti-ATPB Mouse mAb	WB/ICC/IF
PTM-5009	Anti-ATPB Mouse mAb	WB/IHC-P/ICC/IF
PTM-6371	Anti-BAG2 Rabbit mAb	WB
PTM-5482	Anti-BAX Mouse mAb	WB/ICC/IF
PTM-5777	Anti-Bcl-2 Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7008	Anti-Bcl-XL Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6485	Anti-BCS1L Mouse mAb	WB
PTM-5440	Anti-BHMT Rabbit mAb	WB
PTM-5353	Anti-BID Mouse mAb	WB/IHC-P
PTM-6111	Anti-Bid Rabbit mAb	WB/IHC-P/IP
PTM-7145	Anti-Bim Rabbit mAb	WB
PTM-5843	Anti-BMAL1 Rabbit mAb	WB/IHC-P/FC
PTM-6457	Anti-c-IAP1 Mouse mAb	WB/ICC/IF
PTM-6686	Anti-c-IAP2 Rabbit mAb	WB/ICC/IF
PTM-7216	Anti-C/EBP delta Mouse mAb	WB
PTM-6342	Anti-CAB39/MO25 Rabbit mAb	WB
PTM-6911	Anti-CAD Rabbit pAb	WB/IHC-P
PTM-6314	Anti-Calmodulin 1/2/3 Rabbit mAb	WB
PTM-6809	Anti-Calpain 2 Mouse mAb	WB
PTM-5489	Anti-Calpain small subunit 1 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5505	Anti-Calsequestrin 2 Mouse mAb	WB/IHC-P
PTM-6045	Anti-CAMKII alpha Mouse mAb	WB/ICC/IF
PTM-6299	Anti-CAMKII alpha Rabbit mAb	WB/IHC-P
PTM-6807	Anti-Carbonic Anhydrase 9 / CA9 Rabbit pAb	IHC-P
PTM-6661	Anti-Casein Kinase 1 gamma 2 Rabbit mAb	WB/ICC/IF/FC
PTM-5752	Anti-Caspase-3 Rabbit mAb	WB/IHC-P
PTM-6233	Anti-Caspase-6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6781	Anti-Caspase-7 Rabbit mAb	WB/IHC-P
PTM-6782	Anti-Caspase-7 Rabbit pAb	WB/IHC-P
PTM-6085	Anti-Caspase-8 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5630	Anti-Catalase Mouse mAb	WB/IHC-P
PTM-5972	Anti-CFTR Mouse mAb	IHC-P
PTM-5848	Anti-CGA/hCG alpha Rabbit mAb	WB/IHC-P
PTM-7082	Anti-CGA/hCG alpha Rabbit mAb	WB/IHC-P
PTM-6765	Anti-CGA/hCG alpha Rabbit pAb	WB/IHC-P
PTM-6774	Anti-Citrate synthase Mouse mAb	WB/IHC-P/ICC/IF
PTM-5833	Anti-CLPP Rabbit mAb	WB
PTM-5494	Anti-COMMD1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6675	Anti-COX IV (HRP conjugate) Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5937	Anti-COX IV Mouse mAb	WB/IHC-P/ICC/IF
PTM-5908	Anti-COX IV Mouse mAb	WB/IHC-P/ICC/IF
PTM-6528	Anti-COX IV Mouse mAb	WB/FC
PTM-6421	Anti-COX IV Rabbit mAb	WB/IHC-P/FC/IP
PTM-6951	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-7147	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-6891	Anti-CPS1 Rabbit mAb	WB/IHC-P
PTM-6612	Anti-CPT1A Mouse mAb	WB/IHC-P
PTM-7069	Anti-CPT1B Rabbit mAb	WB/IHC-P
PTM-7058	Anti-Creatine kinase B type Rabbit mAb	WB/IHC-P
PTM-6245	Anti-CTPS2 Mouse mAb	WB/ICC/IF
PTM-5771	Anti-CXCR4 Mouse mAb	FC
PTM-6685	Anti-Cyclophilin F Rabbit mAb	WB/IHC-P
PTM-5370	Anti-CYP11A1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7181	Anti-Cytoplasmic Dynein Intermediate chain Mouse mAb	WB
PTM-7104	Anti-DARPP-32 Rabbit mAb	WB/IHC-P
PTM-6171	Anti-DHCR24 Mouse mAb	WB/ICC/IF
PTM-6708	Anti-DHFR Mouse mAb	WB
PTM-6150	Anti-DIO2 Mouse mAb	WB/IHC-P
PTM-6757	Anti-DJ-1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6256	Anti-DOPA Decarboxylase/DDC Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7109	Anti-Dopamine beta Hydroxylase Rabbit mAb	IHC-P
PTM-7057	Anti-DPP4/CD26 Rabbit mAb	WB/IHC-P
PTM-6783	Anti-DPS1 Rabbit pAb	WB
PTM-5428	Anti-ECHS1 Mouse mAb	WB/IHC-P
PTM-6780	Anti-ENO1 Rabbit pAb	WB/IHC-P
PTM-6021	Anti-ENO3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6731	Anti-ENO3 Rabbit pAb	WB/IHC-P
PTM-7124	Anti-ENPP1/PC1 Rabbit mAb	WB/IHC-P
PTM-5430	Anti-EPHX2 Mouse mAb	WB
PTM-6531	Anti-ERp44 Mouse mAb	WB/IHC-P
PTM-6068	Anti-ESD Mouse mAb	WB/ICC/IF
PTM-6816	Anti-Estrogen Related Receptor alpha Rabbit mAb	WB
PTM-6412	Anti-FABP1 Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-6168	Anti-FABP2 Mouse mAb	WB
PTM-6313	Anti-FABP3 Mouse mAb	WB/IHC-P
PTM-6416	Anti-FABP4 Rabbit mAb	WB/IHC-P
PTM-5865	Anti-FABP5 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5236	Anti-Factor XIIIa Rabbit mAb	WB/IHC-P
PTM-5444	Anti-FASN Rabbit mAb	WB/IHC-P/FC/IP
PTM-7439	Anti-Ferritin Heavy Chain Rabbit mAb	WB/IHC-P
PTM-7141	Anti-Ferritin Light Chain Rabbit mAb	WB
PTM-6746	Anti-Ferritin Light Chain Rabbit pAb	WB/IHC-P
PTM-5942	Anti-FIS1 Mouse mAb	WB/ICC/IF
PTM-5545	Anti-FOXD3 Mouse mAb	WB
PTM-6229	Anti-FSH beta Rabbit mAb	WB/IHC-P
PTM-6922	Anti-FTO Rabbit mAb	WB/IHC-P
PTM-5473	Anti-G6PD Mouse mAb	WB/ICC/IF
PTM-6694	Anti-GAD65 Rabbit pAb	WB/IHC-P
PTM-6684	Anti-GAPDH Mouse mAb	WB/IHC-P/ICC/IF/FC
PTM-6737	Anti-GART Mouse mAb	WB/ICC/IF
PTM-7059	Anti-GBA Rabbit mAb	WB/IHC-P
PTM-6959	Anti-GCLC Rabbit mAb	WB
PTM-6178	Anti-GCLM Mouse mAb	WB/IHC-P
PTM-6031	Anti-GLP-1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7112	Anti-Glucagon Rabbit mAb	IHC-P
PTM-5698	Anti-GLUD1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5632	Anti-GLUD1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5203	Anti-GLUT1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6209	Anti-GLUT2 Rabbit mAb	WB/FC
PTM-6992	Anti-GLUT3+GLUT14 Rabbit mAb	WB
PTM-6775	Anti-Glutathione Peroxidase 4 Rabbit mAb	WB
PTM-6297	Anti-Glycogen synthase 1/GYS1 RabbimAb	WB/IHC-P/ICC/IF
PTM-5649	Anti-GNB2 Rabbit mAb	WB
PTM-5959	Anti-GNB2 Rabbit mAb	WB/ICC/IF
PTM-5702	Anti-GNB3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6406	Anti-GNMT Mouse mAb	WB
PTM-6037	Anti-GOT2 Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
	Anti-GPI Rabbit pAb	WB/IHC-P
	Anti-GRP75 Rabbit mAb	WB/ICC/IF
PTM-6767	Anti-GRP94 Mouse mAb	WB
PTM-6180	Anti-GRP94 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5727	Anti-GSS Mouse mAb	WB/IHC-P
PTM-5800	Anti-GST3/GSTP1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5695	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6598	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5992	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5992	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P
PTM-6003	Anti-GST3/GSTP1 Rabbit pAb	WB/IHC-P/ICC/IF
PTM-7140	Anti-GSTK1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6691	Anti-hCG beta Rabbit mAb	IHC-P/FC
PTM-6391	Anti-Heme Oxygenase 1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6196	Anti-Hemopexin Rabbit mAb	WB/IHC-P
PTM-6478	Anti-Hexokinase I Rabbit pAb	WB/IHC-P
PTM-6341	Anti-HGD Rabbit mAb	IHC-P
PTM-6999	Anti-HIF-1 alpha Mouse mAb	WB/IHC-P
PTM-6286	Anti-HK2 Mouse mAb	WB
PTM-5851	Anti-HMGCR Rabbit mAb	WB/ICC/IF/FC/IP
PTM-5371	Anti-HNF-1B Rabbit mAb	WB
PTM-6018	Anti-HSD17B1 Rabbit mAb	WB/IHC-P
PTM-6320	Anti-HTR2C Rabbit mAb	WB
PTM-6672	Anti-Iba1 Rabbit mAb	IHC-P/mIHC-P
PTM-7213	Anti-Iba1/AIF-1 Rabbit mAb	WB/IHC-P
PTM-7108	Anti-IDH1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6749	Anti-IDH2 Rabbit pAb	WB/IHC-P
PTM-5949	Anti-IGFBP3 Rabbit mAb	WB
PTM-6878	Anti-IL-1 beta Mouse mAb	WB
PTM-6475	Anti-IL-1 beta Rabbit mAb	WB
PTM-5116	Anti-IL-18 Mouse mAb	WB
PTM-6799	Anti-IMPDH2 Rabbit mAb	WB/ICC/IF/IP
PTM-6235	Anti-iNOS Mouse mAb	WB/ICC/IF
PTM-6170	Anti-Insulin Rabbit mAb	IHC-P
PTM-6217		
PTM-5642		

Cat.#	Product	Recommended Applications
	Anti-IRG1 Rabbit pAb	WB
	Anti-IVD Mouse mAb	WB/IHC-P
PTM-6769	Anti-LACTB Mouse mAb	WB/IHC-P/ICC/IF
PTM-6484	Anti-LAMP1 Mouse mAb	WB/ICC/IF
PTM-5691	Anti-Langerin Rabbit pAb	WB
PTM-5775	Anti-LC3B Rabbit mAb	WB/IHC-P
PTM-6964	Anti-LDHA Rabbit pAb	WB/IHC-P
PTM-6384	Anti-LDHB Mouse mAb	WB/ICC/IF
PTM-6829	Anti-LDHC Rabbit mAb	WB
PTM-5869	Anti-LDL Receptor Mouse mAb	WB/ICC/IF
PTM-6302	Anti-LDL Receptor Rabbit mAb	WB/IHC-P/FC
PTM-5861	Anti-Legumain Mouse mAb	WB/IHC-P
PTM-6258	Anti-LONP1 Mouse mAb	WB/IHC-P
PTM-6596	Anti-LRP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5867	Anti-LTA4H Rabbit mAb	WB
PTM-6610	Anti-Luteinizing Hormone beta Rabbit mAb	WB/IHC-P
PTM-6024	Anti-Luteinizing Hormone beta Rabbit mAb	IHC-P
PTM-5985	Anti-MAP1LC3A Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5849	Anti-MAT1A Mouse mAb	WB/ICC/IF
PTM-5896	Anti-MAT2B Mouse mAb	WB/ICC/IF
PTM-6022	Anti-MDH1 Mouse mAb	WB
PTM-6450	Anti-MDH2 Rabbit pAb	WB/IHC-P
PTM-5385	Anti-ME2 Mouse mAb	WB/ICC/IF
PTM-7133	Anti-Met Mouse mAb	WB
PTM-6002	Anti-MFF Rabbit mAb	WB
PTM-7177	Anti-MFN1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7185	Anti-MGME1 Mouse mAb	WB
PTM-5475	Anti-Mitofusin 2 Mouse mAb	WB
PTM-5474	Anti-MKS1 Mouse mAb	WB/ICC/IF
PTM-5161	Anti-MRPL45 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5818	Anti-MTAP Rabbit mAb	WB/IHC-P
PTM-5866	Anti-MTCO1 Mouse mAb	WB
PTM-6545	Anti-MTCO2 Mouse mAb	WB/IHC-P
PTM-5109	Anti-MTHFD1 Mouse mAb	WB/ICC/IF
PTM-5167		
PTM-5694		

Cat.#	Product	Recommended Applications
	Anti-mTOR Mouse mAb	WB/IHC-P/ICC/IF
	Anti-mtTFA Rabbit mAb	WB/IHC-P
PTM-6594	Anti-Myoglobin Rabbit mAb	WB/IHC-P
PTM-6993	Anti-Na ⁺ /K ⁺ ATPase alpha-1 Mouse mAb	WB/IHC-P
PTM-7079	Anti-Na ⁺ /K ⁺ ATPase alpha-1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5533	Anti-NaPi2b Rabbit mAb	IHC-P
PTM-6271	Anti-NAPRT1 Mouse mAb	WB
PTM-5348	Anti-NAT1 Mouse mAb	WB
PTM-5774	Anti-NDUFA4L2 Mouse mAb	WB/ICC/IF
PTM-6738	Anti-NDUFB8 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5765	Anti-NDUFS3 Rabbit mAb	WB/IHC-P
PTM-5899	Anti-NFκB p65 Rabbit mAb	WB
PTM-7189	Anti-NRF1 Mouse mAb	WB/ICC/IF
PTM-5808	Anti-NSE Mouse mAb	WB
PTM-5511	Anti-OGDH Mouse mAb	WB/ICC/IF
PTM-6449	Anti-OLR1/LOX1 Rabbit mAb	WB/ICC/IF/FC
PTM-5974	Anti-OMA1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6365	Anti-OTC Mouse mAb	WB/IHC-P
PTM-6335	Anti-OXA1L Mouse mAb	WB/ICC/IF
PTM-6752	Anti-P4HA2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6176	Anti-P4HB Mouse mAb	WB/IHC-P/ICC/IF
PTM-6033	Anti-PADI2 Mouse mAb	WB
PTM-5314	Anti-PC Mouse mAb	WB/IHC-P
PTM-6201	Anti-PCCB Mouse mAb	WB/IHC-P
PTM-5699	Anti-PCK1 Mouse mAb	WB/IHC-P
PTM-5970	Anti-PCK2 Rabbit mAb	WB/ICC/IF
PTM-5939	Anti-PCSK9 Rabbit mAb	WB
PTM-5650	Anti-PDCD2 Rabbit mAb	WB
PTM-7067	Anti-PDHA1 Rabbit mAb	WB
PTM-6410	Anti-PDHE1 Mouse mAb	WB/ICC/IF
PTM-6541	Anti-PDK2 Rabbit mAb	WB
PTM-5386	Anti-PDP1 Rabbit mAb	WB/ICC/IF
PTM-6321	Anti-PDX1 Rabbit pAb	WB/IHC-P
PTM-5967	Anti-Pepsinogen II/PGC Rabbit pAb	WB/IHC-P
PTM-6907		
PTM-6723		

Cat.#	Product	Recommended Applications
PTM-7174	Anti-Perilipin-1 Rabbit mAb	WB
PTM-7027	Anti-Perilipin-2 Rabbit mAb	WB/ICC/IF
PTM-6275	Anti-PGAM1 Mouse mAb	WB/ICC/IF
PTM-5888	Anti-PGAM1 Rabbit mAb	WB/ICC/IF
PTM-6861	Anti-PGAM1 Rabbit pAb	WB/IHC-P
PTM-5975	Anti-PHGDH Mouse mAb	WB/ICC/IF
PTM-6615	Anti-Phospho-(Ser/Thr) PDK1 Docking Motif Mouse mAb	WB
PTM-6515	Anti-Phospho-4E-BP1 (Thr37/46) Rabbit mAb	WB
PTM-6785	Anti-Phospho-AMPK beta 1 (Ser182) Rabbit pAb	WB
PTM-7083	Anti-Phospho-Bcl-2 (Thr56) Rabbit pAb	WB
PTM-6887	Anti-Phospho-CAD (Ser1859) Rabbit pAb	WB
PTM-6855	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit mAb	WB
PTM-6801	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit pAb	WB/ICC/IF
PTM-7178	Anti-phospholamban Mouse mAb	WB/IHC-P
PTM-6639	Anti-PI3 Kinase Class III Rabbit mAb	WB/ICC/IF
PTM-5670	Anti-PI3 Kinase p85 alpha Rabbit mAb	IHC-P/ICC/IF
PTM-7162	Anti-PKA C-alpha Mouse mAb	WB/IHC-P
PTM-7161	Anti-PKC mu Rabbit mAb	WB
PTM-5087	Anti-PKM2 Mouse mAb	WB/ICC/IF/FC
PTM-5820	Anti-PLCB3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7180	Anti-PLCG1 Rabbit mAb	WB
PTM-6859	Anti-PLCG2 Rabbit mAb	WB
PTM-6262	Anti-PPP1CA Mouse mAb	WB/ICC/IF
PTM-6128	Anti-PPP1CB Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6486	Anti-PPP1CB Rabbit mAb	WB/IHC-P
PTM-5934	Anti-PRDX1 Rabbit mAb	WB/IHC-P/FC
PTM-6525	Anti-PRDX3 Mouse mAb	WB/ICC/IF
PTM-6522	Anti-PRDX6 Mouse mAb	WB
PTM-6120	Anti-Prohibitin Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7098	Anti-PSAP Rabbit pAb	WB
PTM-5757	Anti-PSMA/FOLH1 Rabbit mAb	WB/IHC-P/FC
PTM-6207	Anti-PYCR1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7045	Anti-PYK2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6240	Anti-RBP4 Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6246	Anti-REDD1 Mouse mAb	WB/ICC/IF
PTM-5043	Anti-RPE65 Rabbit pAb	WB/IHC-P
PTM-5594	Anti-Rubisco Mouse mAb	WB
PTM-7182	Anti-Rubisco Mouse mAb	WB
PTM-7025	Anti-RXRA Rabbit mAb	WB/ICC/IF
PTM-7164	Anti-Ryanodine Receptor Rabbit mAb	IHC-P
PTM-6988	Anti-Scavenging Receptor SR-BI Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6553	Anti-SCLY Mouse mAb	WB
PTM-5162	Anti-SDHA Mouse mAb	WB/IHC-P
PTM-5164	Anti-SDHB Mouse mAb	WB/IHC-P
PTM-6973	Anti-SLC16A1 Rabbit mAb	WB/IHC-P
PTM-6969	Anti-SLC25A12 Rabbit mAb	WB/IHC-P
PTM-5465	Anti-SLC9A9 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6127	Anti-Smac/Diablo Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5977	Anti-Smac/Diablo Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5681	Anti-SOD1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5703	Anti-SOD2 Mouse mAb	WB/IHC-P
PTM-6080	Anti-Sorbitol Dehydrogenase Mouse mAb	WB
PTM-7053	Anti-SP1 Rabbit mAb	WB/IHC-P
PTM-6434	Anti-SQSTM1/p62 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5483	Anti-SQSTM1/p62 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5508	Anti-SRF Mouse mAb	WB
PTM-7211	Anti-StAR Rabbit mAb	WB/IHC-P
PTM-5651	Anti-STXBP1 Mouse mAb	WB/IHC-P
PTM-3059	Anti-SUCLG1 Rabbit mAb	WB/IHC-P
PTM-3560	Anti-SUCLG2 Rabbit mAb	WB/IHC-P
PTM-7184	Anti-Synapsin-1 Mouse mAb	WB/IHC-P
PTM-6742	Anti-Syntaxin 1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7214	Anti-Syntaxin 3 Rabbit mAb	WB/IHC-P
PTM-6946	Anti-TBLR1/TBL1XR1 Rabbit mAb	WB/IHC-P
PTM-5661	Anti-TFEB Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6897	Anti-TGF beta 1 Rabbit mAb	WB/IHC-P
PTM-7029	Anti-TGF beta 2 Mouse mAb	WB/IHC-P
PTM-5265	Anti-Thymidylate Synthase Rabbit mAb	WB/IHC-P/ICC/IF/IP

Cat.#	Product	Recommended Applications
PTM-6589	Anti-Thymidylate Synthase Rabbit mAb	WB/ICC/IF/FC
PTM-7015	Anti-TOMM22/TOM22 Rabbit mAb	WB/IHC-P
PTM-5693	Anti-TOMM40 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5901	Anti-TPI1 Mouse mAb	WB/ICC/IF
PTM-6165	Anti-TPPII Mouse mAb	WB/IHC-P
PTM-5617	Anti-TPSAB1 Rabbit mAb	IHC-P
PTM-6081	Anti-Transferrin Mouse mAb	WB/ICC/IF
PTM-6311	Anti-Transketolase Mouse mAb	WB/IHC-P/FC
PTM-6325	Anti-TrkA (phospho Y791) Rabbit mAb	WB
PTM-5583	Anti-TSHB Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6009	Anti-TSPAN7 Mouse mAb	ICC/IF/FC
PTM-5689	Anti-TXN Mouse mAb	WB/IHC-P/ICC/IF
PTM-6671	Anti-TXNRD1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5416	Anti-TYMS Mouse mAb	WB/IHC-P
PTM-5652	Anti-UBC9 Mouse mAb	WB
PTM-6053	Anti-UCP1 Rabbit mAb	WB/IHC-P
PTM-5837	Anti-UGDH Mouse mAb	WB/IHC-P/ICC/IF
PTM-5008	Anti-UQCRC2 Mouse mAb	WB/IHC-P
PTM-6847	Anti-UQCRFS1/RISP Rabbit mAb	WB/IHC-P/IP
PTM-6998	Anti-USF1 Rabbit mAb	WB
PTM-6527	Anti-VDAC1/Porin Mouse mAb	WB/IHC-P
PTM-6157	Anti-VDAC1/Porin Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6935	Anti-VDAC1/Porin Rabbit pAb	WB/IHC-P
PTM-5894	Anti-WNT10B Mouse mAb	WB/IHC-P/ICC/IF
PTM-5477	Anti-XIAP Mouse mAb	WB/ICC/IF
PTM-5954	Anti-YME1L1 Mouse mAb	WB/IHC-P/ICC/IF

Developmental & Stem Cell Research Antibodies

PTM-7063	Anti-14-3-3 eta/YWHAH Rabbit pAb	WB/IHC-P
PTM-7040	Anti-14-3-3 gamma Rabbit mAb	WB/IHC-P
PTM-6307	Anti-14-3-3 theta Rabbit mAb	WB
PTM-6392	Anti-5T4 Rabbit mAb	WB/IHC-P/IP

Cat.#	Product	Recommended Applications
PTM-6848	Anti-Adiponectin Rabbit mAb	WB/IHC-P
PTM-6243	Anti-AFP Mouse mAb	WB/ICC/IF
PTM-6397	Anti-AFP Rabbit mAb	WB/IHC-P/FC
PTM-6322	Anti-AFP Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6853	Anti-AGR2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7335	Anti-AID Rabbit mAb	WB/IHC-P
PTM-6758	Anti-AIDA Rabbit pAb	WB/IHC-P
PTM-5794	Anti-Albumin Mouse mAb	WB
PTM-6588	Anti-beta Catenin Rabbit pAb	WB/IHC-P
PTM-5709	Anti-BMI1 Mouse mAb	WB
PTM-5466	Anti-brachyury Mouse mAb	WB
PTM-6677	Anti-c-Kit Rabbit mAb	IHC-P
PTM-5675	Anti-c-Kit Rabbit mAb	IHC-P
PTM-7317	Anti-c-Myc Rabbit mAb	WB/ChIP/IP
PTM-6629	Anti-Calponin 1 Mouse mAb	WB
PTM-6661	Anti-Casein Kinase 1 gamma 2 Rabbit mAb	WB/ICC/IF/FC
PTM-6614	Anti-Casein Kinase 2 beta Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6507	Anti-Casein Kinase 2 beta Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6938	Anti-Caveolin-3 Rabbit mAb	WB/IHC-P
PTM-5058	Anti-CDX2 Rabbit mAb	IHC-P/FC
PTM-6743	Anti-Claudin-6 Mouse mAb	WB
PTM-7242	Anti-COL4A1 Rabbit mAb	WB
PTM-6219	Anti-Collagen I Rabbit mAb	WB
PTM-7147	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-5771	Anti-CXCR4 Mouse mAb	FC
PTM-6444	Anti-CYR61 Mouse mAb	WB/IHC-P
PTM-6795	Anti-DDX3 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6370	Anti-DDX4 Rabbit mAb	IHC-P
PTM-6904	Anti-DDX4 Rabbit pAb	WB/IHC-P
PTM-5770	Anti-Decorin Mouse mAb	WB/IHC-P
PTM-6353	Anti-DKK3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7415	Anti-Dlx5 Mouse mAb	WB
PTM-7057	Anti-DPP4/CD26 Rabbit mAb	WB/IHC-P
PTM-6381	Anti-DVL2 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5060	Anti-E-Cadherin Rabbit mAb	IHC-P
PTM-7374	Anti-EN1/Engrailed 1 Rabbit mAb	WB
PTM-6304	Anti-ESRRB Mouse mAb	WB/ICC/IF
PTM-6168	Anti-FABP2 Mouse mAb	WB
PTM-6313	Anti-FABP3 Mouse mAb	WB/IHC-P
PTM-6416	Anti-FABP4 Rabbit mAb	WB/IHC-P
PTM-6408	Anti-FGFR3 Rabbit mAb	WB/IP
PTM-6790	Anti-Fibronectin Rabbit pAb	WB/IHC-P
PTM-5324	Anti-FOXA1 Rabbit mAb	IHC-P/ICC/IF
PTM-7295	Anti-FOXA2 Rabbit mAb	WB/ChIP
PTM-7002	Anti-FOXC1 Rabbit mAb	WB
PTM-6354	Anti-FRZB Mouse mAb	WB/ICC/IF
PTM-6994	Anti-GATA1 Rabbit mAb	WB/IHC-P
PTM-7444	Anti-GATA1 Rabbit mAb	WB/IP/ChIP
PTM-6465	Anti-GATA2 Mouse mAb	WB
PTM-6077	Anti-GATA3 Mouse mAb	WB/ICC/IF
PTM-7112	Anti-Glucagon Rabbit mAb	IHC-P
PTM-5669	Anti-Glypican 3 Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6751	Anti-Glypican-3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5514	Anti-GNL3 Mouse mAb	WB
PTM-5916	Anti-GPI Mouse mAb	WB/ICC/IF
PTM-7240	Anti-GPR64 Rabbit mAb	WB
PTM-7313	Anti-GSK3 beta Mouse mAb	WB/IHC-P
PTM-6720	Anti-HES5 Rabbit mAb	WB
PTM-6814	Anti-HOXB4 Rabbit mAb	WB
PTM-6415	Anti-ICAM-2 Rabbit mAb	WB/IHC-P/FC
PTM-6698	Anti-IFITM1 Rabbit pAb	WB
PTM-7150	Anti-IGF1 Receptor Rabbit mAb	WB/FC
PTM-5223	Anti-IKAROS Rabbit mAb	WB/IHC-P
PTM-6359	Anti-IL-7R/CD127 Rabbit mAb	ICC/IF/FC
PTM-6340	Anti-Integrin beta 1 Rabbit mAb	WB
PTM-6237	Anti-IRF6 Rabbit mAb	WB/ICC/IF
PTM-5672	Anti-JunB Rabbit mAb	WB
PTM-5641	Anti-JunB Rabbit mAb	IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5238	Anti-KLF4 Rabbit mAb	IHC-P/ICC/IF
PTM-5057	Anti-LMO2 Rabbit mAb	IHC-P
PTM-6726	Anti-LRP6 Rabbit mAb	WB
PTM-5803	Anti-Ly-6A/E Rabbit mAb	WB/IHC-P
PTM-7341	Anti-MEF2C Rabbit mAb	WB/IHC-P
PTM-6066	Anti-MESDC2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5474	Anti-MGME1 Mouse mAb	WB
PTM-6010	Anti-Nanog Mouse mAb	WB
PTM-6670	Anti-NCAM1/CD56 Mouse mAb	FC
PTM-6218	Anti-NCAM1/CD56 Rabbit mAb	WB/FC/IP
PTM-7324	Anti-Notch1 Rabbit mAb	WB/IHC-P
PTM-7225	Anti-Notch4 Rabbit mAb	WB
PTM-6449	Anti-NSE Mouse mAb	WB
PTM-7013	Anti-NUMB Rabbit pAb	WB
PTM-5593	Anti-OCT4 Mouse mAb	WB
PTM-5721	Anti-p75NTR Mouse mAb	ICC/IF/FC
PTM-6395	Anti-p75NTR Rabbit mAb	WB/IHC-P
PTM-5338	Anti-PAX8 Rabbit mAb	IHC-P
PTM-5872	Anti-PAX8 Rabbit mAb	IHC-P/ICC/IF
PTM-7445	Anti-PBX1 Rabbit mAb	WB/IHC-P/IP/ChIP
PTM-7067	Anti-PCSK9 Rabbit mAb	WB
PTM-6907	Anti-PDX1 Rabbit pAb	WB/IHC-P
PTM-7066	Anti-PEG10 Rabbit mAb	WB
PTM-7174	Anti-Perilipin-1 Rabbit mAb	WB
PTM-7423	Anti-Phospho-beta Catenin (Thr41/Ser45) Rabbit mAb	WB
PTM-6921	Anti-Phospho-c-Myc (Ser62) Rabbit mAb	WB
PTM-7429	Anti-Phospho-Cardiac Troponin I (Ser23/24) Rabbit mAb	WB
PTM-7100	Anti-Phospho-CD19 (Tyr531) Rabbit pAb	WB
PTM-6931	Anti-Phospho-CD31 (Tyr713) Rabbit mAb	WB
PTM-6955	Anti-Phospho-CDCP1 (Tyr743) Rabbit pAb	WB
PTM-7024	Anti-Phospho-Stat3 (Ser727) Rabbit mAb	WB
PTM-7041	Anti-Phospho-Vimentin (Ser39) Rabbit mAb	WB
PTM-7432	Anti-Phospho-Vimentin (Ser39) Rabbit mAb	WB/ICC/IF
PTM-7293	Anti-Pim-1 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6173	Anti-PITX2/RGS Rabbit mAb	IHC-P/FC
PTM-5612	Anti-Placental Lactogen Rabbit mAb	WB/IHC-P
PTM-6270	Anti-Prion protein PrP Rabbit mAb	WB/IHC-P
PTM-5509	Anti-RAB3GAP1 Mouse mAb	WB/ICC/IF
PTM-6163	Anti-RBM15 Mouse mAb	WB/IHC-P
PTM-6888	Anti-ROR2 Mouse mAb	WB
PTM-7303	Anti-RUNX1/AML1 Rabbit mAb	WB
PTM-6840	Anti-RUNX1/AML1 Rabbit pAb	WB
PTM-6617	Anti-RUNX2 Rabbit mAb	WB
PTM-7096	Anti-S100B Rabbit pAb	WB/IHC-P
PTM-5323	Anti-SALL4 Rabbit mAb	IHC-P
PTM-5799	Anti-SF1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7031	Anti-SGK1 Rabbit mAb	WB
PTM-6929	Anti-Smad2 (phospho S255) Rabbit mAb	WB
PTM-5153	Anti-Smad2 Mouse mAb	WB
PTM-5054	Anti-Smad4 Rabbit mAb	WB/IHC-P
PTM-5449	Anti-Smad4 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7288	Anti-SOCS3 Rabbit mAb	WB
PTM-5076	Anti-SOX10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7065	Anti-SOX17 Rabbit mAb	WB/IHC-P
PTM-6390	Anti-SOX2 Rabbit mAb	WB/IHC-P
PTM-5072	Anti-SOX2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7053	Anti-SP1 Rabbit mAb	WB/IHC-P
PTM-6736	Anti-SPARC Rat mAb	WB
PTM-5408	Anti-STK3 Mouse mAb	WB
PTM-6143	Anti-TAZ Mouse mAb	WB
PTM-6288	Anti-TCF-3/TCF7L1 Mouse mAb	WB
PTM-6451	Anti-TCF3/E2A Mouse mAb	WB
PTM-6518	Anti-TCF3/E2A Mouse mAb	WB/IHC-P
PTM-7280	Anti-TCF4 Rabbit mAb	WB
PTM-6897	Anti-TGF beta 1 Rabbit mAb	WB/IHC-P
PTM-7029	Anti-TGF beta 2 Mouse mAb	WB/IHC-P
PTM-5400	Anti-TH Mouse mAb	WB
PTM-6295	Anti-Thrombomodulin Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6759	Anti-TIE2 Rabbit pAb	WB/IHC-P
PTM-6978	Anti-TTF1 Rabbit mAb	WB/IHC-P
PTM-5105	Anti-TTF1 Rabbit mAb	IHC-P
PTM-7105	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-6623	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-7086	Anti-Von Willebrand Factor Rabbit mAb	WB/IHC-P
PTM-7298	Anti-WT1 Rabbit mAb	WB/IHC-P
PTM-7035	Anti-WT1 Rabbit pAb	WB

Neuroscience Research Antibodies

PTM-7040	Anti-14-3-3 gamma Rabbit mAb	WB/IHC-P
PTM-5629	Anti-14-3-3 sigma Mouse mAb	WB/IHC-P
PTM-5462	Anti-14-3-3 theta Mouse mAb	WB/ICC/IF
PTM-6942	Anti-AACT Rabbit mAb	WB/IHC-P
PTM-6166	Anti-ABAT Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6794	Anti-Acetyl-SOD2 (Lys68) Rabbit mAb	WB/IHC-P
PTM-7007	Anti-ACOT7 Rabbit pAb	WB/IHC-P
PTM-7056	Anti-ACTH Rabbit mAb	IHC-P
PTM-6848	Anti-Adiponectin Rabbit mAb	WB/IHC-P
PTM-7208	Anti-ADRA1A Rabbit mAb	WB/FC/ICC/IF
PTM-6853	Anti-AGR2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7253	Anti-AIF Mouse mAb	WB/IHC-P
PTM-6148	Anti-AKR7A2 Mouse mAb	WB/ICC/IF
PTM-6735	Anti-AKT (phospho T308) Rabbit mAb	WB
PTM-6643	Anti-AKT1 + AKT2 + AKT3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6476	Anti-AKT1 + AKT2 + AKT3 Rabbit mAb	WB/IHC-P/FC
PTM-7179	Anti-ALDH1A1 Rabbit mAb	WB/IHC-P
PTM-5913	Anti-ALDH1L1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6091	Anti-ALDH5A1/SSADH Mouse mAb	WB/IHC-P
PTM-6666	Anti-ALDH5A1/SSADH Rabbit mAb	WB
PTM-6811	Anti-ALDH9A1 Mouse mAb	WB/IHC-P
PTM-6729	Anti-Aldolase A Rabbit pAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6078	Anti-Aldolase C Mouse mAb	WB/IHC-P/ICC/IF
PTM-5026	Anti-alpha Skeletal Muscle Actin Rabbit mAb	WB/IP
PTM-5886	Anti-alpha SNAP Mouse mAb	WB/ICC/IF
PTM-7176	Anti-alpha tubulin(acetyl K40) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6443	Anti-alpha-Internexin Mouse mAb	WB
PTM-6827	Anti-alpha-Internexin Rabbit pAb	WB/IHC-P/FC
PTM-5439	Anti-alpha-Synuclein Rabbit mAb	WB
PTM-6836	Anti-alpha-Synuclein Rabbit mAb	WB/IHC-P
PTM-6345	Anti-alpha-Synuclein Rabbit mAb	WB/ICC/IF/FC
PTM-7250	Anti-AMPA Receptor 1 Rabbit mAb	WB/IHC-P
PTM-6739	Anti-AMPK alpha 1+AMPK alpha 2 Rabbit mAb	WB
PTM-7149	Anti-APAF1 Rabbit mAb	WB/IHC-P
PTM-6082	Anti-ApoA1 Mouse mAb	WB
PTM-5915	Anti-ApoE Mouse mAb	WB/IHC-P/ICC/IF
PTM-7230	Anti-AQP4 Mouse mAb	WB
PTM-6425	Anti-Arc Mouse mAb	WB/IHC-P/ICC/IF
PTM-5730	Anti-Arginase-1 Rabbit mAb	WB/IHC-P
PTM-7207	Anti-Ataxin-1 Rabbit mAb	WB
PTM-6195	Anti-Ataxin-3 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-7043	Anti-ATG16L1 Mouse mAb	WB
PTM-6514	Anti-ATG9A Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6881	Anti-Axl Rabbit mAb	WB
PTM-7344	Anti-BAF53B Rabbit mAb	WB
PTM-5041	Anti-BCL6 Rabbit mAb	WB/IHC-P
PTM-6793	Anti-Bcr Rabbit mAb	WB/IHC-P
PTM-6485	Anti-BCS1L Mouse mAb	WB
PTM-6588	Anti-beta Catenin Rabbit pAb	WB/IHC-P
PTM-6468	Anti-beta Dystroglycan Mouse mAb	WB/IHC-P/ICC/IF/FC
PTM-6976	Anti-beta III Tubulin Rabbit mAb	WB/IHC-P
PTM-6844	Anti-beta-Amyloid (1-40) Rabbit pAb	WB/IHC-P
PTM-6445	Anti-beta-Arrestin 1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7145	Anti-Bim Rabbit mAb	WB
PTM-5843	Anti-BMAL1 Rabbit mAb	WB/IHC-P/FC
PTM-5826	Anti-BRG1 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-6362	Anti-BRN3A Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6653	Anti-c-Fos Mouse mAb	WB
PTM-6524	Anti-C1q Rabbit mAb	WB
PTM-5140	Anti-Calbindin Mouse mAb	WB/IHC-P
PTM-6417	Anti-Calbindin Rabbit mAb	WB/IHC-P
PTM-6950	Anti-Calnexin Rabbit mAb	WB/IHC-P
PTM-5759	Anti-Calpastatin Rabbit mAb	IHC-P/ICC/IF/FC
PTM-6806	Anti-Calretinin Rabbit pAb	WB/IHC-P
PTM-5505	Anti-Calsequestrin 2 Mouse mAb	WB/IHC-P
PTM-6045	Anti-CAMKII alpha Mouse mAb	WB/ICC/IF
PTM-6299	Anti-CAMKII alpha Rabbit mAb	WB/IHC-P
PTM-6233	Anti-Caspase-6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5793	Anti-Caspase-6 Rabbit mAb	WB/ICC/IF
PTM-6781	Anti-Caspase-7 Rabbit mAb	WB/IHC-P
PTM-6782	Anti-Caspase-7 Rabbit pAb	WB/IHC-P
PTM-6300	Anti-Caspr Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6372	Anti-Cathepsin B Rabbit mAb	WB
PTM-6369	Anti-Cathepsin D Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7436	Anti-Cathepsin K Mouse mAb	WB/IHC-P
PTM-6938	Anti-Caveolin-3 Rabbit mAb	WB/IHC-P
PTM-6521	Anti-CDK1/Cdc2 Mouse mAb	WB/ICC/IF
PTM-7072	Anti-CDK1/Cdc2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5096	Anti-CDK5 Mouse mAb	WB/IHC-P
PTM-6791	Anti-CDR2L Mouse mAb	WB/ICC/IF
PTM-5848	Anti-CGA/hCG alpha Rabbit mAb	WB/IHC-P
PTM-7082	Anti-CGA/hCG alpha Rabbit mAb	WB/IHC-P
PTM-6765	Anti-CGA/hCG alpha Rabbit pAb	WB/IHC-P
PTM-6073	Anti-CHAT Mouse mAb	WB
PTM-6239	Anti-CHMP2B Rabbit mAb	WB/IHC-P
PTM-5148	Anti-Chromogranin A Mouse mAb	WB
PTM-5806	Anti-Chromogranin A Rabbit mAb	WB/IHC-P/FC
PTM-7204	Anti-Chromogranin A Rabbit mAb	IHC-P
PTM-7131	Anti-Claudin-11 Rabbit mAb	WB
PTM-7239	Anti-Claudin-11 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6280	Anti-CIC-2 Rabbit mAb	WB/IHC-P
PTM-7246	Anti-Cleaved caspase-3 Rabbit mAb	WB
PTM-6873	Anti-Clusterin beta Chain Rabbit mAb	WB
PTM-7261	Anti-Clusterin Rabbit mAb	WB/IHC-P
PTM-5435	Anti-CNPase Mouse mAb	WB/ICC/IF
PTM-5729	Anti-Cofilin-1 Rabbit mAb	WB/ICC/IF/FC
PTM-7319	Anti-COMT Rabbit mAb	WB
PTM-5870	Anti-COPS3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6951	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-7147	Anti-COX2/Cyclooxygenase 2 Rabbit mAb	WB/IHC-P
PTM-5682	Anti-CRALBP Rabbit mAb	WB/IHC-P
PTM-7058	Anti-Creatine kinase B type Rabbit mAb	WB/IHC-P
PTM-5595	Anti-CREB Mouse mAb	WB/IHC-P
PTM-6298	Anti-CRMP2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6156	Anti-Ctip2 Rabbit mAb	WB/IHC-P/FC
PTM-6333	Anti-CX3CR1 Rabbit mAb	WB
PTM-5771	Anti-CXCR4 Mouse mAb	FC
PTM-6901	Anti-CXCR5 Rabbit pAb	WB
PTM-7181	Anti-Cytoplasmic Dynein Intermediate chain Mouse mAb	WB
PTM-7104	Anti-DARPP-32 Rabbit mAb	WB/IHC-P
PTM-7165	Anti-DFNA5/GSDME Mouse mAb	WB
PTM-7166	Anti-DFNA5/GSDME Mouse mAb	WB
PTM-6171	Anti-DHCR24 Mouse mAb	WB/ICC/IF
PTM-6708	Anti-DHFR Mouse mAb	WB
PTM-6150	Anti-DIO2 Mouse mAb	WB/IHC-P
PTM-6757	Anti-DJ-1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5589	Anti-DLL3 Rabbit mAb	WB
PTM-7415	Anti-Dlx5 Mouse mAb	WB
PTM-7262	Anti-DNA Polymerase beta Rabbit mAb	WB/IHC-P
PTM-7195	Anti-DNAJC6 Rabbit mAb	WB/IHC-P
PTM-6256	Anti-DOPA Decarboxylase/DDC Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7109	Anti-Dopamine beta Hydroxylase Rabbit mAb	IHC-P
PTM-6140	Anti-Dynamin 1 Mouse mAb	WB/IHC-P
PTM-6161	Anti-Dynamin 1 Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5602	Anti-DYNLL1 Rabbit mAb	WB/IHC-P
PTM-6947	Anti-EAAT1 Rabbit mAb	WB
PTM-7281	Anti-EAAT3 Rabbit mAb	WB/IHC-P
PTM-7217	Anti-EEA1 Mouse mAb	WB/ICC/IF
PTM-7148	Anti-EEA1 Rabbit mAb	WB
PTM-7199	Anti-EGR2 Rabbit mAb	WB
PTM-5804	Anti-EIF2S1 Rabbit mAb	WB/IHC-P
PTM-7374	Anti-EN1/Engrailed 1 Rabbit mAb	WB
PTM-6087	Anti-EPHB4 Rabbit mAb	WB/ICC/IF/FC
PTM-6667	Anti-ErbB3 / HER3 Rabbit mAb	WB
PTM-5904	Anti-ERC1 Mouse mAb	WB/ICC/IF
PTM-6516	Anti-ERK5 Rabbit mAb	WB
PTM-5480	Anti-ERO1L Mouse mAb	WB/ICC/IF
PTM-6967	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6838	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6792	Anti-Estrogen Receptor alpha Rabbit mAb	WB/IHC-P
PTM-6816	Anti-Estrogen Related Receptor alpha Rabbit mAb	WB
PTM-5732	Anti-F4/80 Rabbit mAb	IHC-P
PTM-6703	Anti-FABP7 Rabbit mAb	WB/IHC-P
PTM-7439	Anti-Ferritin Heavy Chain Rabbit mAb	WB/IHC-P
PTM-7141	Anti-Ferritin Light Chain Rabbit mAb	WB
PTM-7266	Anti-ferritin light chain Rabbit mAb	WB
PTM-6746	Anti-Ferritin Light Chain Rabbit pAb	WB/IHC-P
PTM-5387	Anti-FKBP52 Mouse mAb	WB/IHC-P
PTM-5582	Anti-FLI1 Rabbit mAb	WB
PTM-5628	Anti-FOX2 Mouse mAb	WB
PTM-5545	Anti-FOXD3 Mouse mAb	WB
PTM-5688	Anti-FYN Mouse mAb	WB
PTM-6590	Anti-GABAA Receptor alpha 1 Rabbit mAb	WB/IHC-P
PTM-6981	Anti-GABA B Receptor 1 Rabbit mAb	WB/IHC-P
PTM-6694	Anti-GAD65 Rabbit pAb	WB/IHC-P
PTM-6356	Anti-Galectin-3 Mouse mAb	WB
PTM-6818	Anti-Galectin-9 Rabbit pAb	WB/IHC-P
PTM-6124	Anti-GAP43 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6784	Anti-GAP43 Rabbit pAb	WB/IHC-P
PTM-7059	Anti-GBA Rabbit mAb	WB/IHC-P
PTM-6178	Anti-GCLM Mouse mAb	WB/IHC-P
PTM-5957	Anti-GDI2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6305	Anti-Gemin 2 Rabbit mAb	WB
PTM-6104	Anti-GFAP (NT) Rabbit mAb	WB/IHC-P/mIHC-P/ICC/IF
PTM-5044	Anti-GFAP Mouse mAb	WB/IHC-P
PTM-7103	Anti-GIRK2 Rabbit mAb	WB
PTM-3024	Anti-GLO1 Rabbit pAb	WB
PTM-3561	Anti-GLO2 Rabbit pAb	WB
PTM-6031	Anti-GLP-1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7221	Anti-GluA4 Rabbit mAb	WB
PTM-7446	Anti-Glucocorticoid receptor Rabbit mAb	WB/ICC/IF
PTM-5698	Anti-GLUD1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7005	Anti-Glutamate receptor 3 Rabbit mAb	WB
PTM-6920	Anti-Glypican 1 Rabbit mAb	WB
PTM-5543	Anti-GNA13 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5860	Anti-GNAI2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5702	Anti-GNB3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5514	Anti-GNL3 Mouse mAb	WB
PTM-7251	Anti-Golgin-97 Rabbit mAb	WB
PTM-6110	Anti-GOSR1/GS28 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6767	Anti-GPI Rabbit pAb	WB/IHC-P
PTM-7326	Anti-Granulin Mouse mAb	WB/IHC-P
PTM-5965	Anti-GRK2 Rabbit mAb	WB/ICC/IF
PTM-5515	Anti-GRP75 Mouse mAb	WB/IHC-P
PTM-6180	Anti-GRP75 Rabbit mAb	WB/ICC/IF
PTM-7313	Anti-GSK3 beta Mouse mAb	WB/IHC-P
PTM-5695	Anti-GSS Mouse mAb	WB/IHC-P
PTM-5992	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5992	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6003	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7140	Anti-GST3/GSTP1 Rabbit mAb	WB/IHC-P
PTM-6691	Anti-GST3/GSTP1 Rabbit pAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6478	Anti-Heme Oxygenase 1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5258	Anti-Hemoglobin Alpha Rabbit mAb	IHC-P
PTM-6789	Anti-HER4/ErbB4 Rabbit mAb	IHC-P/ICC/IF
PTM-6720	Anti-HES5 Rabbit mAb	WB
PTM-6999	Anti-Hexokinase I Rabbit pAb	WB/IHC-P
PTM-6306	Anti-HIP1 Mouse mAb	WB/ICC/IF
PTM-6259	Anti-HIP2/LIG Rabbit mAb	WB
PTM-7137	Anti-HMGB1 Rat mAb	WB
PTM-7092	Anti-hnRNP A1 Rabbit pAb	WB
PTM-5737	Anti-Homer1 Rabbit mAb	WB/IHC-P
PTM-7322	Anti-HP1 beta Rabbit mAb	WB
PTM-7134	Anti-Hsp70 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7213	Anti-HTR2C Rabbit mAb	WB
PTM-6100	Anti-Huntingtin Mouse mAb	WB/ICC/IF
PTM-7108	Anti-Iba1 Rabbit mAb	IHC-P/mlIHC-P
PTM-6749	Anti-Iba1/AIF-1 Rabbit mAb	WB/IHC-P
PTM-7150	Anti-IGF1 Receptor Rabbit mAb	WB/FC
PTM-6717	Anti-IGF2R/CI-M6PR Rabbit mAb	WB/IHC-P
PTM-6799	Anti-IL-1 beta Rabbit mAb	WB
PTM-6869	Anti-IL-22 Rabbit pAb	WB
PTM-6712	Anti-IL-4 Rabbit mAb	WB/ICC/IF/FC
PTM-6217	Anti-iNOS Mouse mAb	WB/ICC/IF
PTM-6340	Anti-Integrin beta 1 Rabbit mAb	WB
PTM-6696	Anti-IRAK1 Rabbit mAb	WB/IHC-P
PTM-7135	Anti-JAK2 Rabbit mAb	WB
PTM-6948	Anti-JNK2 Rabbit mAb	WB/IHC-P
PTM-5672	Anti-JunB Rabbit mAb	WB
PTM-5641	Anti-JunB Rabbit mAb	IHC-P/ICC/IF
PTM-6899	Anti-Kappa Opioid Receptor Rabbit pAb	WB
PTM-7304	Anti-KEAP1 Rabbit mAb	WB
PTM-7227	Anti-KLC1 Rabbit mAb	WB
PTM-5604	Anti-Laminin Receptor Rabbit mAb	WB/ICC/IF
PTM-5775	Anti-LAMP1 Mouse mAb	WB/ICC/IF
PTM-7342	Anti-LAMP1 Rabbit mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-7382	Anti-LAP3 Mouse mAb	WB/IHC-P
PTM-6384	Anti-LC3B Rabbit mAb	WB/IHC-P
PTM-6949	Anti-Lck Rabbit mAb	WB/IHC-P
PTM-7359	Anti-LIF Rat mAb	WB/IHC-P
PTM-7340	Anti-LIMK1 Rabbit mAb	WB
PTM-5603	Anti-LOX Rabbit mAb	WB/ICC/IF
PTM-6023	Anti-LPCAT1 Mouse mAb	WB
PTM-5853	Anti-LPCAT1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6610	Anti-LRP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5856	Anti-LRPPRC Mouse mAb	WB/IHC-P/ICC/IF
PTM-5985	Anti-Luteinizing Hormone beta Rabbit mAb	WB/IHC-P
PTM-5849	Anti-Luteinizing Hormone beta Rabbit mAb	IHC-P
PTM-6162	Anti-MAG Rabbit mAb	WB/IHC-P
PTM-6856	Anti-MAO-A Rabbit mAb	WB/IHC-P
PTM-6977	Anti-MAPKAPK-2 Rabbit mAb	WB
PTM-7175	Anti-MASH1 Rabbit mAb	WB
PTM-7254	Anti-MED26 Mouse mAb	WB
PTM-7222	Anti-MED7 Rabbit mAb	WB/IHC-P
PTM-7341	Anti-MEF2C Rabbit mAb	WB/IHC-P
PTM-6890	Anti-MEK5 Rabbit mAb	WB
PTM-7333	Anti-Mer Mouse mAb	WB
PTM-7377	Anti-Mer Rabbit mAb	WB/IHC-P
PTM-7177	Anti-Met Mouse mAb	WB
PTM-7004	Anti-Metabotropic Glutamate Receptor 2/MGLUR2 Rabbit pAb	WB/IHC-P
PTM-5475	Anti-MFN1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7151	Anti-Mitofusin2+Mitofusin1 Mouse mAb	WB/IHC-P
PTM-7144	Anti-MMP13 Mouse mAb	WB
PTM-7375	Anti-MMP7 Rabbit mAb	WB/IHC-P
PTM-7259	Anti-MMP9 Rabbit mAb	WB/IHC-P
PTM-7139	Anti-MOG Rabbit pAb	WB/IHC-P
PTM-6800	Anti-Myelin Basic Protein Rabbit pAb	WB/IHC-P
PTM-7079	Anti-Myoglobin Rabbit mAb	WB/IHC-P
PTM-5533	Anti-Na ⁺ /K ⁺ ATPase alpha-1 Mouse mAb	WB/IHC-P
PTM-6271	Anti-Na ⁺ /K ⁺ ATPase alpha-1 Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6670	Anti-NCAM1/CD56 Mouse mAb	FC
PTM-6218	Anti-NCAM1/CD56 Rabbit mAb	WB/FC/IP
PTM-6041	Anti-NCS1 Mouse mAb	WB/ICC/IF
PTM-5402	Anti-NDRG2 Mouse mAb	WB/ICC/IF
PTM-5328	Anti-Nestin Mouse mAb	WB/IHC-P
PTM-6825	Anti-NeuN Rabbit mAb	WB/IHC-P/mIHC-P
PTM-7312	Anti-NeuroD1 Rabbit mAb	WB
PTM-7055	Anti-Neurofilament heavy polypeptide Rabbit mAb	WB/IHC-P
PTM-6874	Anti-Neurofilament-L Rabbit mAb	WB
PTM-6101	Anti-Neuron-specific beta-III Tubulin Mouse mAb	WB/IHC-P/ICC/IF
PTM-6750	Anti-Neuropeptide Y Rabbit pAb	IHC-P
PTM-7042	Anti-NFIB Rabbit mAb	WB/IHC-P
PTM-5808	Anti-NFκB p65 Rabbit mAb	WB
PTM-5591	Anti-NFκB p65 Mouse mAb	WB/ICC/IF
PTM-7324	Anti-Notch1 Rabbit mAb	WB/IHC-P
PTM-5352	Anti-NPM1 Mouse mAb	WB
PTM-6449	Anti-NSE Mouse mAb	WB
PTM-5519	Anti-NSE Rabbit mAb	WB/ICC/IF/FC/IP
PTM-6645	Anti-NSE Rabbit pAb	WB/FC
PTM-7013	Anti-NUMB Rabbit pAb	WB
PTM-7279	Anti-NUP107 Rabbit mAb	WB
PTM-5497	Anti-OGT Mouse mAb	WB/IHC-P/ICC/IF
PTM-6051	Anti-OPCML Rabbit mAb	WB
PTM-6902	Anti-Optineurin Rabbit mAb	WB
PTM-5470	Anti-P2RX4 Mouse mAb	WB/ICC/IF
PTM-6324	Anti-p44/42 MAPK (Erk1/2) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5721	Anti-p75NTR Mouse mAb	ICC/IF/FC
PTM-6395	Anti-p75NTR Rabbit mAb	WB/IHC-P
PTM-6122	Anti-p75NTR Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6201	Anti-PADI2 Mouse mAb	WB
PTM-6154	Anti-PAK3 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5881	Anti-PAK3 Rabbit pAb	IHC-P
PTM-7099	Anti-Pan Trk Rabbit mAb	WB/IHC-P
PTM-5182	Anti-Parvalbumin Mouse mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-7038	Anti-Parvalbumin Rabbit mAb	WB/IHC-P
PTM-7442	Anti-PAX5 Rabbit mAb	WB/ICC/IF/FC
PTM-7067	Anti-PCSK9 Rabbit mAb	WB
PTM-6062	Anti-PDE6D Mouse mAb	WB
PTM-7163	Anti-PEDF Rabbit mAb	WB
PTM-5863	Anti-PELP1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5735	Anti-PELP1 Rabbit mAb	WB/IHC-P
PTM-7174	Anti-Perilipin-1 Rabbit mAb	WB
PTM-6876	Anti-Peripherin Rabbit mAb	WB/IHC-P
PTM-7000	Anti-Peripherin Rabbit pAb	WB/IHC-P
PTM-5535	Anti-PGP 9.5 Mouse mAb	WB/ICC/IF
PTM-5520	Anti-PGP 9.5 Rabbit mAb	WB/IHC-P/IP
PTM-7039	Anti-PHF2 Rabbit mAb	WB
PTM-5975	Anti-PHGDH Mouse mAb	WB/ICC/IF
PTM-6615	Anti-Phospho-(Ser/Thr) PDK1 Docking Motif Mouse mAb	WB
PTM-6785	Anti-Phospho-AMPK beta 1 (Ser182) Rabbit pAb	WB
PTM-7423	Anti-Phospho-beta Catenin (Thr41/Ser45) Rabbit mAb	WB
PTM-6855	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit mAb	WB
PTM-6801	Anti-Phospho-Caveolin-1 (Tyr14) Rabbit pAb	WB/ICC/IF
PTM-6191	Anti-Phospho-CRMP2 (Thr514) Rabbit mAb	WB
PTM-6870	Anti-Phospho-EGFR (Tyr1092) Rabbit pAb	WB/ICC/IF
PTM-7373	Anti-Phospho-EGFR (Tyr1173) Rabbit mAb	WB/ICC/IF
PTM-6065	Anti-Phospho-FGFR1 (Tyr653/654) Mouse mAb	WB
PTM-7427	Anti-Phospho-GAP43 (Ser41) Rabbit mAb	WB
PTM-5230	Anti-Phospho-Girdin (Tyr1765) Rabbit mAb	IHC-P
PTM-7236	Anti-Phospho-IkB alpha (Tyr42) Rabbit mAb	WB
PTM-7425	Anti-Phospho-p53 (Thr377) Rabbit mAb	WB
PTM-7258	Anti-Phospho-Rab8A (Thr72) Rabbit mAb	WB
PTM-7426	Anti-Phospho-Synapsin-1 (Ser9) Rabbit mAb	WB
PTM-7041	Anti-Phospho-Vimentin (Ser39) Rabbit mAb	WB
PTM-7432	Anti-Phospho-Vimentin (Ser39) Rabbit mAb	WB/ICC/IF
PTM-6544	Anti-Pin1 Rabbit mAb	WB/ICC/IF/FC
PTM-6206	Anti-PKC delta Rabbit mAb	WB/ICC/IF
PTM-5433	Anti-PKC gamma Mouse mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7318	Anti-PLA2G6 Rabbit mAb	WB/IHC-P
PTM-7416	Anti-PLK1 Rabbit mAb	WB
PTM-6458	Anti-PLOD3 Mouse mAb	WB
PTM-5784	Anti-PLXNB2 Mouse mAb	WB/ICC/IF
PTM-6155	Anti-PPM1A Rabbit mAb	WB/IHC-P
PTM-5173	Anti-PPP3CB Mouse mAb	WB/IHC-P
PTM-6854	Anti-PPT1 Mouse mAb	WB
PTM-6796	Anti-PRAS40 Rabbit mAb	WB/IHC-P
PTM-5539	Anti-PRDX2 Mouse mAb	WB/ICC/IF
PTM-6525	Anti-PRDX3 Mouse mAb	WB/ICC/IF
PTM-7409	Anti-PRDX5 Mouse mAb	WB/IHC-P
PTM-5552	Anti-Prealbumin/Transthyretin Mouse mAb	WB/IHC-P
PTM-7327	Anti-Prealbumin/Transthyretin Rabbit mAb	WB/IHC-P
PTM-6270	Anti-Prion protein PrP Rabbit mAb	WB/IHC-P
PTM-7297	Anti-PROX1 Mouse mAb	WB/IHC-P
PTM-6242	Anti-PSA/KLK3 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5027	Anti-PSD-95 Mouse mAb	WB
PTM-5738	Anti-PSD95 Rabbit mAb	WB/IHC-P
PTM-5978	Anti-PSMA1 Rabbit mAb	WB/IP
PTM-5895	Anti-PTBP2 Rabbit mAb	WB/IHC-P
PTM-5605	Anti-PTGER2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5608	Anti-PTH Rabbit mAb	WB
PTM-7206	Anti-Rab29 Rabbit mAb	WB
PTM-7405	Anti-Rab3A Rabbit mAb	WB/IHC-P
PTM-5509	Anti-RAB3GAP1 Mouse mAb	WB/ICC/IF
PTM-6722	Anti-Rab5A Mouse mAb	WB/IHC-P
PTM-6958	Anti-Rab5A Rabbit mAb	WB
PTM-5907	Anti-RALYL Mouse mAb	WB/ICC/IF
PTM-5601	Anti-RAN Mouse mAb	WB/IHC-P/ICC/IF
PTM-6657	Anti-Ret Rabbit mAb	WB
PTM-5768	Anti-RIC8A Mouse mAb	WB
PTM-6989	Anti-Robo1 Rabbit mAb	WB
PTM-7187	Anti-RSK2 Rabbit mAb	WB
PTM-7303	Anti-RUNX1/AML1 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
PTM-6840	Anti-RUNX1/AML1 Rabbit pAb	WB
PTM-7025	Anti-RXRA Rabbit mAb	WB/ICC/IF
PTM-7164	Anti-Ryanodine Receptor Rabbit mAb	IHC-P
PTM-6250	Anti-S100A1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6283	Anti-S100A8 Rabbit mAb	WB/IHC-P/FC
PTM-7096	Anti-S100B Rabbit pAb	WB/IHC-P
PTM-5450	Anti-SCGN Rabbit mAb	IHC-P
PTM-6105	Anti-Serum Amyloid P/SAP Rabbit mAb	WB/IHC-P
PTM-7371	Anti-SHIP-2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6534	Anti-SHP2 Mouse mAb	WB/IHC-P/ICC/IF/FC
PTM-6969	Anti-SLC25A12 Rabbit mAb	WB/IHC-P
PTM-7172	Anti-SLK Mouse mAb	WB
PTM-5094	Anti-SNAP25 Mouse mAb	WB/IHC-P
PTM-5703	Anti-SOD2 Mouse mAb	WB/IHC-P
PTM-6113	Anti-SOX1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5876	Anti-SOX11 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5072	Anti-SOX2 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7289	Anti-SPR Rabbit mAb	WB/IHC-P
PTM-5508	Anti-SRF Mouse mAb	WB
PTM-5301	Anti-Stathmin Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5521	Anti-STK4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7210	Anti-STUB1 Rabbit mAb	WB/FC
PTM-5651	Anti-STXBP1 Mouse mAb	WB/IHC-P
PTM-3126	Anti-SUCLA2 Mouse mAb	WB
PTM-7138	Anti-SV2A Rabbit mAb	WB/IHC-P
PTM-7184	Anti-Synapsin-1 Mouse mAb	WB/IHC-P
PTM-6734	Anti-Synaptophysin Rabbit pAb	WB/IHC-P
PTM-6742	Anti-Syntaxin 1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6253	Anti-Syntaxin 10 Mouse mAb	WB
PTM-5999	Anti-Syntaxin 16 Mouse mAb	WB/ICC/IF
PTM-7214	Anti-Syntaxin 3 Rabbit mAb	WB/IHC-P
PTM-6228	Anti-Tau (phospho S262) Rabbit mAb	WB
PTM-6603	Anti-Tau Mouse mAb	WB/IHC-P
PTM-5151	Anti-Tau Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5063	Anti-Tau Rabbit mAb	IHC-P
PTM-7186	Anti-Tau(phospho S396) Rabbit mAb	WB/IHC-P
PTM-7142	Anti-TBK1/ NAK Rabbit mAb	WB
PTM-7146	Anti-TBR1 Rabbit mAb	WB/IHC-P
PTM-7422	Anti-TBX1 Rabbit mAb	WB
PTM-7280	Anti-TCF4 Rabbit mAb	WB
PTM-5469	Anti-TDP43 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6509	Anti-TDP43 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5998	Anti-TDRD3 Rabbit mAb	WB
PTM-5400	Anti-TH Mouse mAb	WB
PTM-5542	Anti-Thrombospondin 1 Mouse mAb	WB
PTM-6759	Anti-TIE2 Rabbit pAb	WB/IHC-P
PTM-6779	Anti-Tissue Factor Rabbit pAb	WB
PTM-5418	Anti-TMEM106B Mouse mAb	WB/IHC-P
PTM-6277	Anti-TMEM119 Mouse mAb	WB
PTM-5994	Anti-TMPRSS2 Rabbit mAb	WB/IHC-P/FC
PTM-6894	Anti-TMS1/ASC Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6834	Anti-TNF receptor II Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5964	Anti-TRIM24 Mouse mAb	WB/IHC-P
PTM-6325	Anti-TrkA (phospho Y791) Rabbit mAb	WB
PTM-6852	Anti-TrkC Mouse mAb	WB
PTM-6728	Anti-TRPV4 Rabbit pAb	WB/IHC-P
PTM-5583	Anti-TSHB Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7346	Anti-UBE3A Rabbit mAb	WB
PTM-7270	Anti-UBXN1 Rabbit mAb	WB/IHC-P
PTM-7238	Anti-ULK1 (phospho S623) Rabbit mAb	WB
PTM-6656	Anti-VAMP3 Mouse mAb	WB/IHC-P
PTM-5884	Anti-VAT1 Mouse mAb	WB/ICC/IF
PTM-6623	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-7105	Anti-VCAM1 Rabbit mAb	WB/IHC-P
PTM-5823	Anti-VCP Mouse mAb	WB/IHC-P/ICC/IF
PTM-7244	Anti-VGLUT1 Rabbit mAb	WB
PTM-5266	Anti-VPAC1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6862	Anti-VPS26 Rabbit pAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-7273	Anti-WFDC2 Rabbit mAb	WB
PTM-5894	Anti-WNT10B Mouse mAb	WB/IHC-P/ICC/IF
PTM-6039	Anti-WWOX Mouse mAb	WB/ICC/IF
PTM-5935	Anti-WWOX Rabbit mAb	WB/IHC-P/ICC/IF/FC

Cell Biology Research Antibodies

PTM-5462	Anti-14-3-3 theta Mouse mAb	WB/ICC/IF
PTM-5213	Anti-4-1BBL/TNFSF9 Mouse mAb	WB
PTM-6119	Anti-A1BG Mouse mAb	WB
PTM-6942	Anti-AACT Rabbit mAb	WB/IHC-P
PTM-6030	Anti-AATF Rabbit mAb	WB/ICC/IF
PTM-7158	Anti-ACE2 Rabbit mAb	WB/IHC-P
PTM-6542	Anti-ACE2 Rabbit pAb	WB/IHC-P
PTM-6794	Anti-Acetyl-SOD2 (Lys68) Rabbit mAb	WB/IHC-P
PTM-5791	Anti-Activin A Receptor Type IB/ALK-4 Rabbit mAb	WB/ICC/IF
PTM-5229	Anti-AFP Rabbit mAb	WB/FC
PTM-7253	Anti-AIF Mouse mAb	WB/IHC-P
PTM-6771	Anti-AIFM2/FSP1 Mouse mAb	WB/IHC-P
PTM-5191	Anti-AKT1 Mouse mAb	WB
PTM-6072	Anti-Albumin Mouse mAb	WB/ICC/IF
PTM-6099	Anti-ALIX Mouse mAb	WB
PTM-6407	Anti-ALIX Rabbit mAb	WB/IHC-P/IP
PTM-5306	Anti-alpha 2 Macroglobulin Mouse mAb	WB
PTM-5303	Anti-Alpha B Crystallin Rabbit mAb	WB/IHC-P
PTM-5216	Anti-alpha smooth muscle Actin Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5671	Anti-alpha smooth muscle Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6472	Anti-alpha Tubulin (acetyl K40) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5654	Anti-alpha tubulin (PE) Mouse mAb	ICC/IF/FC
PTM-5648	Anti-alpha tubulin (PE) Rabbit mAb	ICC/IF/FC
PTM-5282	Anti-alpha Tubulin Mouse mAb	WB/ICC/IF
PTM-5001	Anti-alpha Tubulin Mouse mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5442	Anti-alpha tubulin Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6607	Anti-ALPL Rabbit mAb	WB
PTM-5291	Anti-Androgen Receptor (AR) (C-term) Rabbit mAb	WB/IHC-P
PTM-5891	Anti-Annexin A1 Mouse mAb	WB/IHC-P
PTM-5423	Anti-Annexin A4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5424	Anti-Annexin V Mouse mAb	WB/IHC-P/ICC/IF
PTM-5767	Anti-ANP32A Mouse mAb	WB/IHC-P
PTM-5354	Anti-ANP32B Mouse mAb	WB/ICC/IF
PTM-7149	Anti-APAF1 Rabbit mAb	WB/IHC-P
PTM-5883	Anti-ARPC1B Mouse mAb	WB/ICC/IF
PTM-5988	Anti-ARPC5/p16 ARC Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6093	Anti-ASGR1 Mouse mAb	WB/ICC/IF
PTM-5579	Anti-ASGR1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6060	Anti-ASNA1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7043	Anti-ATG16L1 Mouse mAb	WB
PTM-5606	Anti-ATG4B Rabbit mAb	WB/ICC/IF
PTM-6296	Anti-ATG4B Rabbit mAb	WB/IHC-P/FC
PTM-6142	Anti-ATG5 Mouse mAb	WB/ICC/IF
PTM-6267	Anti-ATG7 Rabbit mAb	WB
PTM-6514	Anti-ATG9A Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6394	Anti-ATM Rabbit mAb	WB/IHC-P
PTM-5260	Anti-ATM Rabbit mAb	IHC-P/ICC/IF
PTM-5170	Anti-ATR Mouse mAb	WB
PTM-5308	Anti-Aurora A Mouse mAb	WB
PTM-5947	Anti-Aurora B Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6125	Anti-Aurora B Rabbit mAb	WB/IHC-P
PTM-5285	Anti-Axl Mouse mAb	WB
PTM-6881	Anti-Axl Rabbit mAb	WB
PTM-5275	Anti-B7-H3/CD276 Rabbit mAb	WB
PTM-6797	Anti-B7-H4 Rabbit mAb	WB
PTM-5420	Anti-BACH1 Mouse mAb	WB
PTM-6084	Anti-BAFF-R/TNFRSF13C Mouse mAb	ICC/IF/FC
PTM-6094	Anti-BAFF/TNFSF13B Mouse mAb	ICC/IF/FC
PTM-6371	Anti-BAG2 Rabbit mAb	WB

Cat.#	Product	Recommended Applications
	Anti-BAG6 Mouse mAb	WB/IHC-P/ICC/IF
	Anti-Bak Rabbit mAb	WB/IHC-P
	Anti-Bak Rabbit pAb	WB/IHC-P
	Anti-BAP31 Rabbit mAb	WB/IHC-P/ICC/IF
	Anti-BAX Rabbit mAb	IHC-P
	Anti-BCAS3 Mouse mAb	WB/ICC/IF
	Anti-Bcl-2 Mouse mAb	WB/ICC/IF
	Anti-Bcl-2 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
	Anti-Bcl-2 Rabbit mAb	WB/ICC/IF
	Anti-Bcl-XL Rabbit mAb	WB/IHC-P/ICC/IF/FC
	Anti-BCLG Mouse mAb	WB/IHC-P/ICC/IF
	Anti-beta Actin (FITC) Rabbit mAb	ICC/IF/FC
	Anti-beta Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
	Anti-beta Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC
	Anti-beta Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC
	Anti-beta Tubulin (HRP) Rabbit mAb	WB/IHC-P
	Anti-beta Tubulin Rabbit mAb	WB/IHC-P/ICC/IF/IP
	Anti-beta-Actin Mouse mAb	WB/IHC-P/ICC/IF
	Anti-BID Mouse mAb	WB/IHC-P
	Anti-Bid Rabbit mAb	WB/IHC-P/IP
	Anti-Bim Rabbit mAb	WB
	Anti-BMI1 Mouse mAb	WB
	Anti-BMP6 Rabbit mAb	WB/ICC/IF
	Anti-BTC Mouse mAb	WB/ICC/IF
	Anti-BTK Mouse mAb	WB
	Anti-BTN2A2/BTF2 Rabbit mAb	WB
	Anti-c-Cbl Mouse mAb	WB/ICC/IF
	Anti-c-IAP1 Mouse mAb	WB/ICC/IF
	Anti-c-IAP2 Rabbit mAb	WB/ICC/IF
	Anti-c-Maf Rabbit mAb	WB
	Anti-c-Myc Mouse mAb	WB
	Anti-Cadherin-17 Rabbit mAb	IHC-P
	Anti-Cadherin-17 Rabbit mAb	IHC-P
	Anti-Caldesmon/CDM Mouse mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5261	Anti-Caldesmon/CDM Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6809	Anti-Calpain 2 Mouse mAb	WB
PTM-5759	Anti-Calpastatin Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5233	Anti-Calponin 1 Rabbit mAb	WB/IHC-P
PTM-5264	Anti-Calponin 1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6389	Anti-Calreticulin Rabbit mAb	WB/IHC-P
PTM-5707	Anti-Calumenin Mouse mAb	WB/ICC/IF
PTM-5625	Anti-CAPRIN1 Mouse mAb	WB/IHC-P
PTM-5467	Anti-CAPZA1 Mouse mAb	WB
PTM-6202	Anti-Carbonic Anhydrase 2/CA2 Rat mAb	WB
PTM-6507	Anti-Casein Kinase 2 beta Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6614	Anti-Casein Kinase 2 beta Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7248	Anti-Caspase-2 Rabbit mAb	WB
PTM-5752	Anti-Caspase-3 Rabbit mAb	WB/IHC-P
PTM-6233	Anti-Caspase-6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6781	Anti-Caspase-7 Rabbit mAb	WB/IHC-P
PTM-5297	Anti-Caspase-7 Rabbit mAb	WB/IHC-P
PTM-6782	Anti-Caspase-7 Rabbit pAb	WB/IHC-P
PTM-6085	Anti-Caspase-8 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5311	Anti-Caspase-9 Mouse mAb	WB
PTM-6369	Anti-Cathepsin D Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7436	Anti-Cathepsin K Mouse mAb	WB/IHC-P
PTM-5356	Anti-Caveolin-1 Mouse mAb	WB
PTM-5025	Anti-Caveolin-1 Mouse mAb	WB
PTM-5696	Anti-CBX2 Rabbit mAb	WB
PTM-5022	Anti-CCR4 Mouse mAb	FC
PTM-5944	Anti-CCT5 Mouse mAb	WB/ICC/IF
PTM-5714	Anti-CCT8 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6489	Anti-Cdc23 Rabbit mAb	WB/ICC/IF
PTM-7410	Anti-Cdc25B Rabbit mAb	WB/IHC-P
PTM-6585	Anti-Cdc25C Mouse mAb	WB
PTM-6850	Anti-Cdc25C Rabbit mAb	WB
PTM-6315	Anti-Cdc27 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5476	Anti-Cdc37 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-6705	Anti-Cdc42 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5919	Anti-CDCP1 Mouse mAb	FC
PTM-6521	Anti-CDK1/Cdc2 Mouse mAb	WB/ICC/IF
PTM-7072	Anti-CDK1/Cdc2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6396	Anti-CDK2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5121	Anti-CDK2 Rabbit mAb	WB/IHC-P
PTM-5636	Anti-CDK4 Mouse mAb	WB
PTM-6131	Anti-CDK4 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5665	Anti-CDK6 Rabbit mAb	WB/ICC/IF/FC
PTM-5741	Anti-CDK8 Rabbit mAb	WB/IHC-P
PTM-5739	Anti-CDK9 Rabbit mAb	WB/IHC-P
PTM-5302	Anti-CDKN2A/p14ARF Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6704	Anti-CDKN2A/p16INK4a Rabbit mAb	WB
PTM-7154	Anti-CDKN2A/p16INK4a Rabbit mAb	WB/IHC-P
PTM-5724	Anti-CEACAM1 Mouse mAb	ICC/IF/FC
PTM-6118	Anti-CEP250 Mouse mAb	WB/ICC/IF
PTM-5893	Anti-CEP55 Mouse mAb	WB/ICC/IF
PTM-6536	Anti-CEP63 Mouse mAb	WB/ICC/IF
PTM-5972	Anti-CFTR Mouse mAb	IHC-P
PTM-7398	Anti-CHCHD2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6469	Anti-CHCHD6 Mouse mAb	WB/IHC-P
PTM-7275	Anti-CIRBP Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5399	Anti-Clathrin heavy chain Mouse mAb	WB/ICC/IF
PTM-5123	Anti-Claudin-1 Rabbit mAb	WB/IHC-P/IP
PTM-5134	Anti-Claudin-1 Rabbit mAb	WB/IHC-P
PTM-5586	Anti-Claudin-1 Rabbit mAb	WB/IHC-P
PTM-5679	Anti-Claudin-7 Rabbit mAb	WB/IHC-P
PTM-7246	Anti-Cleaved caspase-3 Rabbit mAb	WB
PTM-5272	Anti-Cleaved PARP1 p25 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7451	Anti-Cleaved PARP1 p25 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP/CUT&Tag
PTM-6014	Anti-Cleaved PARP1 p85 Rabbit mAb	WB/ICC/IF/FC/IP
PTM-7466	Anti-Cleaved PARP1 p85 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5530	Anti-CNPY2 Mouse mAb	WB/ICC/IF
PTM-5638	Anti-COPE Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5870	Anti-COPS3 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5712	Anti-Cpn10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6008	Anti-CSE1L Mouse mAb	WB/ICC/IF
PTM-6374	Anti-CSNK2A1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5339	Anti-CTAG1B Rabbit mAb	IHC-P
PTM-5653	Anti-CTBP1 Mouse mAb	WB
PTM-5220	Anti-CTEN Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6452	Anti-Cullin 1 Mouse mAb	WB
PTM-5500	Anti-CXCR7 Mouse mAb	WB/ICC/IF
PTM-6659	Anti-Cyclin B1 Rabbit mAb	WB/ICC/IF/IP
PTM-5936	Anti-Cyclin D1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6029	Anti-Cyclin D1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6664	Anti-Cyclin D2 Rabbit mAb	WB/IHC-P
PTM-5950	Anti-Cyclin D3 Mouse mAb	WB/ICC/IF
PTM-5138	Anti-Cyclin D3 Rabbit mAb	WB/IHC-P
PTM-6152	Anti-Cyclin H Mouse mAb	WB/IHC-P/ICC/IF
PTM-6048	Anti-Cyclin Y Mouse mAb	WB
PTM-5411	Anti-Cyclophilin B Mouse mAb	WB/IHC-P/ICC/IF
PTM-6685	Anti-Cyclophilin F Rabbit mAb	WB/IHC-P
PTM-5618	Anti-Cystatin A Rabbit mAb	WB/ICC/IF/FC
PTM-5160	Anti-cytochrome C Mouse mAb	WB/IHC-P
PTM-5351	Anti-Cytochrome C Mouse mAb	WB/IHC-P
PTM-7181	Anti-Cytoplasmic Dynein Intermediate chain Mouse mAb	WB
PTM-5498	Anti-DCPS Mouse mAb	WB/IHC-P
PTM-7183	Anti-DcR1 Rabbit mAb	WB/ICC/IF
PTM-5525	Anti-DDB1 Mouse mAb	WB/IHC-P
PTM-5059	Anti-delta 1 Catenin Rabbit mAb	IHC-P
PTM-5188	Anti-Desmin Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6539	Anti-Desmin Rabbit mAb	WB/IHC-P
PTM-6063	Anti-Desmocollin-1 Rat mAb	WB/IHC-P
PTM-6171	Anti-DHCR24 Mouse mAb	WB/ICC/IF
PTM-6708	Anti-DHFR Mouse mAb	WB
PTM-5905	Anti-DNAJA1 Mouse mAb	WB/ICC/IF
PTM-5673	Anti-DPEP1 Rabbit mAb	IHC-P

Cat.#	Product	Recommended Applications
PTM-6281	Anti-DPEP3 Rabbit mAb	WB/FC
PTM-7057	Anti-DPP4/CD26 Rabbit mAb	WB/IHC-P
PTM-6695	Anti-DR5 Rabbit mAb	WB/IHC-P
PTM-6007	Anti-DRG2 Mouse mAb	WB
PTM-5602	Anti-DYNLL1 Rabbit mAb	WB/IHC-P
PTM-5024	Anti-E Cadherin Mouse mAb	WB/IHC-P/ICC/IF
PTM-6222	Anti-E-Cadherin Rabbit mAb	WB/IP
PTM-7307	Anti-E-Syt1 Mouse mAb	WB
PTM-6745	Anti-EGFR (phospho S1166) Rabbit pAb	WB
PTM-5042	Anti-EGFR Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5438	Anti-EGFR Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7360	Anti-EID1 Rabbit mAb	WB/IHC-P
PTM-5728	Anti-eIF3L Mouse mAb	WB/ICC/IF
PTM-5887	Anti-eIF4G2/DAP5 Mouse mAb	WB/ICC/IF
PTM-6214	Anti-EpCAM Mouse mAb	FC
PTM-5871	Anti-EpCAM Rabbit mAb	WB/IHC-P
PTM-5280	Anti-EphA2 Rabbit mAb	WB/IHC-P
PTM-5430	Anti-EPHX2 Mouse mAb	WB
PTM-5017	Anti-ErbB2 Mouse mAb	WB
PTM-5061	Anti-ERCC1 Rabbit mAb	IHC-P
PTM-6013	Anti-eRF3/GSPT1 Rabbit mAb	WB
PTM-5480	Anti-ERO1L Mouse mAb	WB/ICC/IF
PTM-5504	Anti-ERp29 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6531	Anti-ERp44 Mouse mAb	WB/IHC-P
PTM-6068	Anti-ESD Mouse mAb	WB/ICC/IF
PTM-6290	Anti-EXOC3L2 Mouse mAb	WB
PTM-5574	Anti-Factor XIIIa Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6367	Anti-FADD (phospho S194) Rabbit mAb	WB
PTM-5358	Anti-FAP Rabbit mAb	IHC-P
PTM-5269	Anti-FBXW7 Rabbit mAb	IHC-P/ICC/IF
PTM-5276	Anti-FGFR2 Rabbit mAb	IHC-P
PTM-5958	Anti-Fibrillarin Rabbit mAb	WB/ICC/IF
PTM-6050	Anti-Filamin A Rat mAb	WB/ICC/IF
PTM-6379	Anti-FKBP8 Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5472	Anti-FKBPL Mouse mAb	WB/IHC-P
PTM-5599	Anti-FLII Mouse mAb	WB/IHC-P/ICC/IF
PTM-5369	Anti-Flotillin Mouse mAb	WB/ICC/IF
PTM-7364	Anti-FNDC5 Rabbit mAb	WB
PTM-5878	Anti-FUBP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5273	Anti-Galectin-1 Rabbit mAb	WB/IHC-P
PTM-6356	Anti-Galectin-3 Mouse mAb	WB
PTM-6818	Anti-Galectin-9 Rabbit pAb	WB/IHC-P
PTM-6620	Anti-GAPDH (HRP) Rabbit mAb	WB/IHC-P
PTM-5251	Anti-GAPDH Mouse mAb	WB
PTM-5150	Anti-GAPDH Mouse mAb	WB
PTM-5375	Anti-GAPDH Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5133	Anti-GCDFP-15 Rabbit mAb	IHC-P
PTM-5817	Anti-GCK Mouse mAb	WB/IHC-P/ICC/IF
PTM-6959	Anti-GCLC Rabbit mAb	WB
PTM-6178	Anti-GCLM Mouse mAb	WB/IHC-P
PTM-6104	Anti-GFAP (NT) Rabbit mAb	WB/IHC-P/mIHC-P/ICC/IF
PTM-5180	Anti-GFP Mouse mAb	WB
PTM-6673	Anti-GFP Rabbit mAb	WB
PTM-3024	Anti-GLO1 Rabbit pAb	WB
PTM-3561	Anti-GLO2 Rabbit pAb	WB
PTM-5632	Anti-GLUD1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6992	Anti-GLUT3+GLUT14 Rabbit mAb	WB
PTM-5585	Anti-Glycophorin A Rabbit mAb	IHC-P/FC
PTM-5610	Anti-Glycophorin C Rabbit mAb	WB/IHC-P
PTM-5923	Anti-GOLGA2/GM130 Mouse mAb	WB
PTM-5678	Anti-GOT1 Rabbit mAb	WB/IHC-P
PTM-5246	Anti-Granzyme B Mouse mAb	WB
PTM-5588	Anti-Granzyme B Rabbit mAb	WB
PTM-5779	Anti-GRP78/BiP Mouse mAb	WB/IHC-P/ICC/IF
PTM-6107	Anti-GSK3 alpha Rabbit mAb	WB/ICC/IF/IP
PTM-5855	Anti-GSK3 alpha Rat mAb	WB
PTM-6635	Anti-GSK3 alpha/beta Mouse mAb	WB/IHC-P/ICC/IF
PTM-5695	Anti-GSS Mouse mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5181	Anti-GST Tag Mouse mAb	WB
PTM-5046	Anti-GST Tag Mouse mAb	WB
PTM-6598	Anti-GST3/GSTP1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6391	Anti-GSTK1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6044	Anti-GSTO1 Mouse mAb	WB/ICC/IF
PTM-7212	Anti-HACE1 Rabbit mAb	WB
PTM-7193	Anti-Hamartin Rabbit mAb	WB
PTM-6380	Anti-HEC1 Mouse mAb	WB
PTM-5225	Anti-hENT1 Rabbit mAb	IHC-P
PTM-5452	Anti-HER2/ErbB2 Rabbit mAb	WB/IHC-P/FC
PTM-5458	Anti-HER2/ErbB2 Rabbit mAb	WB/IHC-P
PTM-5623	Anti-HIF-1 beta / ARNT Mouse mAb	WB
PTM-6259	Anti-HIP2/LIG Rabbit mAb	WB
PTM-5262	Anti-HLA-DPA1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5263	Anti-HLA-DPB1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6213	Anti-HP1 alpha Rabbit mAb	WB
PTM-6640	Anti-HPT Mouse mAb	WB
PTM-5431	Anti-HSP20 Mouse mAb	WB/IHC-P
PTM-6116	Anti-HSP60 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5685	Anti-HSP60 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5692	Anti-HSP70 Mouse mAb	WB
PTM-5684	Anti-HSP70 Rabbit mAb	WB/IHC-P/FC
PTM-6934	Anti-HSP70 Rabbit pAb	WB/IHC-P
PTM-5409	Anti-HSP90 alpha Mouse mAb	WB/IHC-P
PTM-5097	Anti-HSP90 beta Mouse mAb	WB/IHC-P
PTM-5734	Anti-HSP90 beta Rabbit mAb	WB/IHC-P/IP
PTM-6970	Anti-HSP90 beta Rabbit pAb	WB/IHC-P
PTM-5208	Anti-Human alpha 1-Acid Glycoprotein Mouse mAb	WB
PTM-7108	Anti-Iba1 Rabbit mAb	IHC-P/mIHC-P
PTM-6749	Anti-Iba1/AIF-1 Rabbit mAb	WB/IHC-P
PTM-7237	Anti-IFI16 Mouse mAb	WB
PTM-6475	Anti-IGFBP3 Rabbit mAb	WB
PTM-5020	Anti-IkB alpha Mouse mAb	WB/IHC-P
PTM-5829	Anti-IKK alpha Rabbit mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-7107	Anti-IKK alpha Rabbit mAb	WB
PTM-5701	Anti-IKK gamma Mouse mAb	WB/IHC-P
PTM-6546	Anti-IKK gamma Rabbit mAb	WB/ICC/IF/FC/IP
PTM-5241	Anti-IL-6 Mouse mAb	WB
PTM-5212	Anti-IL-7R Mouse mAb	FC
PTM-5898	Anti-ILK Mouse mAb	WB/IHC-P/ICC/IF
PTM-5762	Anti-INF2 mouse mAb	WB/IHC-P
PTM-6460	Anti-ING1 Rabbit mAb	WB
PTM-5550	Anti-Integrin alpha 2 / CD49b Mouse mAb	FC
PTM-5429	Anti-Integrin alpha 3 Mouse mAb	WB/IHC-P
PTM-5548	Anti-Integrin alpha 4 / CD49d Mouse mAb	FC
PTM-5597	Anti-INTS3 Mouse mAb	WB
PTM-5890	Anti-IST1 Mouse mAb	WB/IHC-P
PTM-7135	Anti-JAK2 Rabbit mAb	WB
PTM-6456	Anti-KAP1 Mouse mAb	WB/IHC-P
PTM-5760	Anti-Keratin 1 Rabbit pAb	WB/IHC-P
PTM-6678	Anti-Keratin 10 Rabbit mAb	WB/IHC-P
PTM-5048	Anti-Keratin 10 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5943	Anti-Keratin 14 Mouse mAb	IHC-P/ICC/IF
PTM-5391	Anti-Keratin 14 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5049	Anti-Keratin 14 Rabbit mAb	WB/IHC-P
PTM-5986	Anti-Keratin 15 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5341	Anti-Keratin 16 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5340	Anti-Keratin 16 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5067	Anti-Keratin 17 Rabbit mAb	WB/IHC-P/IP
PTM-5062	Anti-Keratin 18 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6508	Anti-Keratin 18 Rabbit mAb	WB/IHC-P
PTM-5796	Anti-Keratin 19 Mouse mAb	WB/ICC/IF
PTM-5645	Anti-Keratin 19 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5611	Anti-Keratin 20 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5454	Anti-Keratin 20 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5050	Anti-Keratin 20 Rabbit mAb	IHC-P
PTM-5012	Anti-Keratin 4 Rabbit mAb	IHC-P/mIHC-P
PTM-6086	Anti-Keratin 4 Rabbit mAb	IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5004	Anti-Keratin 5 Mouse mAb	IHC-P
PTM-5875	Anti-Keratin 5 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5040	Anti-Keratin 5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5080	Anti-Keratin 5 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5492	Anti-Keratin 5/6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5079	Anti-Keratin 6 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5011	Anti-Keratin 7 Rabbit mAb	WB/IHC-P/IP
PTM-5107	Anti-Keratin 8 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6058	Anti-Keratin 8 Rabbit mAb	WB/ICC/IF/FC
PTM-5491	Anti-Keratin HMW Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5490	Anti-Keratin LMW Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5690	Anti-KIF20A Mouse mAb	WB/ICC/IF
PTM-5725	Anti-KLK11 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5425	Anti-KRAS Mouse mAb	WB/ICC/IF
PTM-5287	Anti-KRT15 Mouse mAb	WB/ICC/IF
PTM-7274	Anti-L2HGDH Mouse mAb	WB
PTM-6036	Anti-LACTB2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5746	Anti-Lamin A/C Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6432	Anti-Lamin A/C Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5639	Anti-Lamin A/C Rabbit pAb	WB/IHC-P
PTM-6357	Anti-Lamin B2 Rabbit mAb	WB/IHC-P/IP
PTM-6777	Anti-LATS1 Rabbit mAb	WB
PTM-6829	Anti-LDHA Rabbit pAb	WB/IHC-P
PTM-6596	Anti-Legumain Mouse mAb	WB/IHC-P
PTM-7359	Anti-LIF Rat mAb	WB/IHC-P
PTM-5687	Anti-LPP Mouse mAb	WB
PTM-6610	Anti-LRP1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6481	Anti-LSM12 Mouse mAb	WB
PTM-5255	Anti-LYPD3 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5290	Anti-LYPD3 Rabbit mAb	IHC-P
PTM-6179	Anti-MAD2L1 Mouse mAb	WB/IHC-P
PTM-5516	Anti-MAD2L2 Mouse mAb	WB
PTM-6190	Anti-MALT1/MLT Rabbit mAb	WB/ICC/IF
PTM-5619	Anti-mCherry Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-6491	Anti-MCM2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6164	Anti-MCM2 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5844	Anti-MCM2 Rabbit mAb	WB/IHC-P
PTM-6804	Anti-MCM2 Rabbit pAb	WB
PTM-5733	Anti-MCM3 Rabbit mAb	WB/IHC-P/IP
PTM-6917	Anti-MCM3 Rabbit pAb	WB
PTM-7310	Anti-MDM2 Mouse mAb	WB
PTM-5873	Anti-MEK1 Rabbit mAb	WB/IHC-P/IP
PTM-5397	Anti-MEK2 Mouse mAb	WB
PTM-5963	Anti-MEK4/MKK4 Mouse mAb	WB/ICC/IF
PTM-5219	Anti-Mesothelin Rabbit mAb	WB/IP
PTM-6537	Anti-Mesothelin Rabbit mAb	WB
PTM-7177	Anti-Met Mouse mAb	WB
PTM-5448	Anti-Met Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-7185	Anti-MFF Rabbit mAb	WB
PTM-6215	Anti-MIS12 Mouse mAb	WB
PTM-5155	Anti-MITF Mouse mAb	WB
PTM-7151	Anti-Mitofusin2+Mitofusin1 Mouse mAb	WB/IHC-P
PTM-5818	Anti-MKS1 Mouse mAb	WB/ICC/IF
PTM-7375	Anti-MMP7 Rabbit mAb	WB/IHC-P
PTM-6764	Anti-MMP8 Rabbit mAb	IHC-P
PTM-7259	Anti-MMP9 Rabbit mAb	WB/IHC-P
PTM-5726	Anti-MPP8 rabbit mAb	WB
PTM-5874	Anti-MPTX Rabbit mAb	WB/IHC-P
PTM-5694	Anti-MTHFD1 Mouse mAb	WB/ICC/IF
PTM-6594	Anti-mTOR Mouse mAb	WB/IHC-P/ICC/IF
PTM-5373	Anti-MUC1 Rabbit mAb	WB/IHC-P
PTM-5995	Anti-MUC2 Rabbit mAb	WB/IHC-P/mIHC-P/FC
PTM-5372	Anti-MUC2 Rabbit pAb	IHC-P/ICC/IF
PTM-5271	Anti-MUC4 Rabbit mAb	IHC-P
PTM-5390	Anti-Myc tag Rabbit mAb	WB/ICC/IF/IP
PTM-5487	Anti-Myeloperoxidase Rabbit mAb	WB/IP
PTM-5065	Anti-Myeloperoxidase Rabbit mAb	IHC-P
PTM-5795	Anti-Myosin Light Chain 3 Mouse mAb	WB/IHC-P

Cat.#	Product	Recommended Applications
PTM-5443	Anti-Myosin Light Chain 3 Rabbit mAb	WB/IHC-P/FC
PTM-5221	Anti-N-Cadherin Rabbit mAb	WB/IHC-P/IP
PTM-5415	Anti-NAF1/CISD2 Mouse mAb	WB/ICC/IF
PTM-5774	Anti-NAPRT1 Mouse mAb	WB
PTM-5573	Anti-Napsin A Rabbit mAb	IHC-P
PTM-5305	Anti-NCAM1/CD56 Mouse mAb	WB/ICC/IF/FC
PTM-5765	Anti-NDUFA4L2 Mouse mAb	WB/ICC/IF
PTM-7189	Anti-NDUFS3 Rabbit mAb	WB/IHC-P
PTM-6883	Anti-Nectin-4 Rabbit pAb	WB
PTM-5070	Anti-Nestin Rabbit mAb	IHC-P
PTM-5706	Anti-NF-κB1 p105/p50 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6593	Anti-NF-κB2 p100/p52 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5254	Anti-NFκB p65 Mouse mAb	WB
PTM-5268	Anti-NPM1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5831	Anti-NRBP2 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6216	Anti-NUFIP2 Mouse mAb	WB
PTM-6000	Anti-NUP155 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5716	Anti-NUP50 Mouse mAb	WB/IHC-P
PTM-5547	Anti-NUP62 Mouse mAb	WB/IHC-P
PTM-5885	Anti-NXT1 Mouse mAb	WB
PTM-6335	Anti-OMA1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6176	Anti-OXA1L Mouse mAb	WB/ICC/IF
PTM-5571	Anti-p130 Cas Rabbit mAb	WB/FC
PTM-6510	Anti-p21 Rabbit mAb	WB/IHC-P/IP
PTM-7125	Anti-p21 Rabbit mAb	WB
PTM-6471	Anti-p27 Kip1 Rabbit mAb	WB/IHC-P
PTM-5470	Anti-P2RX4 Mouse mAb	WB/ICC/IF
PTM-7064	Anti-p40 Rabbit mAb	WB
PTM-5850	Anti-p44/42 MAPK (Erk1/2) Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6324	Anti-p44/42 MAPK (Erk1/2) Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5859	Anti-P4HA1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6257	Anti-p63 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6440	Anti-Pan-Actin (FITC) Rabbit mAb	FC
PTM-5447	Anti-Pan-Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC

Cat.#	Product	Recommended Applications
PTM-6015	Anti-PAP Rabbit mAb	WB/IHC-P/FC
PTM-6249	Anti-PAP Rabbit mAb	IHC-P
PTM-5171	Anti-PARP1 Mouse mAb	WB
PTM-5989	Anti-PARP1 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7431	Anti-PARP4 Rabbit mAb	WB
PTM-5033	Anti-PAX5 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6663	Anti-PAX5 Rabbit mAb	WB/IHC-P/IP
PTM-5711	Anti-PCAF Rabbit mAb	WB/IP
PTM-5819	Anti-PCMT1 Mouse mAb	WB/ICC/IF
PTM-6410	Anti-PDCD2 Rabbit mAb	WB
PTM-6482	Anti-PDCD4 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5512	Anti-PDIA6 Mouse mAb	WB/ICC/IF
PTM-5292	Anti-PDK1 Rabbit mAb	IHC-P/ICC/IF
PTM-6846	Anti-Pericentrin Rabbit pAb	WB
PTM-5535	Anti-PGP 9.5 Mouse mAb	WB/ICC/IF
PTM-5520	Anti-PGP 9.5 Rabbit mAb	WB/IHC-P/IP
PTM-5788	Anti-Phospho S6 Ribosomal Protein (Ser235/236) Mouse mAb	WB/ICC/IF
PTM-6615	Anti-Phospho-(Ser/Thr) PDK1 Docking Motif Mouse mAb	WB
PTM-6986	Anti-Phospho-Axl (Tyr702) Rabbit pAb	WB
PTM-7083	Anti-Phospho-Bcl-2 (Thr56) Rabbit pAb	WB
PTM-6887	Anti-Phospho-CAD (Ser1859) Rabbit pAb	WB
PTM-7417	Anti-Phospho-Cdc37 (Ser13) Rabbit mAb	WB
PTM-6786	Anti-Phospho-Cdc37 (Ser13) Rabbit pAb	WB
PTM-6185	Anti-Phospho-Cyclin D1 (Thr286) Rabbit mAb	WB/IHC-P/FC
PTM-6870	Anti-Phospho-EGFR (Tyr1092) Rabbit pAb	WB/ICC/IF
PTM-7373	Anti-Phospho-EGFR (Tyr1173) Rabbit mAb	WB/ICC/IF
PTM-6985	Anti-Phospho-IkB α (Ser36) Rabbit mAb	WB
PTM-7350	Anti-Phospho-NEDD4L (Ser448) Rabbit mAb	WB
PTM-6886	Anti-Phospho-p38 MAPK (Thr180/Tyr182) Rabbit mAb	WB
PTM-6047	Anti-Phospho-p53 (Ser15) Mouse mAb	WB/ICC/IF
PTM-7425	Anti-Phospho-p53 (Thr377) Rabbit mAb	WB
PTM-7126	Anti-Phospho-PAK1 (Thr212) Rabbit mAb	WB
PTM-7381	Anti-Phospho-Rb(Thr826) Rabbit mAb	WB
PTM-6224	Anti-Phospho-RIPK1(Ser161) Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5782	Anti-Phospho-RNA Polymerase II CTD Repeat YSPTSPS (Ser2/5) Mouse mAb	WB
PTM-6601	Anti-Phospho-SIRT1(Ser27) Rabbit mAb	WB/ICC/IF
PTM-6595	Anti-Phospho-SIRT1(Ser47) Rabbit mAb	WB/IHC-P
PTM-5644	Anti-Phospho-Tau (Ser199) Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5228	Anti-PI3 Kinase p110 alpha Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5198	Anti-PI3 Kinase p85 alpha Mouse mAb	WB
PTM-5426	Anti-PI3 Kinase p85 Beta Mouse mAb	WB/IHC-P
PTM-6544	Anti-Pin1 Rabbit mAb	WB/ICC/IF/FC
PTM-7162	Anti-PKA C-alpha Mouse mAb	WB/IHC-P
PTM-5433	Anti-PKC gamma Mouse mAb	WB/IHC-P/ICC/IF
PTM-5852	Anti-PKN1 Rabbit mAb	WB
PTM-5816	Anti-PKR Mouse mAb	WB
PTM-5952	Anti-PLK1 Mouse mAb	WB/ICC/IF
PTM-5659	Anti-PLK1 Rabbit mAb	WB
PTM-7416	Anti-PLK1 Rabbit mAb	WB
PTM-6155	Anti-PPM1A Rabbit mAb	WB/IHC-P
PTM-6796	Anti-PRAS40 Rabbit mAb	WB/IHC-P
PTM-6602	Anti-PRC1 Mouse mAb	WB/IHC-P
PTM-5951	Anti-PRC1 Mouse mAb	WB
PTM-5367	Anti-PRDX1 Mouse mAb	WB/IHC-P
PTM-5934	Anti-PRDX1 Rabbit mAb	WB/IHC-P/FC
PTM-5539	Anti-PRDX2 Mouse mAb	WB/ICC/IF
PTM-6404	Anti-PRDX4 Mouse mAb	WB
PTM-6522	Anti-PRDX6 Mouse mAb	WB
PTM-6538	Anti-PRK2 Rabbit mAb	WB
PTM-6120	Anti-Prohibitin Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6719	Anti-PSMA3 Mouse mAb	WB/IHC-P
PTM-5813	Anti-PSMA6 Mouse mAb	WB/ICC/IF
PTM-5783	Anti-PSMB8 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5864	Anti-PSMD9 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5119	Anti-PTEN Rabbit mAb	WB/IHC-P
PTM-5748	Anti-PTEN Rabbit mAb	WB/IHC-P/FC
PTM-7438	Anti-PWP1 Rabbit mAb	WB/IHC-P
PTM-5210	Anti-PYGO2 Mouse mAb	WB

Cat.#	Product	Recommended Applications
PTM-5575	Anti-QSOX1 Rabbit mAb	WB/FC/IP
PTM-5089	Anti-Rab10 Mouse mAb	WB
PTM-7241	Anti-Rab35 Rabbit mAb	WB
PTM-5815	Anti-Rab9A Rabbit mAb	WB/ICC/IF
PTM-5329	Anti-Rac1 Mouse mAb	WB
PTM-7443	Anti-Rad21 Rabbit mAb	WB/IHC-P/ICC/IF/FC/ChIP
PTM-6069	Anti-Rad23B Mouse mAb	WB/IHC-P
PTM-5601	Anti-RAN Mouse mAb	WB/IHC-P/ICC/IF
PTM-5824	Anti-RARS Mouse mAb	WB
PTM-6773	Anti-Ras Mouse mAb	WB/IHC-P
PTM-5493	Anti-Rb Mouse mAb	WB/IHC-P
PTM-6474	Anti-RBX1 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5460	Anti-RCC1 Mouse mAb	WB/IHC-P
PTM-7009	Anti-RING2/RING1B/RNF2 Rabbit pAb	WB
PTM-7272	Anti-RNF126 Mouse mAb	WB
PTM-7197	Anti-RNF40 Rabbit mAb	WB
PTM-5787	Anti-RNH1 Mouse mAb	WB
PTM-6454	Anti-ROCK1 Mouse mAb	WB/IHC-P
PTM-5379	Anti-ROCK2 Mouse mAb	WB
PTM-6169	Anti-ROCK2 Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6926	Anti-ROCK2+ROCK1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6828	Anti-RON Rabbit mAb	WB
PTM-7167	Anti-RON Rabbit mAb	WB
PTM-5933	Anti-RPA32/RPA2 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5621	Anti-RPL17 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6433	Anti-RPL26L1 Rabbit mAb	WB/IHC-P
PTM-5763	Anti-RPS12 Mouse mAb	WB/ICC/IF
PTM-5835	Anti-RPS14 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5834	Anti-RPS20 Rabbit mAb	WB/IHC-P
PTM-7068	Anti-RPS3 Rabbit mAb	WB
PTM-5146	Anti-RRM2 Mouse mAb	WB/ICC/IF
PTM-5594	Anti-Rubisco Mouse mAb	WB
PTM-7182	Anti-Rubisco Mouse mAb	WB
PTM-6250	Anti-S100A1 Rabbit mAb	WB/IHC-P/ICC/IF/FC

Cat.#	Product	Recommended Applications
PTM-7233	Anti-S100A11 Mouse mAb	WB
PTM-5677	Anti-S100A13 Rabbit mAb	IHC-P/ICC/IF/FC
PTM-5809	Anti-S100A2 Rabbit mAb	WB/IHC-P
PTM-5122	Anti-S100B Rabbit mAb	WB/IHC-P
PTM-5113	Anti-S100B Rabbit mAb	WB/IHC-P/ICC/IF
PTM-7096	Anti-S100B Rabbit pAb	WB/IHC-P
PTM-5451	Anti-S100P Rabbit mAb	IHC-P/FC
PTM-5592	Anti-S6 Ribosomal Protein Mouse mAb	WB/ICC/IF
PTM-5862	Anti-S6 Ribosomal Protein Mouse mAb	WB
PTM-6933	Anti-SA2 Rabbit mAb	WB/IHC-P
PTM-7308	Anti-SAPK/JNK (phospho T183/Y185) Mouse mAb	WB
PTM-5805	Anti-SATB1 Rabbit mAb	WB/IHC-P/FC
PTM-5131	Anti-SATB1 Rabbit mAb	IHC-P
PTM-5294	Anti-SATB1 Rabbit mAb	IHC-P/ICC/IF
PTM-6461	Anti-SEC11A Mouse mAb	WB
PTM-7220	Anti-SENP1 Rabbit mAb	WB
PTM-5434	Anti-SEPT2 Mouse mAb	WB/IHC-P
PTM-6147	Anti-Septin 7 Mouse mAb	WB/ICC/IF
PTM-6740	Anti-SERPINB9 Mouse mAb	WB
PTM-7234	Anti-SET/TAF-I Rabbit mAb	WB/IHC-P
PTM-7229	Anti-SFT Rabbit mAb	WB
PTM-5825	Anti-SHARPIN Mouse mAb	WB
PTM-6813	Anti-SKP2 Mouse mAb	WB
PTM-5554	Anti-SLAMF6/CD352 Mouse mAb	ICC/IF/FC
PTM-6198	Anti-SLC25A17 Mouse mAb	WB/IHC-P
PTM-6753	Anti-SLFN11 Rabbit mAb	WB
PTM-6127	Anti-Smac/Diablo Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5977	Anti-Smac/Diablo Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5842	Anti-SMARCA2/BRM Rabbit mAb	WB/IHC-P
PTM-7160	Anti-SMC1A Rabbit mAb	WB
PTM-5845	Anti-smooth muscle Myosin heavy chain 11 Rabbit mAb	WB/IHC-P
PTM-5406	Anti-SND1 Mouse mAb	WB
PTM-6097	Anti-SNX5 Mouse mAb	WB/ICC/IF
PTM-5363	Anti-SOX10 Rabbit mAb	WB/IHC-P/ICC/IF

Cat.#	Product	Recommended Applications
PTM-6736	Anti-SPARC Rat mAb	WB
PTM-5350	Anti-SPINK1 Rabbit mAb	IHC-P
PTM-5457	Anti-ST13 Rabbit mAb	WB/IHC-P/FC
PTM-7211	Anti-StAR Rabbit mAb	WB/IHC-P
PTM-5090	Anti-STAT3 Mouse mAb	WB/IHC-P
PTM-6025	Anti-STBD1 Mouse mAb	WB/IHC-P
PTM-5360	Anti-STING Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5359	Anti-STING Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6114	Anti-STIP1/ST11 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5398	Anti-STK11 Mouse mAb	WB
PTM-5408	Anti-STK3 Mouse mAb	WB
PTM-5708	Anti-STRAP Mouse mAb	WB/IHC-P/ICC/IF
PTM-5892	Anti-STT3A Mouse mAb	WB/IHC-P/ICC/IF
PTM-5427	Anti-Syntaxin 6 Mouse mAb	WB
PTM-6212	Anti-SYVN1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-7269	Anti-TACC3 Rabbit mAb	WB/IHC-P
PTM-5274	Anti-TCL1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5821	Anti-TDRKH Mouse mAb	WB
PTM-5745	Anti-Tenascin X Rabbit mAb	IHC-P
PTM-5661	Anti-TFEB Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-6064	Anti-TFIP11 Mouse mAb	WB
PTM-5773	Anti-TGF beta Receptor II Mouse mAb	WB/IHC-P/ICC/IF
PTM-5822	Anti-TGN46 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5357	Anti-Thrombomodulin Rabbit mAb	IHC-P/ICC/IF
PTM-5518	Anti-TIA1 Rabbit mAb	WB/FC/IP
PTM-7285	Anti-Tissue Plasminogen Activator Rabbit mAb	WB/IHC-P
PTM-5192	Anti-TLR4 Mouse mAb	WB
PTM-5463	Anti-TMEM111 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6641	Anti-TMIGD2/IGPR1 Rabbit mAb	WB/IP
PTM-6894	Anti-TMS1/ASC Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6834	Anti-TNF receptor II Rabbit mAb	WB/IHC-P/ICC/IF
PTM-6200	Anti-TNFAIP1 Mouse mAb	WB/ICC/IF
PTM-6165	Anti-TPPII Mouse mAb	WB/IHC-P
PTM-6405	Anti-TPRKB Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5414	Anti-TPT1 Mouse mAb	WB/IHC-P
PTM-6056	Anti-TPX2 Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-7354	Anti-TRADD Mouse mAb	WB/IHC-P
PTM-5534	Anti-TRAF2 Mouse mAb	WB/ICC/IF
PTM-5047	Anti-TRAF4 Mouse mAb	WB
PTM-5346	Anti-TRBC1 Rabbit mAb	WB
PTM-7339	Anti-TRIAD1 Rabbit mAb	WB
PTM-5576	Anti-TROP2 Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-5344	Anti-TROP2 Rabbit mAb	WB/IHC-P
PTM-5108	Anti-TSG101 Mouse mAb	WB
PTM-7030	Anti-TTK Mouse mAb	WB
PTM-6671	Anti-TXNRD1 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6441	Anti-UBAC1 Mouse mAb	WB
PTM-5513	Anti-UBC12 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5652	Anti-UBC9 Mouse mAb	WB
PTM-6376	Anti-UBE2C Mouse mAb	WB/IHC-P/ICC/IF
PTM-7235	Anti-UBE2L3 Rabbit mAb	WB/IHC-P
PTM-7346	Anti-UBE3A Rabbit mAb	WB
PTM-7270	Anti-UBXN1 Rabbit mAb	WB/IHC-P
PTM-5837	Anti-UGDH Mouse mAb	WB/IHC-P/ICC/IF
PTM-6847	Anti-UQCRFS1/RISP Rabbit mAb	WB/IHC-P/IP
PTM-5897	Anti-USP15 Mouse mAb	WB/IHC-P/ICC/IF
PTM-5830	Anti-USP5 Mouse mAb	WB/ICC/IF
PTM-5761	Anti-USP8 Mouse mAb	WB/IHC-P
PTM-5278	Anti-VEGF Receptor2 Rabbit mAb	WB
PTM-5392	Anti-Villin Rabbit mAb	WB/IHC-P
PTM-5029	Anti-Villin Rabbit mAb	IHC-P
PTM-5376	Anti-Vimentin Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP/mIHC-P
PTM-5168	Anti-Vinculin Mouse mAb	WB/IHC-P
PTM-5224	Anti-Vinculin Rabbit mAb	WB/IHC-P
PTM-5840	Anti-VPS18 Mouse mAb	WB/IHC-P/ICC/IF
PTM-6039	Anti-WWOX Mouse mAb	WB/ICC/IF
PTM-5935	Anti-WWOX Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5477	Anti-XIAP Mouse mAb	WB/ICC/IF

Cat.#	Product	Recommended Applications
PTM-5405	Anti-XRCC5 Mouse mAb	WB/ICC/IF
PTM-5165	Anti-ZO-1 Mouse mAb	IHC-P

Loading Control Antibodies

Cytoplasm

PTM-6702	Anti-Actin (Plant) Mouse mAb	WB
PTM-5001	Anti-alpha Tubulin Mouse mAb	WB/IHC-P/ICC/IF
PTM-5282	Anti-alpha Tubulin Mouse mAb	WB/ICC/IF
PTM-5442	Anti-alpha Tubulin Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5436	Anti-beta Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC
PTM-5028	Anti-beta Actin Rabbit mAb	WB/IHC-P/ICC/IF/FC/IP
PTM-6622	Anti-beta Tubulin (HRP) Rabbit mAb	WB/IHC-P
PTM-6414	Anti-beta Tubulin Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-6620	Anti-GAPDH (HRP) Rabbit mAb	WB/IHC-P
PTM-6684	Anti-GAPDH Mouse mAb	WB/IHC-P/ICC/IF/FC
PTM-5375	Anti-GAPDH Rabbit mAb	WB/IHC-P/ICC/IF/IP
PTM-5097	Anti-HSP90 beta Mouse mAb	WB/IHC-P
PTM-5734	Anti-HSP90 beta Rabbit mAb	WB/IHC-P/IP

Nucleus

PTM-6621	Anti-Histone H3 (HRP) Rabbit mAb (NT)	WB/IHC-P
PTM-1002RM	Anti-Histone H3 Rabbit mAb (CT)	WB/IHC-P/IP/ICC/IF
PTM-1001RM	Anti-Histone H3 Rabbit mAb (NT)	WB/IHC-P/ICC/IF/FC/ChIP
PTM-1002	Anti-Histone H3 Rabbit pAb (CT)	WB/IHC-P
PTM-1009	Anti-Histone H4 Mouse mAb (NT)	WB/IHC-P
PTM-1015RM	Anti-Histone H4 Rabbit mAb (CT)	WB/IHC-P
PTM-5187	Anti-PCNA Mouse mAb	WB/IHC-P
PTM-6409	Anti-HDAC1 Rabbit mAb	WB/IHC-P/ICC/IF/IP

Nuclear Membrane

PTM-5495	Anti-Lamin B1 Rabbit mAb	WB/IHC-P/mIHC-P
----------	--------------------------	-----------------

Whole Blood / Serum / Plasma

PTM-5794	Anti-Albumin Mouse mAb	WB
----------	------------------------	----

Cat.#

Product

Recommended Applications

Mitochondria

PTM-6675	Anti-COX IV (HRP conjugate) Rabbit mAb	WB/IHC-P
PTM-6421	Anti-COX IV Rabbit mAb	WB/IHC-P/FC/IP
PTM-6157	Anti-VDAC1/Porin Rabbit mAb	WB/IHC-P/ICC/IF/FC

Membrane Protein

PTM-6271	Anti-Na ⁺ /K ⁺ ATPase alpha-1 Rabbit mAb	WB/IHC-P/ICC/IF
PTM-5009	Anti-ATPB Mouse mAb	WB/IHC-P/ICC/IF

Cell Junction

PTM-6939	Anti-Vinculin Rabbit mAb	WB/IHC-P
----------	--------------------------	----------

Isotype Controls

PTM-5284	Mouse mAb IgG1 Isotype Control	FC
PTM-5283	Mouse mAb IgG1 Isotype Control	IHC-P/ICC/IF/IP/ChIP
PTM-6984	Mouse mAb IgG1 Isotype Control (APC Conjugate)	ICC/IF/FC
PTM-6136	Mouse mAb IgG1 Isotype Control (FITC Conjugate)	ICC/IF/FC
PTM-6135	Mouse mAb IgG1 Isotype Control (PE Conjugate)	ICC/IF/FC
PTM-5073	Rabbit mAb IgG Isotype Control	IHC-P/ICC/IF/IP/ChIP
PTM-6628	Rabbit mAb IgG Isotype Control	IHC-P/ICC/IF/FC
PTM-6968	Rabbit mAb IgG Isotype Control (APC)	ICC/IF/FC

Secondary Antibodies

PTM-6261	Goat Anti-Rabbit IgG H&L (HRP) pAb	ELISA/WB/IHC-P
PTM-6252	Goat Anti-Rabbit IgG H&L pAb	WB/ELISA/ICC/IF
PTM-7168	Goat Anti-Mouse IgG H&L (HRP) pAb	WB/IHC-P


Tag Antibodies

Cat.#	Product	Recommended Applications
PTM-6075	Anti-DYKDDDDK Tag Mouse mAb (Binds to FLAG® Tag Sequence)	WB/IP/ChIP
PTM-5577	Anti-DYKDDDDK Tag Rabbit mAb (Binds to FLAG® Tag Sequence)	WB/ChIP
PTM-5180	Anti-GFP Mouse mAb	WB
PTM-6673	Anti-GFP Rabbit mAb	WB
PTM-5181	Anti-GST Tag Mouse mAb	WB
PTM-5046	Anti-GST Tag Mouse mAb	WB
PTM-5177	Anti-HA tag Mouse mAb	WB
PTM-5389	Anti-HA tag Rabbit mAb	WB/IP/ICC/IF/ChIP
PTM-5614	Anti-His tag Rabbit mAb	WB/ICC/IF/IP
PTM-5619	Anti-mCherry Mouse mAb	WB
PTM-5390	Anti-Myc tag Rabbit mAb	WB/ICC/IF/IP

Tag Antibodies

Anti-GFP Rabbit mAb

PTM-6673




1.293T(Neg)

2.293T transfected with GFP

Anti-Flag Tag Rabbit mAb

PTM-5577



1.293T(Neg)


2.293T transfected with PFKFB3

Published Case Studies

Cell Metabolism

IF=29.0

The diapause-like colorectal cancer cells induced by SMC4 attenuation are characterized by low proliferation and chemotherapy insensitivity




Cat.# Product

PTM-1401	Anti-L-Lactyl Lysine Rabbit pAb
PTM-1409	Anti-L-Lactyl-Histone H4 (Lys5) Mouse mAb
PTM-1414RM	Anti-L-Lactyl-Histone H3 (Lys14) Rabbit mAb
PTM-1411RM	Anti-L-Lactyl-Histone H4 (Lys12) Rabbit mAb
PTM-1415RM	Anti-L-Lactyl-Histone H4 (Lys8) Rabbit mAb
PTM-1406RM	Anti-L-Lactyl-Histone H3 (Lys18) Rabbit mAb

Cell Research

IF=44.1

α -myosin heavy chain lactylation maintains sarcomeric structure and function and alleviates the development of heart failure




Cat.# Product

Custom Antibody	anti- α -MHC K1897 lac antibody
PTM-1404	Anti-L-lactyllysine antibody conjugated agarose beads
PTM-1401RM	Anti-L-Lactyl Lysine Rabbit mAb

nature cell biology

IF=21.3

De novo pyrimidine biosynthetic complexes support cancer cell proliferation and ferroptosis defence



Cat.#


Product

Custom Antibody	Anti-phospho-UMPS S214 antibody
-----------------	---------------------------------

nature metabolism

IF=20.8

Glutarate regulates T cell metabolism and anti-tumour immunity



Cat.#


Product

PTM-1151	Anti-Glutaryllysine Rabbit pAb
----------	--------------------------------

nature neuroscience

IF=25.0

Emergence of consciousness from anesthesia through ubiquitin degradation of KCC2 in the ventral posteromedial nucleus of the thalamus



Cat.#


Product

PTM-1106	Anti-Ubiquitin Rabbit pAb
----------	---------------------------

Published Case Studies

Cell IF=64.5


Metabolic regulation of homologous recombination repair by MRE11 lactylation



Cat.#	Product
PTM-1401RM	Anti-L-Lactyl Lysine Rabbit mAb
PTM-101	Anti-Acetyllysine Mouse mAb
Custom Antibody	Anti-MRE11-K673la antibody

nature IF=64.8


Diverse modes of H3K36me3-guided nucleosomal deacetylation by Rpd3S



Cat.#	Product
PTM-157	Anti-Acetyl-Histone H3 (Lys14) Mouse mAb
PTM-158	Anti-Acetyl-Histone H3 (Lys18) Mouse mAb
PTM-115	Anti-Acetyl-Histone H3 (Lys23) Rabbit pAb
PTM-160	Anti-Acetyl-Histone H3 (Lys27) Mouse mAb
PTM-162	Anti-Acetyl-Histone H3 (Lys56) Mouse mAb
PTM-625	Anti-Tri-Methyl-Histone H3 (Lys36) Rabbit pAb
PTM-1009	Anti-Histone H4 Mouse mAb (NT)

Cell IF=64.5


Linking chromatin acylation mark-defined proteome and genome in living cells



Cat.#	Product
PTM-545RM	Anti-Crotonyl-Histone H3 (Lys27) Rabbit mAb
PTM-1293	Anti-β-Hydroxybutyryl-Histone H3 (Lys27) Mouse mAb
PTM-502	Anti-Crotonyllysine Mouse mAb
PTM-801	Anti-2-Hydroxyisobutyryllysine Rabbit pAb
PTM-1201	Anti-β-Hydroxybutyryllysine Rabbit pAb

nature IF=64.8


Lysine catabolism reprograms tumour immunity through histone crotonylation



Cat.#	Product
PTM-501	Anti-Crotonyllysine Rabbit pAb

nature IF=64.8

Lysine catabolism reprograms tumour immunity through histone crotonylation



Cat.#	Product
PTM-6055	Anti-SUMO-1 Rabbit mAb (C-terminal)

Published Case Studies

PMID	Title	Journal
37468628	Diverse modes of H3K36me3-guided nucleosomal deacetylation by Rpd3S	<i>Nature</i>
37198486	Lysine catabolism reprograms tumour immunity through histone crotonylation	<i>Nature</i>
36921622	Lactate regulates cell cycle by remodelling the anaphase promoting complex	<i>Nature</i>
31645732	Metabolic regulation of gene expression by histone lactylation	<i>Nature</i>
32968279	Modulating TRADD to restore cellular homeostasis and inhibit apoptosis	<i>Nature</i>
23535595	Succinate is an inflammatory signal that induces IL-1 β through HIF-1 α	<i>Nature</i>
36099882	Acquired semi-squamization during chemotherapy suggests differentiation as a therapeutic strategy for bladder cancer	<i>Cancer Cell</i>
30008322	SET1a-mediated mono-methylation at K342 regulates YAP activation by blocking its nuclear export and promotes tumorigenesis.	<i>Cancer Cell</i>
38653238	Alanyl-tRNA synthetase, AARS1, is a lactate sensor and lactyltransferase that lactylates p53 and contributes to tumorigenesis	<i>Cell</i>
38128537	Metabolic regulation of homologous recombination repair by MRE11 lactylation	<i>Cell</i>
36868209	Linking chromatin acylation mark-defined proteome and genome in living cells	<i>Cell</i>
34384544	Erythroid mitochondrial retention triggers myeloid-dependent type I interferon in human SLE	<i>Cell</i>
34174187	The single-cell epigenomic and transcriptional landscape of immunity to influenza vaccination	<i>Cell</i>
21925322	Identification of 67 histone marks and histone lysine crotonylation as a new type of histone modification	<i>Cell</i>
29706550	Single-cell chromatin modification profiling reveals increased epigenetic variations with aging	<i>Cell</i>
37777506	STAT5 promotes PD-L1 expression by facilitating histone lactylation to drive immunosuppression in acute myeloid leukemia	<i>STTT</i>
36872368	Phosphorylation of PHF2 by AMPK releases the repressive H3K9me2 and inhibits cancer metastasis	<i>STTT</i>
31637018	Global profiling of O-GlcNAcylated and/or phosphorylated proteins in hepatoblastoma	<i>STTT</i>
32200651	Branched-chain amino acid catabolism promotes thrombosis risk by enhancing tropomodulin-3 propionylation in platelets.	<i>Circulation</i>
37443257	α -myosin heavy chain lactylation maintains sarcomeric structure and function and alleviates the development of heart failure	<i>Cell Res</i>
36882514	Lysine 2-hydroxyisobutyrylation of NAT10 promotes cancer metastasis in an ac4C-dependent manner	<i>Cell Res</i>
29192674	Landscape of the regulatory elements for lysine 2-hydroxyisobutyrylation pathway	<i>Cell Res</i>
28497810	Class I histone deacetylases are major histone decrotonylases: evidence for critical and broad function of histone crotonylation in transcription	<i>Cell Res</i>
37789393	LINC00922 decoys SIRT3 to facilitate the metastasis of colorectal cancer through up-regulation the H3K27 crotonylation of ETS1 promoter	<i>Mol Cancer</i>
37684641	CircXRN2 suppresses tumor progression driven by histone lactylation through activating the Hippo pathway in human bladder cancer	<i>Mol Cancer</i>
35303422	Positive feedback regulation of microglial glucose metabolism by histone H4 lysine 12 lactylation in Alzheimer's disease	<i>Cell Metab</i>

PMID	Title	Journal
38232735	Augmentation of scleral glycolysis promotes myopia through histone lactylation	<i>Cell Metab</i>
37979583	The tRNA-GCN2-FBXO22-axis-mediated mTOR ubiquitination senses amino acid insufficiency	<i>Cell Metab</i>
37797623	High dietary fructose promotes hepatocellular carcinoma progression by enhancing O-GlcNAcylation via microbiota-derived acetate	<i>Cell Metab</i>
37543034	The diapause-like colorectal cancer cells induced by SMC4 attenuation are characterized by low proliferation and chemotherapy insensitivity	<i>Cell Metab</i>
37463576	Hexokinase 2-mediated gene expression via histone lactylation is required for hepatic stellate cell activation and liver fibrosis	<i>Cell Metab</i>
34879239	NEAT1 is essential for metabolic changes that promote breast cancer growth and metastasis	<i>Cell Metab</i>
32275862	Pimozide alleviates hyperglycemia in diet-induced obesity by inhibiting skeletal muscle ketone oxidation.	<i>Cell Metab</i>
30146488	Quiescent endothelial cells upregulate fatty acid β -oxidation for vasculoprotection via redox homeostasis.	<i>Cell Metab</i>
30146486	Impairment of angiogenesis by fatty acid synthase inhibition involves mTOR malonylation.	<i>Cell Metab</i>
29719225	A potent and specific CD38 inhibitor ameliorates age-related metabolic dysfunction by reversing tissue NAD ⁺ decline.	<i>Cell Metab</i>
28380375	A class of reactive acyl-CoA species reveals the non-enzymatic origins of protein acylation.	<i>Cell Metab</i>
24703693	Lysine glutarylation is a protein posttranslational modification regulated by SIRT5.	<i>Cell Metab</i>
38679071	Nucleolin lactylation contributes to intrahepatic cholangiocarcinoma pathogenesis via RNA splicing regulation of MADD	<i>J Hepatol</i>
38703775	Glucose-driven histone lactylation promotes the immunosuppressive activity of monocyte-derived macrophages in glioblastoma	<i>Immunity</i>
36973513	Emergence of consciousness from anesthesia through ubiquitin degradation of KCC2 in the ventral posteromedial nucleus of the thalamus	<i>Nat Neurosci</i>
35395187	A small-molecule cocktail promotes mammalian cardiomyocyte proliferation and heart regeneration	<i>Cell Stem Cell</i>
33450185	Histone crotonylation promotes mesoendodermal commitment of human embryonic stem cells	<i>Cell Stem Cell</i>
38531969	Deciphering functional roles of protein succinylation and glutarylation using genetic code expansion	<i>Nat Chem</i>
38448615	Regulation of urea cycle by reversible high-stoichiometry lysine succinylation	<i>Nat Metab</i>
38413806	Histone butyrylation in the mouse intestine is mediated by the microbiota and associated with regulation of gene expression	<i>Nat Metab</i>
37605057	Glutarate regulates T cell metabolism and anti-tumour immunity	<i>Nat Metab</i>
32839595	Glis1 facilitates induction of pluripotency via an epigenome–metabolome–epigenome signalling cascade	<i>Nat Metab</i>
32694864	SIRT7 couples light-driven body temperature cues to hepatic circadian phase coherence and gluconeogenesis	<i>Nat Metab</i>
37481436	Martynoside rescues 5-fluorouracil-impaired ribosome biogenesis by stabilizing RPL27A	<i>Sci Bull</i>
37291265	De novo pyrimidine biosynthetic complexes support cancer cell proliferation and ferroptosis defence	<i>Nat Cell Biol</i>
31871320	Ketogenesis-generated β -hydroxybutyrate is an epigenetic regulator of CD8 ⁺ T-cell memory development	<i>Nat Cell Biol</i>

PMID	Title	Journal
36646829	A natural allele of proteasome maturation factor improves rice resistance to multiple pathogens	<i>Nat Plants</i>
36513638	CENP-F-dependent DRP1 function regulates APC/C activity during oocyte meiosis I	<i>Nat Commun</i>
36414671	SEN1 prevents steatohepatitis by suppressing RIPK1-driven apoptosis and inflammation	<i>Nat Commun</i>
36333310	YiaC and CobB regulate lysine lactylation in <i>Escherichia coli</i>	<i>Nat Commun</i>
36253417	Sirtuin5 protects colorectal cancer from DNA damage by keeping nucleotide availability	<i>Nat Commun</i>
35927268	Deficiency of WTAP in hepatocytes induces lipotrophy and non-alcoholic steatohepatitis (NASH)	<i>Nat Commun</i>
35449131	KMT5A-methylated SNIP1 promotes triple-negative breast cancer metastasis by activating YAP signaling	<i>Nat Commun</i>
35296687	Global profiling of regulatory elements in the histone benzoylation pathway	<i>Nat Commun</i>
35013160	Altered succinylation of mitochondrial proteins, APP and tau in Alzheimers disease	<i>Nat Commun</i>
38670996	HBO1 catalyzes lysine lactylation and mediates histone H3K9la to regulate gene transcription	<i>Nat Commun</i>
38615060	Lipopolysaccharide binding protein resists hepatic oxidative stress by regulating lipid droplet homeostasis	<i>Nat Commun</i>
38438368	PsAF5 functions as an essential adapter for PsPHB2-mediated mitophagy under ROS stress in <i>Phytophthora sojae</i>	<i>Nat Commun</i>
38413566	Deaggregation of mutant <i>Plasmodium yoelii</i> de-ubiquitinase UBP1 alters MDR1 localization to confer multidrug resistance	<i>Nat Commun</i>
38216551	Mutant KRAS-activated circATXN7 fosters tumor immunoescape by sensitizing tumor-specific T cells to activation-induced cell death	<i>Nat Commun</i>
38167340	Histone lactylation couples cellular metabolism with developmental gene regulatory networks	<i>Nat Commun</i>
37863889	Lactylation of METTL16 promotes cuproptosis via m6A-modification on FDX1 mRNA in gastric cancer	<i>Nat Commun</i>
37803016	Decreased liver B vitamin-related enzymes as a metabolic hallmark of cancer cachexia	<i>Nat Commun</i>
37737243	A chemical catalyst enabling histone acylation with endogenous acyl-CoA	<i>Nat Commun</i>
37277331	ACL and HAT1 form a nuclear module to acetylate histone H4K5 and promote cell proliferation	<i>Nat Commun</i>
37198164	Small-sample learning reveals propionylation in determining global protein homeostasis	<i>Nat Commun</i>
36914626	De novo lipogenesis fuels adipocyte autophagosome and lysosome membrane dynamics	<i>Nat Commun</i>
34654833	The δ subunit of F1Fo-ATP synthase is required for pathogenicity of <i>Candida albicans</i>	<i>Nat Commun</i>
34545082	Isonicotinylation is a histone mark induced by the anti-tuberculosis first-line drug isoniazid	<i>Nat Commun</i>
33741957	The essential role of PRAK in tumor metastasis and its therapeutic potential	<i>Nat Commun</i>
30154464	Lysine benzoylation is a histone mark regulated by SIRT2	<i>Nat Commun</i>
30026585	SIRT7 has a critical role in bone formation by regulating lysine acylation of SP7/Osterix	<i>Nat Commun</i>

PMID	Title	Journal
29416026	Sirtuin5 contributes to colorectal carcinogenesis by enhancing glutaminolysis in a deglutarylation-dependent manner	<i>Nat Commun</i>
29317660	Microbiota derived short chain fatty acids promote histone crotonylation in the colon through histone deacetylases	<i>Nat Commun</i>
29158494	Mrg15 stimulates Ash1 H3K36 methyltransferase activity and facilitates Ash1 Trithorax group protein function in Drosophila	<i>Nat Commun</i>
29070843	Characterization of histone acylations links chromatin modifications with metabolism	<i>Nat Commun</i>
29030587	NatD promotes lung cancer progression by preventing histone H4 serine phosphorylation to activate Slug expression	<i>Nat Commun</i>
33311474	Ubiquitination of RIPK1 regulates its activation mediated by TNFR1 and TLRs signaling in distinct manners	<i>Nat Commun</i>
31113955	Ubiquitylome study identifies increased histone 2A ubiquitylation as an evolutionarily conserved aging biomarker	<i>Nat Commun</i>
30659183	Malonylation of GAPDH is an inflammatory signal in macrophages	<i>Nat Commun</i>
26158509	Biotin starvation causes mitochondrial protein hyperacetylation and partial rescue by the SIRT3-like deacetylase Hst4p	<i>Nat Commun</i>
27436229	SIRT7 is a histone desuccinylase that functionally links to chromatin compaction and genome stability	<i>Nat Commun</i>
35738262	TCR activation directly stimulates PYGB-dependent glycogenolysis to fuel the early recall response in CD8+ memory T cells	<i>Mol Cell</i>
35447080	CDYL1-dependent decrease in lysine crotonylation at DNA double-strand break sites functionally uncouples transcriptional silencing and repair	<i>Mol Cell</i>
35320754	Lactylation-driven METTL3-mediated RNA m6A modification promotes immunosuppression of tumor-infiltrating myeloid cells	<i>Mol Cell</i>
34856123	Quantitative subcellular acyl-CoA analysis reveals distinct nuclear metabolism and isoleucine-dependent histone propionylation	<i>Mol Cell</i>
38423014	A di-acetyl-decorated chromatin signature couples liquid condensation to suppress DNA end synapsis	<i>Mol Cell</i>
38176415	OXCT1 functions as a succinyltransferase, contributing to hepatocellular carcinoma via succinylating LACTB	<i>Mol Cell</i>
37311463	Histone H3 lysine 27 crotonylation mediates gene transcriptional repression in chromatin	<i>Mol Cell</i>
37244254	IL-1 β -associated NNT acetylation orchestrates iron-sulfur cluster maintenance and cancer immunotherapy resistance	<i>Mol Cell</i>
32937100	Histone deacetylase 3 couples mitochondria to drive IL-1 β -dependent inflammation by configuring fatty acid oxidation.	<i>Mol Cell</i>
23746352	SIRT4 coordinates the balance between lipid synthesis and catabolism by repressing malonyl CoA decarboxylase.	<i>Mol Cell</i>
28803779	Chromodomain protein CDYL acts as a crotonyl-CoA hydratase to regulate histone crotonylation and spermatogenesis.	<i>Mol Cell</i>
31676231	Recognition of histone crotonylation by Taf14 links metabolic state to gene expression.	<i>Mol Cell</i>
31000437	Regulation of UCP1 and mitochondrial metabolism in brown adipose tissue by reversible succinylation.	<i>Mol Cell</i>
29775581	p300-mediated lysine 2-hydroxyisobutyrylation regulates glycolysis.	<i>Mol Cell</i>
25818647	Intracellular crotonyl-CoA stimulates transcription through p300-catalyzed histone crotonylation.	<i>Mol Cell</i>
27105113	Dynamic competing histone H4 K5K8 acetylation and butyrylation are hallmarks of highly active gene promoters.	<i>Mol Cell</i>


Company Profile

PTM BIO is a platform-driven biotechnology company integrating R&D, production, and technical services. With mass spectrometry-based proteomics as its core technology and antibody reagents as a key product pillar, PTM BIO provides integrated solutions covering protein discovery, validation, and detection.


PTMab® is PTM BIO's premium antibody brand, offering high-quality antibody reagents and custom antibody development services for research, pharmaceutical R&D, and IVD applications. Its product portfolio includes pan-PTM antibodies, histone modification antibodies, site-specific antibodies, and recombinant proteins, supporting studies in epigenetics, tumor immunology, metabolism, neuroscience, and diagnostics.

CUT&Tag-Validated Antibodies — Newly Launched


PTM-1427RM Anti-L-Lactyl-Histone H3 (Lys18) Rabbit mAb-ChIP Grade




PTM-1415RM Anti-L-Lactyl-Histone H4 (Lys8) Rabbit mAb




Strict Validation. Exceptional Specificity



① Dot Blot Validation



② Drug Treatment Validation



③ Blocking Peptide

